



SIMPLIFIED DIABETES MANAGEMENT.

EVENCARE® G2™
Blood Glucose Monitoring System

EVENCARE® G2™

Blood Glucose Monitoring System with Voice Guidance

Advanced technology makes it easier to care
for residents with diabetes.

Easy-to-use

Voice assistance guides you
through the testing process and
provides you with your results.

Auto code system

Minimizes errors and
inaccurate results.

Fast test time

Results in 6 seconds.

Small sample size

Small blood requirement, only
1.0uL, means less discomfort.

Alternate site testing

Option to test on the forearm,
palm or fingertips.

Memory option

Stores 300 results and can provide
7-, 14- and 30-day averaging.

with Voice Guidance

"103 mg/dL"



Item No.	Description	Pkg
MPH1540	Glucose Meter	1 Each
MPH1550	Test Strips	50 ea/bx, 12 bx/cs
MPH1550Z	Test Strips	50 ea/bx
MPH1560	Control Solution High/Low	6 bx/cs
MPH1560Z	Control Solution High/Low	1 bx
MPHBAAAZ	Replacement Battery, AAA	24 ea/bx

Contact your Medline representative to request
a product trial.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)