



## Portable Suction Machine AC and battery power

**Item#**  
**SU100DC**

Includes internal rechargeable 12-volt lithium battery, which can be used for up to 1 hour of operation per charge.

## Suction Machine AC power

**Item#**  
**SU100AC**

### Both machines include

- Machine
- Carry Bag
- 800cc Suction Canister
- AC Adapter & Power Cord
- 6ft Suction Tubing
- 3.5in Suction Tube Connector
- 4.5in Suction Tube Connector
- Bacteria Filter

### Both machines feature

- Effective removal of respiratory fluids
- Durability
- Variable vacuum control
- Quiet and easy to use
- Two-year warranty\*
- HCPCS# E0600

\*The SU100DC internal rechargeable battery is covered by a 90-day warranty

### Specifications

Specialty	Small/Lightweight
Carry Case	Yes
Dimensions (L x H x W in)	12.83" x 6.37" x 8.12"
Weight (lbs.)	SU100DC: 5.29 lbs / SU100AC: 4.08 lbs
Vacuum Range (mm Hg)	50-560
LPM	25
Sound Level (dBA)	60 dBA
Canister Type w/ purchase of machine	Disposable
Uses Suction Canister Size (cc)	800 /1000/1200
Machine Warranty	2 years
Can be operated inside carry case?	Yes

### Accessories and replacement parts

<b>RES024-MKIT</b>	Replacement Filter Kit (connectors and filter)
<b>RES026S-SS</b>	Replacement Suction Kit (canister, tubing, connectors, and filter)
<b>SU100DC-Adapt</b>	Replacement Power Adapter (without cord)
<b>SU100DC-Cord</b>	Replacement Power Cord (without adapter)
<b>SU100DC-BAG</b>	Carrying/Storage Bag
<b>SU100DC-BATTERY</b>	Replacement internal rechargeable battery for SU100DC



SCAN TO LEARN MORE



## Suction Kit for SU100DC

Item# **RES026S-SS**

RES023AW Suction Canister  
RES024 Suction Tubing Connector 3.5in  
RES024M Suction Tubing Connector 4.5in  
BF1438 Bacteria Filter  
RES025 Suction Tubing 1/4in x 6ft

HCPCS#  
CANISTER - A7000  
TUBING - A7002



## Suction Canister

Item# **RES023AW**

800cc • Single Use  
48/case  
HCPCS# A7000

Item# **RES025AW1200**

1200cc • Single Use  
48/case


Float valve shutoff



## Suction Filter Component Kit

Item# **RES024XL-3KIT**

RES024XL Suction Tubing Connector 16in  
BF3812 Filter  
RES028 Elbow Connector



## Suction Filter Component Kit for SU100DC

Item# **RES024-MKIT**

RES024 Suction Tubing Connector 3.5in  
RES024M Suction Tubing Connector 4.5in  
BF1438S Bacteria Filter



## Suction Kit

Item# **RES026S-E**

RES023AW Suction Canister  
RES024M Suction Tubing Connector 4.5in  
RES025 Suction Tubing 6ft  
RES028 Elbow Connector  
BF3812 Bacteria Filter

HCPCS#  
CANISTER - A7000  
TUBING - A7002



**To order or for more information**

P 312.997.9980 E [sales@sunsethcs.com](mailto:sales@sunsethcs.com) W [SUNSETHCS.COM](http://SUNSETHCS.COM)

