

WORLD'S LIGHTEST POWERCHAIR!

Feather POWER CHAIR

At only 33 lbs.*, the Feather Power Chair is the lightest power wheelchair on the planet.

This lightweight, airline-compliant, **one-piece folding** power chair makes life on the move a breeze. The Feather Power Chair folds to 13" x 28" x 27", making it super easy for storage and travel. At a top speed of 4 MPH it can keep up with a brisk walk, while an 11 mile battery range keeps going all day. **Flip-back armrests** make sliding in and out easy, and an optional seat belt and anti-tippers are available for added safety.

FEATURES:

- just 33 lbs*
- 11 mile drive range
- folding frame
- airline compliant battery^




SPECIFICATIONS

Weight Capacity	250 lbs.	Ground Clearance	3"
Product Weight w/out batteries	33 lbs.	Front Wheel Size	7"
Product Weight w/batteries	37 lbs.	Rear Wheel Size	10"
Seat Width	18"	Turning Radius	35"
Seat Depth	16"	Max Climbing Angle	10 degrees
Seat-to-Floor Height	19"	Max Speed	4 mph
Overall Width	24"	Driving Range	11 Miles
Overall Length	34"	Battery Type	24V 8Ah
Overall Height	34"	Warranty on Frame	One Year Limited
Folded Depth	27"	Warranty on Wearables	Six Months Limited
Folded Height	28"	Warranty on Electronics	Six Months Limited
Folded Width	13"	Warranty on Batteries	Three-Months Limited

Lightweight



Foldable



Locking Brakes



Arrives Assembled



*weight without battery. ^removable lithium battery is airline compliant. Check with individual carrier for guidelines.