## WORLD'S LIGHTEST WHEELCHAIR!

## Feather CHAIR

You have to lift it to believe it - it weighs only 13.5 pounds without the rear wheels!



Folds easily so it can go anywhere. Easy to load into the trunk of a car for quick trips around town. The Feather Chair has quick-release wheels and anti-tippers available for safety. Caregiver brakes on the push handles increase safety and convenience.

## **FEATURES:**

- only 13.5 lbs without wheels
- folding frame
- anti-tippers
- quick-release wheels
- optional elevating leg rests

## SPECIFICATIONS

2.20170,4110140			
Weight Capacity	250 lbs.		
Product Weight	19 lbs.	Front Wheel Size	6"
Weight Without Wheels	13.5 lbs.	Rear Wheel Size	22"
Seat Width	18"	Axle Type	Quick Release
Seat Depth	16"	Armrest Options	Full Length
Seat-to-Floor Height	19"	Armrest Height	9"
Overall Width	25"	Back Height	16"
Overall Length	30"	Upholstery Type	Nylon
Overall Height	35"	Upholstery Color	Black w/ Orange Trim
Folded Depth	29"	Footrests	Swingaway
Folded Height	28"	Warranty on Frame	Lifetime Limited
Folded Width	15"	Warranty on Wearables	Six Months Limited

Lightweight



Locking Brakes



**Foldable** 

Arrives Assembled



