

FREESTYLE COMFORT PORTABLE OXYGEN CONCENTRATOR

The FreeStyle Comfort ensures you get the therapy you need with convenience, mobility and peace of mind.



Key clinical features ...

- 5 pulse dose settings deliver up to 1050 mL to address broad range of prescriptions
- **UltraSense®** offers sensitive breath detection ensuring you get oxygen even if your breath becomes shallow
- **autoDOSE** is a safety feature to ensure automatic delivery of oxygen even if no breath is detected to avoid any gaps in your therapy
- **myCAIRE™** is a secure, telehealth solution offering you the ability to take control of your oxygen therapy

... for your comfort

- Large, easy-to-read and glow-in-the-dark LCD display with user-friendly interface and operations
- Ergonomic, curved design allows you to move freely
- Small and lightweight (5 lbs/2.3 kg) so you can go further
- Durable, robust design is built to last

FREQUENTLY ASKED QUESTIONS

Why is a portable oxygen concentrator better than a traditional oxygen cylinder?

A portable oxygen concentrator provides you with more freedom and time away from home based on its ability to operate on both battery or electrical power (AC and DC), while cylinders offer a limited supply of oxygen, requiring you to carry multiple tanks.

Can I travel with my oxygen concentrator on an airplane?

Yes. The Freestyle Comfort meets FAA guidelines for use on commercial flights to anywhere in the world.

How long will my battery last?

Your FreeStyle Comfort battery lasts up to 8 hours with an 8-cell battery and up to 16 hours with an 16-cell battery at setting 1. Battery life depends on the flow rate.

How will I know the battery needs recharging?

Audio and visual notifications will alert you when the device requires charging.

Will I receive a warranty on my FreeStyle Comfort?

Yes, you will receive an industry leading 3-year warranty that covers the device and sieve beds because we stand behind our product and it's reliability.





For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy. CAIRE products are trusted globally by patients, providers, clinicians, and the military. Let us put our experience to work for you.

www.caireinc.com

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.