



HEEL RAISER[®]

Effectively offloads heels.



Every heel needs protection.



PROTECT AND SUSPEND THE HEEL.

Clinical Study Shows Effective Results¹

Using a pressure-sensing mat, 19 patients in a long-term care facility were tested for heel pressure on a mattress and in a Heel Raiser. The Heel Raiser completely offloaded the vulnerable heel area. Here are the pressure mapping conclusions:




Red dotted line indicates heel



Average pressure on contact area **115.23**
Mattress had the highest heel pressure.



Average pressure on contact area **43.83**
Heel Raiser completely offloaded the vulnerable heel area.



Item No.	Description	Pkg
MDT823296	Heel Raiser, Heel Protector	1/ea

BENEFITS

- » Completely offloads the heel.
- » Top folds down to accommodate small feet or for more airflow to the toes.
- » Anatomical strapping pattern for best fit on the patient.
- » Soft comfortable material to help prevent friction and shear.
- » Open heel design promotes air circulation and allows for easy monitoring.
- » Includes mesh laundry bag to simplify washing and sorting.
- » One size fits most.

Ask your Medline representative for a FREE trial in your facility today!



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

REFERENCES 1. Data on file and available upon request.