

## JEANIE RUB® VARIABLE SPEED MASSAGER

- Deep, penetrating massage with a twist of the switch
- Helps reduce pain, muscle tension, and stiffness
- Variable speeds provide soothing massage or more invigorating massage
- An industry standard among professional massage therapists
- Helps stimulate the circulatory system
- Fosters healing of strained muscles and sprained ligaments



### THERAPEUTIC RELIEF

The Jeanie Rub® Variable Speed Massager delivers 1100 - 4600 rpm for a deep, penetrating massage. Help reduce pain, muscle tension, and stiffness at the twist of the switch. The orbital action of the massager provides a softer and deeper massage than "percussion" type massagers. Lower speeds deliver a soothing massage while faster speeds are more invigorating.

### HIGH QUALITY CONSTRUCTION

This product is built with a cast aluminum construction and has a red and black powder coat finish that looks great and resists abrasions, ensuring years of use. Jeanie Rub massagers are built in the USA to ANSI Standards and are CSA certified, ensuring they have been evaluated to meet or exceed consumer safety standards. Pad measures 4" x 9" (10cm x 23cm). Jeanie Rub accessories sold separately. 12' grounded power cord, 120V AC Nominal power source required.

### SPECIFICATIONS

**Product Selection:**  
PRO-3401 Jeanie Rub®  
Variable Speed Massager

**Product Weight:** 7 lbs  
**Dimensions:** 10.25" x 8.25"

**Case Dimensions:**  
19.125" x 10.875" x 8.5"  
**Case Weight:** 28 lbs  
**Case Quantity:** 4

**Latex:** This product is not made with natural rubber latex.



Available From:



## JEANIE RUB® VARIABLE SPEED MASSAGER ACCESSORIES

Professional package includes:

- Jeanie Rub Variable Speed Massager
- Paraspinal Accessory
- Extremity Accessory
- Carrying Bag
- Fleece Cover

*Professional Package  
PRO-3405*



Available From:



### SPECIFICATIONS

Product Selection:	Dimensions:	Wt:	Case Qty:	Case Wt:
PRO-3405: Professional Package	13.5" x 10.25" x 12.5"	11 #	1	11 #
ACC-896 Paraspinal Accessory	10.25" x 5" x 2"	.15 #	3	2 #
ACC-897 Extremity Accessory	10.25" x 5" x 4.5"	.2 #	6	2 #
ACC-884 Wall Hanger	14" x 5" x 5"	1.5 #	1	3 #
ACC-882 Fleece Cover	9.5" x 5.5" x .5"	.1 #	12	4 #
ACC-886 Sheepskin Cover	10" x 8" x 1"	.3 #	6	3 #
ACC-885 Carrying Bag	12" x 8" x 5"	1 #	6	2 #
ACC-893 Terry Cover	10" x 7"	.2 #	3	2 #
ACC-892 Sure-Chek Cover	9" x 4"	.15 #	3	2 #