



## ELK Lifting Cushion

The ELK is a lifting cushion designed to lift from the floor.

### Easy to fit

- Lifting cushion is inflated using an Airflo Compressor with a simple hand control
- Simple to use, requiring very little training
- Suitable for use inside or out
- Can be used by one carer

### Simple to maintain

- Clean using a proprietary non-abrasive mild liquid cleaner.
- Easy to deflate, fold away and store
- The material is designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

### Safe to use

- Lifts the frailest person up to 980 lbs from the floor
- Exceptionally rigid and stable
- Ideal for carers of individuals with neurological disorders who fall regularly

### Improves well-being

- Provides a safe, dignified lift for a fallen person
- Reduces the risk of musculoskeletal injury for the carer
- Can be used in confined spaces where traditional lifting equipment may not fit

### Technical Specifications

Max user weight	980 lbs
Max height	22"
Min height	Flat
Seat depth	20"
Width	22 1/2"
Total product weight	8 lbs
Airflo Compressor	24
Lifting performance per charge	Up to 9 lifts



Product code: MPCA070400

## Features

1. **4-Way Hand Control.**
2. **Made from very strong polyurethane** coated nylon material.
3. **Remove or keep stability bar** to provide extra comfort for the user.
4. **Can be used to lift a stretcher.**
5. **Easy to setup** - connect the four hoses from the hand control to the corresponding numbers and coloured connectors on the Camel.
6. **Powered by an Airflo 24 Compressor.**

## Accessories

- **ELK Stretcher Bars**
- **Slidesheet**
- **ELK End Cap Set**



**For more information visit**

[www.mangarhealth.com](http://www.mangarhealth.com)

Email: [sales@mangarhealth.com](mailto:sales@mangarhealth.com)

Phone: 0800 2800 485

