

# Proactive Medical Air Mattress Selection Guide

Protekt Aire	3000	3500	3600	3600AB	4000DX	4600DX	4600AB	5000DX	6000	6400	6450	6000AB	6500	7000	8000	8600AB	8900	9900	9900AB	Self-Adjusting	
<b>Pump Style</b>	Analog	Analog	Analog	Analog	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Digital	Blower	Blower	(O) Analog	
<b>Low Air Loss</b>	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	
<b>Alternating Pressure</b>	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	O
<b>Pulsation</b>														X			X	X	X		
<b>Lateral Rotation</b>														X							
<b>Self-Adjusting</b>																				X	
<b>3" Fiber Safety Base</b>		X						X			X										
<b>Cell-on-Cell Safety Base</b>			X	X		X	X		X			X	X	X	X	X	X	X	X	X	
<b>Raised Side Rails (foam)</b>	O	O	O		O	O		O	O	O	O		O		O			O		O	
<b>Raised Side Rails (air)</b>				X			X					X		X	O	X	X		X		
<b>Air Flow Per Minute</b>	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	9 LPM	18 LPM	18 LPM	18 LPM	18 LPM	18 LPM	1275 LPM	1275 LPM	9 LPM
<b>Alternating Time(s) min</b>	10	10	10	10	10	10	10	10	10-25	10-25	10-25	10-25	10-25	10-30	10-25	10-25	10-30	5-30	5-30	10	
<b>Static Mode</b>	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
<b>Fowler (Seat Inflate)</b>									X	X	X	X	X	X	X	X	X	X	X	X	
<b>Weight Settings-Dial/LED</b>	Dial	Dial	Dial	Dial	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	LED	Dial
<b>Cover Material</b>	Nylon	Nylon	Nylon	Stretch	Nylon	Nylon	Stretch	Nylon	Nylon	Nylon	Nylon	Stretch	Nylon	Nylon	Nylon	Stretch	Nylon	Nylon	Nylon	Nylon	Stretch
<b>Mattress Height</b>	8"	8"	8"	8"	8"	8"	8"	8"	8"	8"	8"	8"	8"	10"	8"	10"	8"	8"	8"	8"	6"
<b>Mattress Width</b>	36"	36"	36"	36"	36"	36"	36"	36"	36"	36"	36"	36"	36"	36"	36",42"	42"-60"	42",48"	36"	36"-60"	36"-48"	36"-60"
<b>Weight Capacity</b>	350	350	350	350	400	400	400	400	450	450	450	450	450	660	660	660	660	660	1,000	1,000	500
<b>Warranty</b>	18 Month	18 Month	18 Month	18 Month	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	2 Year	5 Year

\*Custom sizes are available \*\*Orders received by 3pm EST ship same day

Private label program available

For more information, please call 855-237-7622

[www.ProactiveMedical.com](http://www.ProactiveMedical.com)