

Nuvo^{Mark 5} *Lite³*



Nuvo^{Mark 5} *Lite*

Nuvo Lite 3

230 Volt
50 Hz

--

935

Standard Model Part Number

OCSI Model Part Number
(Oxygen Concentration Status Indicator)

Nuvo Lite

115 Volt 230 Volt
60 Hz 50 Hz

520

920

525

925

180 Watts

90% (+6.5% / -3%)
at 3 LPM

0.125 to 3 LPM

40 dBA typical

7.0 PSI

32 lbs / 14.5 kg

14"W x 9"D x 23"H
36 cm x 23 cm x 58.5 cm

Average Power Consumption

Oxygen Concentration

Liter Flow

Sound Level

Outlet Pressure

Weight

Dimensions

330 Watts

90% (+6.5% / -3%)
at 5 LPM

0.125 to 5 LPM

45 dBA typical

7.0 PSI

32 lbs / 14.5 kg

14"W x 9"D x 23"H
36 cm x 23 cm x 58.5 cm

Accessories available upon request -- Please contact your sales rep!

Additional Features



Lockable Flow Control Valve



Metal Outlet Fitting



Resettable Circuit Breaker

Nuvo⁸ Mark 5



Nuvo¹⁰ Mark 5

Nuvo 8

115 Volt 230 Volt
60 Hz 50/60 Hz

580

980

Standard Model Part Number

--

--

585

985

OCSI Model Part Number
(Oxygen Concentration Status Indicator)

595

1005

Nuvo 8		Nuvo 10	
115 Volt 60 Hz	230 Volt 50/60 Hz	115 Volt 60 Hz	230 Volt 50 Hz
580	980	--	--
585	985	595	1005
Average Power Consumption		700 Watts	600 Watts
500 Watts	490/625 Watts	Oxygen Concentration	
90% (+5% / -3%) at 8 LPM		90% (+5% / -3%) at 10 LPM	
2 to 8 LPM		2 to 10 LPM	
15 PSIG (1.2 BAR)		20 PSIG (1.38 BAR)	
54 lbs / 24 kg		54 lbs / 24 kg	
16"W x 16"D x 28"H 39 cm x 40 cm x 71 cm		16"W x 16"D x 28"H 39 cm x 40 cm x 71 cm	
Dimensions		Dimensions	

Accessories available upon request -- Please contact your sales rep!

Additional Features

- Built on the Reliable Nuvo Platform
- Easy Access to filter and hour meter
- Superior Grade of Molecular Sieve
- Spring Mounted High Efficiency Compressor



CE
0413

2010-8406 REV C