



PAIN MANAGEMENT  
TECHNOLOGIES



PRODUCT CATALOG

---



In 2000, Pain Management Technologies started with electrotherapy equipment. Over the years we have diversified our product offerings and greatly improved technologically. Today our unmatched levels of customer service, product quality, and innovation makes us a leader in the manufacturing of medical devices and health products.

# 2019 2020

---

## PRODUCT CATALOG

### OUR GUARANTEE

- 1 Provide quality products at unbeatable value
- 2 Give you the highest level of care and concern
- 3 Provide immediate shipping and prompt delivery

### PRODUCT OFFERINGS

#### **ELECTROTHERAPY**

Wide range of portable devices made with state-of-the-art technology

#### **ULTRASOUND**

Wide range of portable ultrasound devices including our new TENS/ Ultrasound combo device

#### **ELECTRODES AND CONDUCTIVE GARMENTS FOR ELECTROTHERAPY**

Unique, versatile body wrap system that can target a specific treatment site with great accuracy

#### **LIGHT THERAPY**

From Skin to pain relief, light therapy provides results

#### **WATER THERAPY SYSTEMS**

Hot/cold water therapy system, and a cold-only system with therapy wraps to treat all body parts

#### **BRACING AND SUPPORTS**

Diverse range of LSO and TSLO bracing as well as ROM knee braces, and soft goods

#### **INFRARED HEATING PADS**

True infrared carbon fiber fabric heating pads, and our Thermotech automatic moist heating pads

#### **SEAT CUSHIONS**

Mobility assistive and comfortable seat cushions

#### **KINESIOLOGY TAPE**

Speeds healing by slightly lifting skin away from sore or injured tissues

#### **TOPICAL ANALGESICS**

Highly effective, available in generic and name-brand formulations

#### **TRACTION**

Portable and at-home traction devices

#### **MASSAGE**

Self massage tools and powered massage tools

# iTENS

iTENS™ works by sending tiny electrical signals through the skin to intercept pain signals from reaching the brain and helps release endorphins, the body's natural pain-fighting chemical.



**LARGE WINGS**  
for back, shoulder and large treatment areas

AVAILABLE IN:  
**BLUE**  
ITENSLGBLUE



**GREY**  
ITENSLGGREY



**SMALL WINGS**  
for ankle, knee, elbow, wrist and other small treatment areas

AVAILABLE IN:  
**RED**  
ITENSSMRED



**GREEN**  
ITENSSMGREEN



**LONG STRIP**  
for muscle joints and long treatment areas

AVAILABLE IN:  
**PURPLE**  
ITSTRIPPR



**WHITE**  
ITSTRIPWH



DESCRIPTION	SIZE	QUANTITY	PRODUCT CODE
SMALL WINGS	1.75" x 5.5"	3 PAIRS / BOX	ITRODESM
LARGE WINGS	2.5" x 6.5"	3 PAIRS / BOX	ITRODELG
LONG STRIP	1.5" x 7"	3 PAIRS / BOX	ITRODESTR

\*Available in bundle packs.



#### ACCESSORIES

1. Charger, charging cord
2. AC adaptor (sold separately)



Body part specific settings



Track your results and chart your progress




Ultra-conductive laser printed silver

#### WHAT IS TENS?

T.E.N.S stands for Transcutaneous Electrical Nerve Stimulation. For the past Thirty years, doctors have been prescribing TENS as an effective way to reduce or even eliminate pain. The iTENS communicates with your nerve endings to suppress pain via two theories:

- Gate Control Theory - overriding the nerve pain signals before they reach the brain
- Endorphin release theory - releasing the body's own natural pain fighting chemicals called endorphins.

#### Now with rotating wing technology.





### Ultima 3t

#### THE ULTIMA 3T TENS

A dual-channel device with three modes of operation (B, N, and M). Its pulse width and pulse frequency are adjustable. It is also equipped with a 30, 60, 90 minute timer, and has a protective cover over the controls. The U3t is ultra slim and runs on AAA batteries that provide 2 TIMES longer life span than 9volt batteries. Over 850 continuous minutes of battery life.

PRODUCT CODE: U3TN



### Ultima 5

#### THE ULTIMA 5 TENS

A dual-channel five-mode device with two wave form adjustments and a large LCD screen. This device is remarkable for its function, and traditional looks mixed with new age digital style. Also includes a soft carrying case. The Ultima Five is a quality choice at a valued price.

PRODUCT CODE: U5



### Ultima 20

#### THE ULTIMA 20 TENS

A 20-mode device with dual channels and three wave forms. The ULTIMA 20 offers 12 preset modes and 8 manual modes. Also includes a body part diagram for plug-in-play settings.

PRODUCT CODE: U20B



#### THE USTIM

A dual-channel device for muscle exercise and treatment of physical injury. It has 12 modes of operation that are user friendly and written in practical terms. It has two wave forms and a large LCD screen. It is by far one of the best portable EMS devices on the market.

PRODUCT CODE: USTIM



#### ULTIMA COMBO

The Ultima Combo is our premier combo stimulator. It combines the technology of the TENS U20, and the EMS Ustim, both operating with a detailed and comprehensive body-part diagram.

- Dual channel device
- LED Backlight
- Patient compliance timer
- Digital body part diagram

PRODUCT CODE: UCOMBO

# NerveSpa

## ADVANCED NEUROPATHY STIMULATOR

The Nerve Spa is an Advanced Neuropathy stimulator that utilizes enhanced micro-current specifications coupled with a monophasic wave form to target pain associated with early onset peripheral neuropathy in the hands and feet.

The system is designed to deliver an aquatic electrotherapy stimulation that surrounds the feet with a uniform level of current.

System includes device, electrodes, leadwires, foot bath, conductive salt, and effervescent tablets.

PRODUCT CODE: NSFB100





#### NERVE SPA FOOT BATH SUPPLY KIT

Contains 1 salt container and 2 tubes of effervescent tubes.

PRODUCT CODES:  
NSKIT



**STRETCH** by NerveSpa

#### PERFORMANCE DIABETIC SOCK

Not until now, has there been a high quality diabetic sock.

PRODUCT CODES:  
Knee high: **SDSK**  
Ankle high: **SDSA**

#### NERVE REBUILDER FORMULA BY NERVE SPA PERIPHERAL NEUROPATHY SUPPLEMENT

Manage pain, inflammation and the loss of sensation due to peripheral neuropathy. Consisting of a blend of amino acids, nutrients, and botanicals that work in concert to increase and restore balance in the nervous system.

- 2 capsules per serving
- 120 capsules per bottle

PRODUCT CODES:  
NRB100





### ADVANCED MULTI-MODE STIMULATOR

The Ultima Neo is an multi-mode electrotherapy stimulator that combines enhanced TENS, EMS, Interferential, and Microcurrent. It is a dual channel device with a plug-in-play body part diagram.

The Ultima Neo operates on an high capacity lithium-ion battery which allows the user to use the Interferential function without being tethered to the wall, as well as long battery runs on the TENS/EMS and Micro modes. A charging plug and universal leadwires are included.

PRODUCT CODE: **UNEO**



### TOUCH STIM

#### TOUCH SCREEN TENS

An OTC dual-channel 24-mode device that has an easy to select settings – 12 TENS and 4 muscle release modes, and a large touch screen LCD screen.

PRODUCT CODE: **TSTE10**



#### IFC4000 DIGITAL

The IFC4000 is a digital sine-wave Interferential that runs on a Li-ion ultrathin battery for untethered treatments. The IFC digital is a dual channel full featured digital device that is compact and portable.

PRODUCT CODE: **IFC4000D**



#### MICRO EQC

A noninvasive means of pain relief with broad applications for acute pain. It can be used for the relief and management of post surgical and post traumatic acute pain problems. Its pulse amplitude & pulse frequency are adjustable and the timer can be set to 10, 20, or 40 minutes. specs are 0.5, 8.0 & 80.0 hz and will switch to 0.3, 8.0 & 80.0 hz.

PRODUCT CODE: **EQC**



## UltimaOTC

### THE ULTIMA OTC TENS

A dual-channel five-mode device that has an easy to select body-part diagram on a large LCD screen. This device is remarkable for its function, and traditional looks mixed with new-age digital style. Also includes a soft carrying case, and is approved for an over-the-counter no prescription necessary sale.

PRODUCT CODE:  
**UOTCW, UOTCB**



## QUAD STIM PLUS

### QUAD STIM-PLUS

The Quad-Stim is a combination Transcutaneous Electrical Nerve Stimulation (TENS) and Electrical Muscle Stimulation (EMS) digital device used for pain relief, relaxation, and stimulating muscles. The device has four controllable and independent output channels, each of which can be connected to an electrode pair. The device has adjustable pulse rate and width, 5 TENS modes, 3 EMS modes, 2 adjustable timers, 24 preset programs, and can store up to 60 sets of operation records. Now available as an Over-The-Counter device.

PRODUCT CODE: **QUADSTIM**



### THE IF 4000 INTERFERENTIAL

An analogue Interferential stimulation device which is used for symptomatic relief and management of post surgical and post traumatic acute and chronic pain. The IF-4000 has four frequency shifts including continuous, 1/1 abrupt, 8/8 abrupt, and 10/10 ramped options.

PRODUCT CODE: **IF4000**



### ADVANCED RUSSIAN STIM WITH INTERFERENTIAL

A dual-channel device for muscle exercise and strength training. Powered by an AC/DC adaptor or 9volt battery. The 2500Hz provides a deeper and more comfortable stimulation, packaged with interferential therapy for deep tissue treatments.

PRODUCT CODE: **R3500**



### DIGITAL HIGH VOLT PULSED GALVANIC STIMULATOR

Digital Galvanic Stimulation is most useful in acute injuries associated with major tissue trauma with bleeding or swelling. In contrast to TENS or IFC units, which apply alternating current, Galvanic Stimulators apply direct current. Direct current creates an electrical field over the treated area that, theoretically, changes blood flow and moves fluid. The positive pad behaves like ice, causing reduced circulation to the area under the pad and reduction in swelling. The negative pad behaves like heat, causing increased circulation, reportedly speeding healing.

PRODUCT CODE: **DIGITAL**





## Ultima Neuro

### NEUROPATHY STIMULATOR

The Ultima Neuro is a Neuropathy stimulator that utilizes enhanced TENS specifications coupled with a monophasic wave form to target pain associated with early onset peripheral neuropathy in the hands and feet.

PRODUCT CODE: **UNEURO-H, UNEURO-F**

Available conductive garments sold separately:



GLOVE



SOCK



CUFF



## JumpStim

### JUMP STIM

A simple and reliable OTC TENS stimulator that is perfect for travelling. The wired connection between active and ground allows for a "true" current loop. Offers 5 push-of-a-button programs indicated by an LED color. Includes long-lasting Li-ion rechargeable battery.

PRODUCT CODE: **JS100**

# softcycle

## PELVIC FLOOR STIMULATOR

Softcycle is an electrotherapy medical device kit used to train the muscles of pelvic floor. Pelvic floor muscle training can help strengthen the muscles under the uterus, bladder, and bowel (large intestine). They can help women who have problems with urine leakage or bowel control.

It is a convenient, affordable, reliable and effective solution that helps patients improve their pelvic health and gain greater control over incontinence.

PRODUCT CODE: SC100



## SOFTCYCLE PROBES

Softcycle probes are specially designed for treating Urinary Incontinence.

PRODUCT CODE: UPROBE-S, UPROBE-L



# OH!STIM

## ORGASMS BACKED BY CLINICAL RESEARCH

8 out of 9 participants reported more intense arousal, better vaginal lubrication, or increased ability to reestablish orgasm.

University of Michigan scientists uncovered an effective and non-invasive treatment for female sexual dysfunction and a method to enhance an Orgasm: Electrotherapy. Light and targeted electro stimulation can improve and enhance sexual response, including orgasm.

Oh! Stim is built on the foundation of statistically significant results shown during scientific research.

PRODUCT CODE: OHSTIM





### MEDICAL GRADE ULTRASOUND/TENS COMBO

Pain affects the quality and enjoyment of life, especially for those who suffer chronic pain. The Ultrasound Combo combines electrical therapy with ultrasound therapy in one device.

- 14 TENS programs with 1 MHz ultrasound
- Easy application LCD backlight indication
- Power Supply: medical adapter
- Professional grade

PRODUCT CODE: UTC1000



ULTRASOUND SPECS	
ACOUSTIC FREQUENCY	1 MHz ± 10%
GENERATOR OUTPUT	2.0W ± 10%
MODULATION WAVE SHAPE	20 Hz ± 10%
5 LEVEL OUTPUT INTENSITY WITH DUTY FACTOR	5%, 20%, 50%, 80%, 100%
TIMER	10 min
ACTUAL RADIATING AREA (AER)	4.0 cm <sup>2</sup> ± 20%
ACTUAL INTENSITY	0.5 W/cm <sup>2</sup> ± 20%
RBN	5.0 ± 30%
BEAM TYPE	collimated
MATERIAL OF APPLICATOR	stainless steel
TENS SPECS	
OUTPUT CHARACTERISTICS	constant voltage (CV)
READING RESOLUTION	1V
PROGRAM	each 7 programs for C and U Mode
SINGLE DIRECTIONAL PULSE	43 x 43 x 10 cm / 17" x 17" x 4"
FREQUENCY	2-150 Hz
PULSE DURATION	60-250uS
MAXIMUM OUTPUT	80Vpp (1 K ohm)



### PM 2000 PORTABLE ULTRASOUND MACHINE

This hand-held portable ultrasound device provides therapeutic heat. Ultrasound penetrates far beneath the skin to soothe and comfort chronic or acute pain.

PRODUCT CODE: PM2000


PM 2000	
FREQUENCY	1MHz +/-10%
PULSE WIDTH	2ms +/-10%
REPETITION RATE	150Hz +/-10%
WAVE SHAPE	Pulsed
TEMPORAL MAX. POWER	4W +/-10%
TEMPORAL MAX. EFFECTIVE INTENSITY	1.57W/cm <sup>2</sup> +/-10%
EFFECTIVE (EFF.) AVERAGE EFF. INTENSITY	0.47W/cm <sup>2</sup> +/-10%
EFF. MAX TEMPORAL INTENSITY/EFF. AVG. EFF INTENSITY= TEMPORAL PEAK AVG. INTENSITY RATIO	1.57/0.47=3.34
ERA (EFF. RADIATING AREA)	6.16 square centimeter +/-5%
BNR (BEAM NON-UNIFORMITY RATIO)	Max. 5.6
AUTO-TIME SETTING	30 minutes +/-1min.
OUTPUT INTENSITY	Low, Medium, High
SIZE OF MAIN UNIT	172(L)x54(w)x42(h)mm
WEIGHT OF MAIN UNIT	120g (w/o adaptor)



### WINSTIM- ELECTROTHERAPY + ULTRASOUND

Winstim 4 Channel Electrotherapy modalities in a single unit including combination of ultrasound and electrotherapy.

PRODUCT CODE: WINSTIMU



### CLINICAL COMBO IN, TENS, GALVANIC RUSSIAN

Analgesic Pulsar AP439, the cutting edge multi waveform, multiple channels device presents the maximum applications of pain relief.

- 4 Channel, 8 Electrodes
- Symmetrical Biphasic Square Wave
- MKB Panel with Digital Display
- Sine Wave, Symm Biphasic, DC
- Carrier Freq. 4000Hz and 2500Hz
- Pulse Rate: 1-150Hz
- Pulse Width: 50µs-250µs
- Output: 0 – 120V Peak @ 500 ohm
- Timer: 1-60 min

PRODUCT CODE: AP 439



### SOUND PRO

The SoundPro is our most economical therapeutic ultrasound series which can provide 1MHz as well as 3Mhz treatments. The SoundPro has a very versatile 5cm soundhead that works well for all body parts.

PRODUCT CODE: SPRO

# soft-touch

## SOFT-TOUCH ELECTRODES

A self adhesive reusable disposable electrode with a pig tail pin connector. Available with a carbon and a blue cloth backing. Oval corners help minimize edge curl. Our brand of electrodes delivers great dispersion, durability, and is value priced. Average use is 10-15 uses, with resealable pouch.



### SOFT-TOUCH SILVERS

The Soft-Touch Silvers use highly conductive silver backing with an ultra-low impedance and a clean blue gel.

SIZE	ITEM #
1.25" x 1.25" (4-pack)	SPS1250
2" x 2" (4-pack)	SPS2020
2" ROUND (4-pack)	SPS2000
2" x 4" (4-pack)	SPS2040
3" ROUND (4-pack)	SPS3000



### SOFT-TOUCH CLINICAL GRADE

The Soft-Touch Clinical grade uses a carbon backing with clean blue gel with a low impedance gel.

SIZE	ITEM #
1" ROUND (4-pack)	SP1000
1.5" x 1.5" (4-pack)	SP1515
1.5" x 3.1" (4-pack)	SP1531
1.5" x 4" (4-pack)	SP1540
2" ROUND (4-pack)	SP2000
2" SQUARE (4-pack)	SP2020
1.5" x 6" (2-pack)	SP1560
3" ROUND (4-pack)	SP3000
3.3" x 6" BUTTERFLY (1-pack)	SP3360



### SOFT-TOUCH MEDICAL GRADE

The Soft-Touch Medical grade uses a carbon backing with a standard grade gel.

SIZE	ITEM #
2" ROUND (4-pack)	FA2000
2" SQUARE (4-pack)	FA2020



### SOFT-TOUCH BASIC GRADE

The Soft-Touch Basic grade uses a carbon backing with a standard grade gel at a thinner durometer.

SIZE	ITEM #
2" ROUND (4-pack)	BG2000
2" SQUARE (4-pack)	BG2020



### STARBURST HYPOALLERGENIC ELECTRODES

The Starburst electrodes are very unique with a conductive pattern that helps radiate the current outward and dispersively. A very skin sensitive gel is used as well.

SIZE	ITEM #
2" ROUND (4-pack)	EP85340
2" SQUARE (4-pack)	EP85345
SUPERIOR SILVER WITH POLYHESIVE BLUE GEL ELECTRODES	
	EP85305

### SOFT TOUCH FOAMS WITH PMT GEL

SIZE	ITEM #
2" ROUND (4-pack)	FAF2000
2" SQUARE (4-pack)	FAF2020



TEN CARE + ACCESSORIES	
PRETENS GEL	M57170
PULSE WIDTH	MS71130
REPETITION RATE	M571130
WAVE SHAPE	UP220
PRETENS WIPES	PTW
LEAD WIRE (PAIR)	ET99090101
ELECTROLYTE SPRAY 4OZ.	CS4
ELECTROLYTE SPRAY 8OZ.	CS8
SPLITTER CABLE	SCABLE



# EB ENERGY BRACE

## ELECTROTHERAPY GARMENTS

Now electrotherapy can be easily applied to hard to reach spots on the body with incredible accuracy. This unique conductive mesh material can stimulate large or multiple areas of the body.

The EB design works in conjunction with a wide variety of Electro-Medical devices.



	CODE	DESCRIPTION	SIZES (INCHES)
	EKD47	ELECTRIC KNEE W/1 - 4X7 DUAL ELECTRODE	KNEE
	EGLOVE	ELECTROTHERAPY GLOVE - ONE SIZE FITS ALL	HAND
	ESOCK	ELECTROTHERAPY SOCK - ONE SIZE FITS ALL	FOOT
	EASLEEVE	ELECTROTHERAPY SLEEVE MEDIUM ARM*	ARM
	ELSLEEVE	ELECTROTHERAPY SLEEVE LARGE LEG	LEG
	ECUFF	ELECTROTHERAPY CUFF - ONE SIZE FITS ALL	N/A
	EV	ELECTRIC VEST W/1 - 4X10 DUAL ELECTRODE	VEST
	EBFOOT	ELECTROTHERAPY FOOT - ONE SIZE FITS ALL	N/A

	CODE	DESCRIPTION	SIZES (INCHES)
	EBQ	QUAD GARMENT ELECTRODE - LOW, MID OR UPPER BACK	BACK
	EB523	ELECTROTHERAPY SINGLE CONDUCTIVE PAD	2 x 3
	EB535	ELECTROTHERAPY SINGLE CONDUCTIVE PAD	3 x 5
	EB547	ELECTROTHERAPY SINGLE CONDUCTIVE PAD	4 x 7
	EBD35	ELECTROTHERAPY DUAL CONDUCTIVE PAD	3 x 5
	EBD47	ELECTROTHERAPY DUAL CONDUCTIVE PAD	4 x 7
	EBD410	ELECTROTHERAPY DUAL CONDUCTIVE PAD	4 x 10
	EBD35/S	ENERGY BRACE KIT - SHOULDER - DUAL CONDUCTIVE 3X5 PAD W/ WRAP & VERSITILE JOINT WRAP	3 x 5
	EBD35/C	ENERGY BRACE KIT - CERVICAL - DUAL CONDUCTIVE 3X5 PAD W/WRAP	3 x 5
	EBD47/K	ENERGY BRACE KIT - KNEE - DUAL CONDUCTIVE 4X7 PAD W/WRAP	4 x 7

	CODE	DESCRIPTION	SIZES (INCHES)
	EBD410/B	ENERGY BRACE KIT - BACK - DUAL CONDUCTIVE 4X10 PAD W/WRAP	4 x 10
	EBDB840B	ENERGY BRACE KIT - DELUXE BACK - WITH-2 DUAL CONDUCTIVE 4X10 PADS W/WRAP	8 x 10
	EB1-410	ELECTRIC BACK WITH 1 - 4X10 DUAL ELECTRODE	BACK
	EB2-410	ELECTRIC BACK WITH 2 - 4X10 DUAL ELECTRODE	BACK
	VJW	VERSATILE JOINT WRAP	N/A
	EBW440	ELASTIC BANDAGE 4X40	N/A
	EBW840	ELASTIC BANDAGE 8X40	N/A
	CS4	CONDUCTIVE SPRAY - REGULAR SIZE	4 OZ.
	CS8	CONDUCTIVE SPRAY - LARGE SIZE	8 OZ.



**X-FLEX KINESIOLOGY TAPE**

4-way stretch kinesiology tape means maximum flexibility. The durable, high-density fabric relieves pain and pressure from swelling or injured muscles and joints. The water-repellent but breathable tape stays on 5-7 days through showers, swimming, and workouts.

PRODUCT CODE: XFLEXBEIGE, XFLEXBLACK, XFLEXBLUE, XFLEXPINK



**JSTIM** 

**JOINT STIMULATION SLEEP THERAPY SYSTEM**

The JStim is a Proprietary break through in the treatment of Osteoarthritis or Rheumatoid Arthritis. It is FDA and Medicare cleared, and through combinational therapy it deploys the use of advanced electrotherapy, and true compression therapy.



**3-STRAP KNEE SYSTEM (ELECTRODE)**

The JStim Knee System utilizes a high-quality knee compression wrap, and specific conductive silver fabric patches that connect to the JStim stimulation unit.

PRODUCT CODE: J1000KN

**JSTIM HAND SYSTEM**

The JStim Hand System is similar to the knee system, but is worn like a glove; this particular treatment set-up consists of silver garment fabric in the shape of a glove along with a special compression glove to be fitted over the conductive components in order to block the outside current and to give effective treatments to the affected hand.

PRODUCT CODE: J1000H



**JSTIM DEVICE**  
PRODUCT CODE: JSTIM1000



**ABOUT THE JSTIM SYSTEM**

- One size fits most
- Conductive Spray and 36" Lead Wires included with all systems
- Patent pending
- Medicare coded
- SADMERC approved
- FDA cleared

All systems include a custom carrying case.



## brightwand

UVB PHOTOTHERAPY

### UVB PHOTOTHERAPY

The brightwand Compact UV Light Therapy Lamp is for the treatment of Psoriasis, vitiligo, and other skin conditions. Narrow band UVB light treatment is a targeted UVB light spectrum which minimizes the development of psoriasis, lesions or re-pigmentation of one's normal skin coloring.

PRODUCT CODE: BWUV100



## INFÄREX

### INFARED LIGHT WAND

This hand-held infrared light wand provides therapeutic massage and deep relief. Infrared heat penetrates far below the surface of the skin to soothe and comfort the nerves.



PRODUCT CODE: PMTLT



**Theralamp Hand Held**  
PRODUCT CODE: TLHH

# LASER HERO

## LASER THERAPY HELMET

The Laser Hero is an advanced, FDA-cleared system for battling hair loss and preventing receding, thinning hair in both men and women using clinically proven low-level laser therapy (LLLT) and red light diodes.

The hair regrowth helmet delivers non-invasive, painless, and safe, low-level power and light energy to the scalp to energize follicles and promote natural hair growth for improved thickness, fullness, and pigment.

The hands-free design is easy to use, small in size, and lightweight, making it ideal for home use. The Laser Hero is comfortable, secure, and effective for all hair types.



Available in two models:

### STANDARD

- 80 laser diodes and 30 LED red lights
- 7-9 months of recommended usage

PRODUCT CODE: **LTH100B**

### PREMIUM

- 204 laser diodes
- 6 months of recommended usage

PRODUCT CODE: **LTH200G**

Intended for seborrheic alopecia, alopecia areata and other hair loss treatment
Medical-grade semiconductor laser diodes (LD), and red light LEDs deliver uniform treatment
Dot matrix arrangement of laser light source to promote rapid hair growth
Partition irradiation technology and LLLT low energy soft laser irradiation technology
Intelligent sensor to ensure safety and effectivity during treatment
Real-time detection of temperature and humidity



# THERALAMP

## RELIEVE JOINT AND MUSCLE PAIN RELIEF INCREASE BLOOD VESSEL RESILIENCE AND PROMOTES CELLULAR FUNCTIONING

The Theralamp uses a special filter to ensure UV proof and high transmission of infrared heat to penetrate deep into the skin. The radiant warmth dilates blood vessels, speeds oxygen flow and increases blood circulation, helping the body to detoxify and transport necessary substances for rebuilding and nourishing body tissue.

- Specially designed ceramic glass cover to safely block all harmful UV
- 200 watt power lamp with 5000 hours of life
- Adjustable tilt of 35 degrees upward and 10 degrees downward for convenient use
- One button operation with 1-15 minute timer settings

PRODUCT CODE: TL100H



# THERALAMP RED LIGHT TORCH

## THERALAMP – RED LIGHT TORCH

Targeted pain relief through high output red light therapy. 3000mW, 660nm and 850nm wavelength

PRODUCT CODE: TLRLT10







## **POLAR** VORTEX

Provides controlled cold water therapy for extended periods, allowing this thermal energy to meet a variety of indications helping patients and practitioners manage pain or swelling.

- Reduces pain and swelling to help speed up rehabilitation
- Portable and lightweight
- Clear reservoir allows for easy to see ice-level

PRODUCT CODE: **PV105**



## **arctic** ice system

The Arctic Ice System delivers localized cold therapy to patients either in the home or during their hospital visit. With a stylish new design, the AIS provides four to five hours of continuous cold therapy; lightweight and fully portable.

PRODUCT CODE: **CTU2**



## **aqua** relief system

### **NEW HOSPITAL GRADE SYSTEM FOR USE IN THE HOME**

The Aqua Relief System delivers pain relief to achy feet and other body parts due to arthritic pain, carpal tunnel syndrome, back pain, and many other pain conditions. The ARS delivers localized cold or hot water fluid therapy to patients either in the home or during their hospital visit.

- Hospital-grade plug
- 60 min. timer, plus auto shut off
- Hot and Cold heat settings
- Anti-interference device

*\* All units for use with a universal cold therapy pad (sold separately) that fits many different areas of the body, including: knees, shoulders, ankles, elbows, calves and limbs.*

PRODUCT CODE: **ARS2000**



# ARCTICICE CLEAR

## ARCTIC ICE CLEAR


The Arctic Ice Clear delivers localized cold therapy to patients either in the home or during their hospital visit. Offers four thermal therapy settings with auto-shutoff. Simple to use and offers 4-5 hours of continuous cold therapy.

**Portable:** Designed for easy portability and storage

**Convenient:** Clear cooler design for easy refill indication

**Simple:** Easy to use by LCD & button. LCD shows the treating time, working time and temperature. Comfortable: adjustable sleeves are made of soft breathable fabric.

PRODUCT CODE: CTU2C



## COLD CUBES

Conveniently swap out your ice and enhance the duration of continued cold therapy with our all new cold cubes. Enhance the cold therapy duration from 6-8 hours, instead of 3-4 hours when using regular ice that needs emptied.

PRODUCT CODE: CCB2



## VERSATILE JOINT WRAP "VJW"

Conveniently wrap the universal water therapy bladder around an elbow, shoulder, hand, knee, or ankle/foot using the new and improved versatile Joint wrap.

PRODUCT CODE: VJW



	THERAPY PADS	ITEM #
	UNIVERSAL THERAPY PAD	TPU
	WATER THERAPY - STANDARD SINGLE BOOTIE	TPSNG-SM
	WATER THERAPY - STANDARD DOUBLE BOOTIE	TPDBL-SM
	WATER THERAPY - LARGE SINGLE BOOTIE	TPSNG-LG
	WATER THERAPY - LARGE DOUBLE BOOTIE	TPDBL-LG
	FACE COLD PAD	FC
	HEAD COLD PAD	IC

	THERAPY PADS	ITEM #
	LARGE BACK PAD - 13" x 24"	TP1324
	THERAPY MITT	TPMITT
	VERSATILE JOINT WRAP	VJW
	CERVICAL THERAPY PAD	TCP
	SHOULDER COLD PAD	SCP
	HAND COLD PAD	HCP
	ANKLE COLD PAD	ACP
	WATER THERAPY KNEE PAD	TPKNEE

## ICY WRAP

### CRYO-COMPRESSION THERAPY

The Icy Wrap delivers simultaneous cold and compression therapy to control pain and swelling caused by sprains, pulled muscles, injured ligaments, or strained tendons.

PRODUCT CODE: LR100

DESCRIPTION	ITEM #
KNEE WRAP	IWKNEE
ANKLE WRAP	IWANKLE
SHOULDER WRAP	IWSHOULDER
BACK WRAP	IWBACK
ELBOW WRAP	IWELBOW
HAND WRAP	IWHAND



**RAPID KNEE  
FRONT WRAP-ON KNEEWRAP  
COMFORT FIT ELASTIC**

The Rapid Knee Double upright adjustable hinge knee brace allows for effectively setting flexion or extension settings needed for the given treatment needs. This Rapid Knee is a wrap on the front type of brace with a very comfortable high tech flat knitting.

PRODUCT CODE: **RK150**

SIZE	CIRCUMFERENCE OF THE KNEE
SMALL	30CM - 34CM
MEDIUM	34CM - 38CM
LARGE	38CM - 42CM
X-LARGE	42CM - 46CM
XX-LARGE	46CM - 50CM
XXX-LARGE	50CM - 55CM



**RAPID KNEE  
FRONT WRAP-ON KNEEWRAP  
COMFORT FIT ELASTIC**

The Rapid Knee Double upright adjustable hinge knee brace allows for effectively setting flexion or extension settings needed for the given treatment needs. This Rapid Knee is a slip-on type of brace made of high tech comfort fit elastic.

PRODUCT CODE: **RK150PLUS**



**RAPID KNEE  
RIGID WRAP-ON KNEE BRACE**  
The Rapid Knee Double upright adjustable hinge rigid knee brace allows for effectively setting flexion or extension settings needed for the given treatment needs. This Rapid Knee brace is a slip or wrap on rigid type of brace, also has quick release snaps.

PRODUCT CODE: **RK200**



**PMT SOFT KNEE WRAP**  
Soft lightweight neoprene knee wrap provides compression and can accommodate the EB Wraps.

PRODUCT CODE: **SOFTKNE**



**RAPID KNEE OA BRACE**  
The Rapid Knee single upright adjustable OA unloader knee brace allows for adjustable flexion or extension settings needed for the given treatment. Adjustable straps allow for a one-size-fits-all comfortable and snug fit.

LEFT OR RIGHT  
**ONE-SIZE-FITS-ALL**

PRODUCT CODE: **OARKR, OARKL**



Rapid Wrist



**RAPID WRIST**

The Rapid Wrist is a dual function wrist splint and finger traction tool. Treat wrist pain, stiffness, and discomfort and build strength with digital exercise through adjustable traction.

PRODUCT CODE: RW200



**TR RapidThumb**

**RAPID THUMB**

The Rapid Thumb is a rigid thumb brace for pain relief through immobilization of the carpometacarpal (CMC) and metacarpophalangeal (MCP) joints.

PRODUCT CODE: RT100

**THE stringback**

NEW ONE-SIZE-FITS-ALL



**THE STRING BACK LSO SPINAL BRACING SYSTEM**

The String Back LSO utilizes a fixed wheel and string system to cinch the user into a supine position and immobilize the back region. The String Back system provides circumferential compression, and is essential in stabilizing the spine, and promoting the healing process following post-operative back surgery.

- Through our String Back technology, the brace is cinched to the spine to create an anti-gravity lift
- Breathable mesh allows for active and long term use
- Provides an anatomically correct fit and adjustable lumbar compression

**INDICATED FOR USE FOR:**

- Chronic back pain
- Stenosis
- Spondylosis
- Degenerative disc disease
- And more!

PRODUCT CODE: SB150, SB100



**RA RapidAnkle**

**RAPID ANKLE**

Provides ankle support and stability with an ambidextrous design (fits left and right foot). Indicated for: Tibialis tendonitis, ankle instability, sprains/strains, osteoarthritis, ankle fracture, and non-traumatic ruptures of tendon.

- Adjustable velcro
- Calf and lateral support
- Full flexion hinge
- Moldable foot plate

PRODUCT CODE:

RA110SM - Size small/medium

RA110LGXL - Size Large/XL

# THERMOTECH.

## MEDICAL GRADE INFRARED AUTOMATIC MOIST HEAT PAD

Provides moist heat treatments by a process of fomentation without adding water for the relief of all types of pain conditions. During the healing process, moist heat therapy has proven most effective to ease pain, reduce stiffness and increase circulation. Includes flannel cover.



ITEM #	TECHNOLOGY	SIZES (INCHES)
S766	ANALOGUE	KING – 26" x 14"
S767	ANALOGUE	MEDIUM – 18" x 14"
S768	ANALOGUE	MINI – 19" x 7"

ITEM #	TECHNOLOGY	SIZES (INCHES)
S766D	DIGITAL	KING – 26" x 14"
S767D	DIGITAL	MEDIUM – 18" x 14"
S768D	DIGITAL	MINI – 19" x 7"
S769D	DIGITAL	SHOULDER – 20" x 20"



## THERMOTECH CONVENTIONAL MOIST HEATING PAD STANDARD

An economical choice for creating moist or dry heat. Just add water to the sponge layer to create moist heating.

ITEM #	TECHNOLOGY	SIZES (INCHES)
S708M	ANALOGUE	KING – 24" x 11.5"
S707M	ANALOGUE	MED. – 14" x 11.5"

# THERMORELIEF



Available with red flannel cover

## THERMORELIEF MEDICAL GRADE AUTOMATIC MOIST DIGITAL HEAT PAD

Provides moist heat treatments by a process of fomentation without adding water for the relief of all types of pain conditions. During the healing process, moist heat therapy has proven most effective to ease pain, reduce stiffness and increase circulation. Includes flannel cover.

PRODUCT CODE: S767-DC (Medium), S766-DC (King size)

## THERMORELIEF BASIC MOIST OR DRY HEATING PAD

Helps relieve pains of inflamed joints caused by arthritis and rheumatism. High-heat automatic moist heat therapy increases blood flow to sore muscles. Ideal for treatment of back, abdomen, shoulder, chest, and more.



PRODUCT CODE: TTE100



### HIGH PERFORMANCE INFRARED HEAT THERAPY

This far infrared heat therapy system provides soothing relief from pain, muscle tension and stress. Far infrared is the safe form of light energy from the sun that we depend upon to remain healthy.

- Speed recovery
- Relax muscles
- Relieve pain
- Fits up to size 42" waist
- Improve circulation

PRODUCT CODE:  
QBW, QWW



### ARCTICXL

#### EXTRA-LARGE COLD PACK

Extra-large insulated cold pack that is ultra-flexible with superior coverage for the high performance athlete. Includes a wrap strap for affixing it to the treatment site.

Size does matter.

Available in two sizes: 10" x 6" and 18" x 12"

PRODUCT CODE: AXL106, AXL1812



### MECHANICAL LIFTING SEAT

Maintain an active independent lifestyle with help from Seat Boost. Seat Boost is the portable mechanical lifting seat that helps people stand up or sit down without assistance.

- Safe, durable and reliable
- Easy care with washable covers and cushions that can be cleaned with mild soap and water
- Adjustable lift for those from 94lbs to 350lbs
- No batteries or electricity needed. Seat Boost is self-powered; take it anywhere

PRODUCT CODE: SB1000S, SB1000L



### PRESSURE RELIEF CUSHION SYSTEM

This advanced pressure relief cushion system is especially designed for the patients to use on various mobile devices. The MobiCushion is fully portable, with capabilities for both battery power and wall outlets. With the advanced pressure sensor technology and unique one piece air cell pad design, it provides alternating low pressure to the patients consistently.

- Suitable for most wheelchairs. Portable pump with carry bag and cushion design for easy mobility.
- Adjustable patient comfort and alternating cycle time. Features visual and audible operating status and fault indicators.
- Unique seat cushion design.
- Stretched, vapor permeable, water resistant and fire retardant coverlet.
- Maximum patient weight capacity up to 120 kg or 264 lb. Therapy Mode: Static / Alternate.

PRODUCT CODE: MC1000





## STRETCH

### TOE TRACTION SOCKS

Give your toes a traction like stretch that not only relaxes but realigns for proper foot health. Helps with bunions, hammer toes, and plantar fasciitis.

PRODUCT CODE: TTS



## Therchands

### COMPRESSION ARTHRITIC GLOVES

Keep your hands comfortable while you engage in daily activities. They provide mild compression to aid in improving circulation and relieving arthritis symptoms.

See online for sizing details.

PRODUCT CODE: TH100B, TH100P



### DR. ARCHY FOOT MASSAGER

The design of Dr.Archy takes advantage of the muscle structure of the foot and lower leg to provide a natural low-impact exercising motion that supports, stretches, and massages your arches.

Dr.Archy can provide temporary relief from pain associated with plantar fasciitis, spurs and tired feet.

The all-in-one foot tool can be used to therapeutically treat your arches at your desk, while traveling, and just about anywhere you need.

PRODUCT CODE: DA1



### DR. NECKY NECK MASSAGER

Dr. Necky is the wonderful little tool that will work your neck into relief like a professional masseuse. The special rolling wheels have a knuckle – like pattern that works like no other tool on the market, to loosen up those tight muscle knots.

Dr. Necky can be used as a light massage or a deep tissue rub. Best of all, it is small, portable, and convenient to use anywhere!

PRODUCT CODE: DN1



### BODY DRUMMER HOT/COLD VIBRATING MASSAGER

The Body Drummer hot/cold massager is a unique massager that not only provides a soothing vibrating massage but also allows for either hot or cold treatment simultaneously. The body drummer is handheld and rechargeable.

PRODUCT CODE: BD100





# BODY DRUMMER PRO

## BODY DRUMMER PRO HIGH RPM VIBRATION MASSAGER

The Body Drummer Pro is an intense, painless vibrational massager that sends vibrations deep into the muscles for pain relief.


Powerful deep tissue massager enhances performance, promotes circulation, increases range of motion, and relieves pain.

A lightweight design and built with a powerful motor which reaches up to 3900 RPM, includes 4 accessories, while battery lasts up to 2 hours on a full charge.

Light, powerful, and quiet.

PRODUCT CODE: **BDPRO100**

BODY DRUMMER PRO	
BATTERY CAPACITY	7.4V 2600 mAh
POWER BY	DC Adapter Input 100-240V/50-60Hz
OUTPUT VOLTAGE	9V DC 1A
MAX CURRENT	4000 mA
POWER CONSUMPTION	30 W
1ST MODE SPEED	2800 RPM (46 Hz)
2ND MODE SPEED	3400 RPM (56 Hz)
3RD MODE SPEED	3900 RPM (65 Hz)
CHARGING TIME	3-4 Hours
WORK FOR WHEN FULL CHARGE	2 Hours
WEIGHT (ONLY MASSAGER)	682 g (1 lb)
WEIGHT (FULL SET)	960 g (2.12 lbs)
OVERALL SIZE OF MASSAGER	23.5 x 6 x 19 cm



# BODY DRUMMER PRO X

## BODY DRUMMER PRO X PERCUSSION MASSAGER

The Body Drummer Pro X is an intense, painless percussion massager that sends percussion waves deep into the muscles for pain relief.

Powerful deep tissue massager enhances performance, promotes circulation, increases range of motion, and relieves pain.

A heavy duty design and built with a powerful motor. Includes 5 accessories and the battery lasts up to 3 hours on a full charge.

Heavy duty, powerful, and with a brushless motor.

PRODUCT CODE: **BDPRO-X**

BODY DRUMMER PRO	
POWER BY	DC Adapter Input 100-240V/50-60Hz
OUTPUT VOLTAGE	24V DC 1A
MAX CURRENT	4000 mA
POWER CONSUMPTION	68 W
1ST MODE SPEED	2700 RPM (46 Hz)
2ND MODE SPEED	2900 RPM (56 Hz)
3RD MODE SPEED	3100 RPM (65 Hz)
CHARGING TIME	3-4 Hours
WORK FOR WHEN FULL CHARGE	3 Hours
WEIGHT (ONLY MASSAGER)	682 g (1 lb)
WEIGHT (FULL SET)	960 g (2.12 lbs)
OVERALL SIZE OF MASSAGER	23.5 x 6 x 19 cm





### PAIN GRENADE

Military Grade pain relief through hot and cool therapy. The low-odor formula is nourishing and moisturizing to the skin while warming and cooling relief activates without feeling sticky or greasy. Available in 3 formulations – Gel, Spray, and Roll-On.

SIZE	ITEM #
4 oz. GEL	PGGEL
4 oz. SPRAY	PGSPRAY
3.4 oz. ROLL-ON	PGROLL

**ACTIVE INGREDIENTS:**  
Menthol 3.0%  
Camphor 3.2%

**INACTIVE INGREDIENTS:**  
Eucalyptus Globulus Leaf Oil  
Arnica Montana Flower  
Hamamelis Virginiana (Witch Hazel)  
Camellia Sinensis Leaf  
Aloe Barbadensis Leaf  
Glucosamine  
Petersonii Oil  
Magnesium



## FLIP STRETCH

### FLIP STRETCH

Inversion therapy has been used for centuries for stretching, strength building, and pain relief. An inversion table can help you counter the effects of gravity on your body by placing you at a downward angle for just a few minutes. Inversion tables stretch and release tension around the spine, helping to reduce stress, improve blood circulation, and increase flexibility.

With an inversion table, you can secure yourself in the downward position and adjust the angle of the table for ultimate benefits and stress relief.

For your convenience, support, and comfort, the backrest is hard and allows for accessories to be added.

PRODUCT CODE: **FS100**



## THE BACK BAT

### SPINAL TRACTION LUMBAR SUPPORT PILLOW

Discover what complete back support feels like with the vertical spinal traction cushion that relieves pressure on your neck and back. Originally part of an inversion table for reducing lower back pain, this cushion can be used seated or lying down.

Use it in your favorite office chair, or in a lying down position to elongate the spine and improve posture.

The lightweight, flexible design can be easily moved from place to place, making it ideal for travel or in your favorite chair at home. With durable grippers and an elastic strap, it can be fastened into any seat.

PRODUCT CODE: **BBT100**



## neckCloud

### NECK CLOUD

Neck Cloud is a portable, easy to use cervical traction device created to help relieve neck pain. Relax and feel like you're floating away on a cloud of relief.

PRODUCT CODE: **NC100**



## Acupillow

### NECK STRETCH AND MASSAGE

Brings relief through acupressure massage and cervical traction. Firm but flexible design gives convenient comfort and relaxation.

PRODUCT CODE: **AP100R**





## DYNAMIC WEDGE

AUTOMATIC LUMBAR TRACTION DEVICE

### AUTOMATIC LUMBAR TRACTION DEVICE

The Dynamic Wedge Automatic Lumbar Traction Device uses heat, light, vibration, and pneumatic air to create traction in the lower lumbar region and relieve back pain.

The device is specially designed to incorporate a variety of treatments in one. The dynamic functions give the choice of a pre-programmed automatic setting, which runs through levels of traction and vibration, or manual settings, which control the intensity of traction and the option of vibration massage at 3,000 RPM. The temperature is easily adjustable, while infrared light penetrates beneath the skin, deep into the tissue where healing occurs.

PRODUCT CODE: **LR100**



## DYNAMIC WEDGE

CERVICAL

### AUTOMATIC NECK TRACTION DEVICE

Heat therapy, electrotherapy, and multifunction air traction in an automatic device for relief from neck pain.

PRODUCT CODE: **DWC100**



## Air Collar

### NECK TRACTION DEVICE COLLAR

This neck traction device collar gently stretches neck muscles allowing the vertebral discs to realign, freeing the nerve root tissue from pressure.

PRODUCT CODE: **NTD3000**



CVT1000

CVT2000

### NECK SUPPORT COLLAR

- Inflatable cervical traction unit used to treat vertebrae symptoms, disease and discomfort
- Portable and convenient for use when traveling
- Air filled neck stretching device, ideal for home or office use



**PNEUMATIC CERVICAL TRACTION UNIT  
INDICATED FOR THE USE, TREATMENT AND RELIEF  
OF CERVICAL MUSCULOSKELETAL OR NEUROLOGICAL IMPAIRMENTS**

The Theratrak gently stretches neck muscles allowing the vertebral discs to realign, freeing the nerve root tissue from the pressure of the discs. The Theratrak helps “locked” neck muscles to slowly stretch and relax, allowing misaligned vertebrae to resume their normal supportive position.

The Theratrak cradles the head and neck via a curvature in the neck portion of the device, as well as a special head rest section that is vertically pneumatic. Both the headrest and neck region, inflate with over 30 lbs of adjustable pneumatic pressure to create traction in the cervical curve.

- Clinical grade, designed for the home
- Innovative and effective design

PRODUCT CODE: PCT100S, PCT100R, PCT100W



**THERATRAC – GLIDE**

Indicated to help treat and relieve herniated discs, radiculopathy, osteoarthritis, and pain or headaches due to soft tissue or joint stiffness. The Theratrak Glide is a clinical grade cervical traction device, utilizing an easy to use dial and cable ratchet system that applies traction to while in the supine position.

PRODUCT CODE: TTG10





**THERATRAC – LSO**  
**INDICATED FOR THE USE, TREATMENT AND RELIEF OF CERVICAL MUSCULOSKELETAL OR NEUROLOGICAL IMPAIRMENTS**

The Theratrak LSO (model TT100) utilizes a decompression system to cinch the user into a supine position and immobilize the back region. The Theratrak LSO system provides circumferential compression and decompression, and is essential in stabilizing the spine. Incremental stretching of the spinal column helps stimulate the healing process following post-operative back surgery.

PRODUCT CODE: **TT100**



**THERATRAC – LUMBAR TRACTION BRACE**


A longitudinal force to the axis of the spinal column is the definition of spinal traction. In other words, pulling parts of the spinal column in opposing directions to change or stabilize damaged spine positions.

The Theratrak is an easy to use, highly functional, and very effective lumbar traction device, giving the user the ability to treat their condition at the most comfortable clinic on earth...their home.

Theratrak utilizes pneumatic bladder technology giving it the unique ability to administer a moderate lumbar decompression stretch, to gently contour the innate spinal path.

PRODUCT CODE: **LT100M, LT100L**

SIZE	WAIST MEASUREMENT
SM /MEDIUM	29" – 37"
LARGE / XL	37" – 45"



**INCLUDES BONUS ARCH**

UltimateArch

**SPINAL TRACTION DEVICE**

Ultimate Arch Spinal Traction Device is a lower back tool to help back or sciatica pain.

Chronic back pain can be hard to address. The unique design of the Ultimate Arch gives the relief that you have been looking for, in the comfort and convenience of your own home. It's travel friendly, and easy to set up, use, and store. With 3 available settings, the Ultimate arch gives you the stretch you need to release tight muscles and low backpain, while improving posture. Includes: 3-level base, flexible support arch with foam comfort pad, chair attachment strap.

PRODUCT CODE: **UA100**



AIRE

**AIR COMPRESSION NECK MASSAGER**

Experience the ultimate neck massage with AIRE; it combines air pressure technology and vibration massage with soothing heat.

PRODUCT CODE: **ACN100**





**CONDUCTIVE COPPER METAL EAR CLIP ELECTRODES**

Conductive copper metal ear clip electrodes for stimulating acupuncture points on the ear lobe or for CES application. Plastic base handles of the clip accepts 2mm male pin from your lead wire.

Apply self-adhesive felt tip pads to the inside surface of the ear clips, then wet the felt tips with water or conductive gel. Has just enough tension for it to clip to the base of the ear lobe with a comfortable fit. 25 felt tip pads included.

MODEL	ITEM #
Conductive copper metal ear clip electrodes, approx. 1.75" long	EARCLP-M



**DUAL POINT PROBE ELECTRODE**

Stimulate acupuncture or pain points with this hand-held remote dual-point stimulation probe. Base of the handle accepts the positive and negative lead wire pins from your stimulator. Carbon rubber tips at the end of flexible stems offer safe and easy point stimulation. Works best with conductive gel or wet felt-tip pads applied to the carbon rubber tips.


SIZE	ITEM #
ENERGIZER 9 VOLT	EN22
ENERGIZER AA	EN91
9 VOLT RECHARGEABLE	9VR
9 VOLT CHARGER	9VCRG
AA RECHARGEABLE	AAR
AA CHARGER	AACRG

MODEL	ITEM #
CARBON RUBBER-TIPPED DUAL POINT PROBE ELECTRODE, APPROX. 6.75" LONG	DUOPOINT



**TIGER LEAD WIRES**  
Ultra-strong and flexible lead wires for TENS and EMS machines.

PRODUCT CODE: TLW



**2" CARBON ELECTRODES**  
Reusable 2" round carbon rubber pad.

PRODUCT CODE: CR20



**URO-PROBE AND DUO-PROBE UNIVERSAL VAGINAL SENSORS**

MODEL	ITEM #
URO-PROBE VAGINAL SENSOR LARGE	UProbe-L
URO-PROBE VAGINAL SENSOR SMALL	UProbe-S
DUO-POINT PROBE VAGINAL SENSOR LARGE	DUOPOINT



PAIN MANAGEMENT  
TECHNOLOGIES 

The logo for Pain Management Technologies, featuring the company name in a clean, sans-serif font and a stylized icon of three overlapping, upward-pointing rectangular shapes to the right.





PAINTECHNOLOGY.COM  
INFO@PAINTECHNOLOGY.COM  
1.800.239.7880