

PATHWAY® Aluminum Wheelchair Ramp Kit Technical Specifications

	RAMPS		
	4' RAMP	6' RAMP	
MADE OF ALUMINUM	✓	✓	
SLIP-RESISTANT SURFACE	✓	✓	
USABLE DIMENSIONS (L x W)	4' x 36"	6' x 36"	
OVERALL DIMENSIONS (L x W)	4′ x 40.1875″	6' x 40.1875"	
OVERALL DIMENSIONS – INCLUDING BASE FEET (L x W)	4′ x 41.1875″	6′ x 41.1875″	
SIDE RAIL CURB HEIGHT	2"	2"	
UPPER HANDRAIL HEIGHT	34.5"	34.5"	
LOWER HANDRAIL HEIGHT	24"	24"	
HANDRAIL DIAMETER	1.5"	1.5"	
WEIGHT W/O HANDRAILS	29.2 lbs.	42.3 lbs.	
WEIGHT W/HANDRAILS	45.2 lbs.	61.8 lbs.	
WEIGHT CAPACITY (DO NOT EXCEED)	1,000 lbs.	1,000 lbs.	
WARRANTY	Lifetime Limited	Lifetime Limited	
COUNTRY OF ORIGIN	USA	USA	

PLATFORMS			
4' x 4' PLATFORM W/ 4' HANDRAILS	5' x 5' PLATFORM W/ 5' HANDRAILS		
✓	✓		
✓	✓		
48" x 48"	60" x 60"		
53.5" x 53.5"	65.5" x 65.5"		
53.75" x 53.75"	65.75" x 65.75"		
Platform handrails include 1.5" x 2" curb	Platform handrails include 1.5" x 2" curb		
34.5"	34.5"		
24"	24"		
1.5"	1.5"		
57.2 lbs.	71.9 lbs.		
73.3 lbs.	95.5 lbs.		
1,000 lbs.	1,000 lbs.		
Lifetime Limited	Lifetime Limited		
USA	USA		

	GROUND TRANSITION	UPPER TRANSITION	BRIDGE PLATE
MADE OF ALUMINUM	✓	✓	✓
SLIP-RESISTANT SURFACE	✓	✓	✓
USABLE DIMENSIONS (L x W)	18" x 35.375"	3.5" x 35.375"	8" x 35.375"
OVERALL DIMENSIONS (L x W)	18" x 35.375"	3.5" x 35.375"	8" x 35.375"

	LANDING PAD W/O LOWER TRANSITION	LOWER TRANSITION	LANDING PAD W/LOWER TRANSITION STRAIGHT CONFIGURATION	LANDING PAD W/LOWER TRANSITION TURN CONFIGURATION
MADE OF ALUMINUM	✓	✓	✓	✓
SLIP-RESISTANT SURFACE	✓	✓	✓	✓
USABLE DIMENSIONS (L x W)	53.25" x 54.5"	6" x 35.375"	59.48" x 54.375"	53.25" x 60.17"
OVERALL DIMENSIONS (L x W)	53.64" x 54.5"	6" x 35.375"	59.48" x 54.375"	53.64" x 60.17"

	BASE FOOT
MADE OF POLYMER	✓
OVERALL DIMENSIONS (L x W)	7.875" x 7.875"

- Weights and dimensions are approximate and subject to change.
- If using with mobility equipment, refer to your mobility equipment's owner guide for acceptable usage, including proper slope/incline, chair direction, etc.; never exceed its recommendations.