

Sklar Disinfectant

BACTERICIDAL • VIRUCIDAL • FUNGICIDAL • TUBERCULOCIDAL
2 MINUTE KILL TIME

Sklar Disinfectant is a ready to use, highly effective cleaner and disinfectant for use in intensive care units and high-risk patient care areas, surgical suites, infant care units, plasmapheresis facilities and other areas where strict environmental control of cross contamination is critical.

Features	Benefits
Ready to Use	May be applied directly to surfaces or used as an immersion solution. See Instructions for Use for more details.
Broad Spectrum, Non-Caustic Formula	Ideal for Medical, Dental, Ophthalmic and surgical equipment and surfaces.
2 Minute Kill Time	Completely kills infectious disease-causing microorganisms including H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA and A. Baumannii. See Efficacy Report for a complete listing.
Safe for Various Surfaces	Ideal for use on surgical equipment and surfaces. Safely use on inanimate surfaces including plastics, stainless steel, chrome, glass, Plexiglas® and baked-on painted surfaces, without leaving residue.
Multi-Application Use	Highly effective as a surface disinfectant and cleaner, instrument decontaminant solution and ultrasonic disinfection cleaning solution. Always follow label instructions.
OSHA Compliant	Meets OSHA blood-borne pathogen requirements for disinfection in healthcare facilities and is formulated following strict FDA Good Manufacturing Practices.

Portfolio

SOLD SEPARATELY

Product Number	Features	Quantities
10-1643	24oz. Spray Bottles	Case of 6
10-1642	24oz. Spray Bottle	1 Bottle
10-1653	1 Gallon Bottles	Case of 4
10-1654	1 Gallon Bottle	1 Bottle
10-1644	2.5 Gallon Bottles	Case of 2

