

QUICK  
REFERENCE  
GUIDE

# SOLUTIONS INSPIRED BY YOU

Shiley™ Neonatal and Pediatric  
Tracheostomy Products



**Medtronic**  
Further. Together

## FEATURES AND SPECIFICATIONS



All Shiley™ neonatal and pediatric tracheostomy tubes have a recessed area behind the connector to provide enhanced comfort.

1. Tubes beginning with the letter N (NEF, NCF) feature a lower flange angle for neonatal patients.
2. Tubes beginning with the letter P (PEF, PCF, PELF, and PLCF) feature a higher flange angle for pediatric patients.
3. Products with the letter C feature the TaperGuard™ low-profile, low-pressure cuff.
4. Products with the letter L offer a longer cannula.

### Shiley™ Custom Tracheostomy Tube Service Online Ordering Tool

When a standard tracheostomy tube doesn't offer the right fit, turn to the Shiley™ custom tracheostomy tube service. Modified lengths, additional sizes, customized cuffs, and modified flange angles are just some of the available options. Use our online ordering tool at [ShileyCustomTrach.com](https://ShileyCustomTrach.com) or call 800-635-5267 for assistance.

### Shiley™ Tracheostomy Home Care Education

Visit [LivingwithaVent.com](https://LivingwithaVent.com) for tracheostomy home care videos and guides.

## Shiley™ Neonatal Tracheostomy Tubes

PRODUCT SKU	DESCRIPTION	I.D.	O.D.	LENGTH
2.5NEF	Neonatal tracheostomy tube cuffless	2.5 MM	4.2 MM	28 MM
3.0NEF		3.0 MM	4.8 MM	30 MM
3.5NEF		3.5 MM	5.4 MM	32 MM
4.0NEF		4.0 MM	6.0 MM	34 MM
4.5NEF		4.5 MM	6.7 MM	36 MM
2.5NCF	Neonatal tracheostomy tube with Taperguard™ cuff	2.5 MM	4.2 MM	28 MM
3.0NCF		3.0 MM	4.8 MM	30 MM
3.5NCF		3.5 MM	5.4 MM	32 MM
4.0NCF		4.0 MM	6.0 MM	34 MM
4.5NCF		4.5 MM	6.7 MM	36 MM

## Shiley™ Pediatric Tracheostomy Tubes

PRODUCT SKU	DESCRIPTION	I.D.	O.D.	LENGTH
2.5PEF	Pediatric tracheostomy tube cuffless	2.5 MM	4.2 MM	38 MM
3.0PEF		3.0 MM	4.8 MM	39 MM
3.5PEF		3.5 MM	5.4 MM	40 MM
4.0PEF		4.0 MM	6.0 MM	41 MM
4.5PEF		4.5 MM	6.7 MM	42 MM
5.0PEF		5.0 MM	7.3 MM	44 MM
5.5PEF		5.5 MM	7.9 MM	46 MM
5.0PELF	Extra long pediatric tracheostomy tube cuffless	5.0 MM	7.3 MM	50 MM
5.5PELF		5.5 MM	7.9 MM	52 MM
6.0PELF		6.0 MM	8.5 MM	54 MM
6.5PELF		6.5 MM	9.0 MM	56 MM
2.5PCF	Pediatric tracheostomy tube with Taperguard™ cuff	2.5 MM	4.2 MM	38 MM
3.0PCF		3.0 MM	4.8 MM	39 MM
3.5PCF		3.5 MM	5.4 MM	40 MM
4.0PCF		4.0 MM	6.0 MM	41 MM
4.5PCF		4.5 MM	6.7 MM	42 MM
5.0PCF		5.0 MM	7.3 MM	44 MM
5.5PCF		5.5 MM	7.9 MM	46 MM
5.0PLCF	Extra long pediatric tracheostomy tube with Taperguard™ cuff	5.0 MM	7.3 MM	50 MM
5.5PLCF		5.5 MM	7.9 MM	52 MM
6.0PLCF		6.0 MM	8.5 MM	54 MM
6.5PLCF		6.5 MM	9.0 MM	56 MM

© 2019, 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 01/19-12-AW-0045b(3)-[WF#2931265]

6135 Gunbarrel Avenue  
Boulder, CO 80301

800.635.5267

**medtronic.com**

**Medtronic**