

## MATERIAL SAFETY DATA SHEET

### Section 1: Product and Company Information

Product Name: Ostomy Appliance Deodorant  
Product Code: SNS40208, SNS40202  
Company: Safe n Simple LLC  
Phone number: Toll Free: 844-767-6334  
Website: [www.sns-medical.com](http://www.sns-medical.com)

### Section 2: Composition and Ingredient Information

Chemical Name: Proprietary  
Synonyms: None  
Cas No: N/A  
Chemical Formula: Proprietary  
Chemical Family: Liquid  
Molecular Weight: N/A

Chemical Name(s)	%Wt	Hazardous
	No	No

### Section 3: Physical Data

Appearance: Blue aqueous dispersion  
Odor: Mild  
Solubility: Soluble in water, insoluble in oils, alcohol and usual solvents  
Specific Gravity: 0.9900-1.200  
Ph @ 25C: 5.50-6.50  
Vapor Density (Air=1): N/A  
Decomposition: N/A  
Boiling Point: N/A  
Melting Point: N/A

### Sections 4: Fire and Exposition Hazard Data

Flash Point: N/A  
Flammable Limits: N/A  
Extinguishing Media: Water, Foam and carbonic anhydride  
Special Fire Fighting Procedure: If vapors or mist are generated, use a self-contained Respirator  
Unusual Fire & Explosion Hazards: N/A

### Section 5: Health Hazard Data and Safe Handling Information

Toxicity: Harmful by ingestion, irritating to eyes  
Threshold Limit Value: N/A  
Effect of Overexposure: N/A

Eye contact: Avoid contact with eyes  
Skin contact: N/A  
Inhalation: N/A  
Ingestion: Harmful by ingestion  
Mutagenicity: N/A

### **Section 6: Emergency and First Aid Procedures**

Eye contact: Wash with plenty of water  
If irritation persists, seek medical advice  
Skin Contact: Wash with plenty of water  
If irritation persists, seek medical advice  
Inhalation: Remove to fresh air  
Seek medical advice  
Ingestion: Drink water and induce vomiting  
Seek medical advice  
Note to Physician: N/A

### **Section 7: Reactivity Data**

Stability: Stable  
Incompatibilities: None  
Hazardous Decomposition Products: Combustion may produce carbon dioxide and reducing agents  
Harzardous Polymerization: None  
Conditions to Avoid: Incompatible with strong oxidizing and reducing agents

### **Sections 8: Spill or Leak Procedure**

In case material is released or spilled: Dispose of through an authorized waste disposal Company  
Waste Disposal Method: Prevent environmental contamination  
Dispose of quantities greater than 10 gallons through an authorized waste disposal company

### **Section 9: Special Protection Information**

Ventilation Requirements: None  
Protective Equipment: None

### **Section 10 - Disposal Considerations**

Waste Disposal Methods: DISPOSE OF IN ACCORDANCE W/LOCAL, STATE & FEDERAL REGULATIONS

### **Section 11 – Transport Information**

DOT Shipping Name: Not Restricted in accordance with CFR 49 Section 173.125(1)

DOT UN Number: Not Restricted in accordance with CFR 49 Section 173.125(1)

Notes: The data provided in this section is for information only. Please apply the appropriate regulations to properly classify your shipment.

## **Section 12 – Regulatory Information**

Active ingredients are not listed by ROTECS, OSHA, IARC, NTP, or EPA. No evidence to date implicates as a carcinogen or toxin.

## **Section 13: Other Information**

Safe n Simple LLC. is a registered Medical Devices manufacturer as designated by the FDA.

Disclaimer: This product is exempt from Safety Data Sheet regulations as the product is for consumer use. (Provided with this information by the compiling agencies): This information contained in this SDS is offered as a guide to the handling of this specific material. The information contained in this Safety Data Sheet (SDS) is offered as a guide to the use and handling of this material. It has been prepared in good faith by technically knowledgeable Personnel and it believed to be correct as of the effective date listed. All safety aspects of all Safe n Simple's products are thoroughly evaluated prior to commercialization Safe n Simple shall not be held liable for any damages, losses or injuries of any nature which may result from the use of or reliance upon any information contained in this SDS. Each individual should make a determination as to the suitability of the information for his or her particular purpose(s). Safe n Simple, LLC and the United States of America in no manner whatsoever, expressly or implied, warrants this information to be accurate and disclaims all liability for its use. Any person utilizing this document should seek competent professional advice to verify and assume responsibility for the suitability of this information to their particular situation.