

# REUSABLE QUILTED UNDERPADS





• Surface  
• Soaker  
• Barrier

**Our underpads are designed with premium components for multi-wash performance:**

- 3-layer quilted lock stitch construction
- Quilted surface fabric
- Heavyweight soaker
- Waterproof barrier
- Rounded corners w/bar tack security

**Four Popular Surface Fabrics:**

- Birdseye 100% cotton for ultimate comfort
- Sheeting 50% cotton-50% polyester for economy
- Brush Polyester for comfort and absorbency.
- Fashion printed woven, an alternative to white, to hide stains


**And Two Barrier Fabrics:**

- Vinyl-supported with knitted backing for economy  
(The green barrier fabric is our standard. Other colors upon request)
- Vinyl-supported woven broadcloth for durability (upon request)

Face Style	Supersize 42" x 35"	X-large 35"x35"	Large 29"x35"	Medium 23"x35"	Chairpad 17"x24"	Chairpad
<b>16"x18" Case Packing</b>	<b>2 DOZEN</b>	<b>2 DOZEN</b>	<b>2 DOZEN</b>	<b>3 DOZEN</b>	<b>5 DOZEN</b>	<b>5 DOZEN</b>
100% cotton birdseye	455-B428G	454-B358G	452-B298G	451-B238G	450-B178G	n/a
50/50 cotton/poly sheeting	548-W428G	544-W358G	542-W298G	541-W238G	540-W178G	n/a
100% brushed poly	470-BP428G	470-BP358G	470-BP298G	470-BP238G	470-BP178G	n/a
50/50 Printed woven	435-P428G	435-P358G	434-P298G	433-P238G	432-P178G	440-P168G

Note: All pads also available in a woven barrier. Call for info.

Many products are also offered with extension tucks for minimal pad shifting



**NOTE:**  
Extension tucks are a durable polycotton sheeting 18" long on both sides of the underpad. To order, add a "T" to the end of style numbers for X-large, large and medium sizes. (e.g.: B358GT for X-large Birdseye)



1-800-935-2241

1

501 Broad Avenue • Ridgely, NJ 07657 | [www.comfortconcepts.com](http://www.comfortconcepts.com) | [info@comfortconcepts.com](mailto:info@comfortconcepts.com) | fax 201-941-6770