



## PROTECTION EVERY STEP OF THE WAY!

The TCC-EZ<sup>®</sup> Total Contact System is an easy-to-use system that supports patient offloading compliance. TCC-EZ is a specialized offloading device designed to provide a protective wound healing environment.

The TCC-EZ is indicated for non-infected wounds with adequate blood supply:

- Diabetic Foot Ulcers located on the plantar, mid-foot or heel
- Post-Operative Care
- Charcot Neuroarthropathy

### Features and Benefits

#### Higher Patient-Reported Preference Over Conventional Total Contact Cast (TCC)

- The lightweight, woven design allows better airflow while wearing the device
- Enforces patient compliance between office visits and supports wound healing

#### Easy to Use and Time-Saving Design

- Prefabricated, woven design has been applied in under 10 minutes<sup>1</sup>
- Requires minimal training time<sup>1</sup>, and up to 4 times as many patients were casted compared with clinics using traditional casting<sup>2</sup>

#### Confidence in Wound Closure

- TCC has demonstrated closure rates of 88% of DFUs in 43 days<sup>3</sup>
- TCC is recognized as the Gold Standard and Preferred method for offloading DFUs<sup>2,4</sup>

## TCC-EZ Healing Chamber®



TCC has been shown to reduce  
amputation rates by more than half  
(Non-TCC 5.2%, TCC 2.2%)<sup>2</sup>



1. In-house data.
2. Fife CE, Carter MJ, Walker D, Thomson B, Eckert KA. Diabetic foot ulcer off-loading: The gap between evidence and practice. Data from the US Wound Registry. *Adv Skin Wound Care.* 2014;27(7):310-316.
3. Bloomgarden ZT. American Diabetes Association 60th Scientific Sessions, 2000. *Diabetes Care.* 2001;24(5):946-951.
4. Armstrong DG, Nguyen HC, Lavery LA, van Schie CH, Boulton AJ, Harkless LB. Off-loading the diabetic foot wound: a randomized clinical trial. *Diabetes Care.* 2001;24(6):1019-1022.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com](http://integralife.com)

INTEGRA®  
LIMIT UNCERTAINTY