

# NITRILECARE EXAMINATION GLOVES

## Specifications & Certifications



### Multi-Purpose




### Product Technical Data Sheet

Type	Powder-Free, Examination Glove
Specification	Non-Sterile/Disposable
Cuff	Beaded
Weight	Average 3.5g (0.12 fl.oz) +/-0.3g
Colour	Indigo Blue
Primary Material	Nitrile
Surface	External: Finger textured
Powder Control	<-2µmg/glove
Packaging	100 pcs/box
Size & Product Code	Small: NG-NTGR-001 Medium: NG-NTGR-002 Large: NG-NTGR-003 Extra-Large: NG-NTGR-004
Product Name	Nitrile care Examination Gloves
Brand	Careplus

Application Setting  
Medium risk medical, dental, procedures, pathology labs, patient care and food handling.

### Certifications



Product Properties			
Test Method	Characteristics	Requirement	Median
ASTM D6319	Width (M)	95 ± 10	95
	Length (M)	Min 230	240 (9.5")
	Thickness single (mm)	Min 0.05	0.10 (4.0mils)
ASTM D5151	Freedom from holes	No leakage	NA
	Before accelerated aging		
ASTM D412	Tensile Strength (MPa)	Min 14	30
	Ultimate elongation (%)	Min 500	650
ASTM D573	After accelerated aging at (70 ±2° C 166±2 hr)		
	Tensile Strength (MPa)	Min 14	14
	Ultimate elongation (%)	Min 400	400
Palm Width	Small: 80mm +/- 10	Weight	3.5g +/- 0.3g (0.12 fl oz)
	Medium: 95mm +/- 10	Length	9.5 inch glove
	Large: 110mm +/- 10	Thickness	3.5 mil glove
	Extra-Large: 120mm +/- 10		
US Standards	ASTM D6319, ASTM F1671 (Viral)		
FDA Information	FDA 510(k) - K172015 Click here for <a href="#">FDA Website</a>		
Quality Standards	ISO 9001:2015, EN ISO 13485:2016		