

EZ-ACCESS Pathway U-Shaped Ramps | Assembled Ramp Specifications

PS20U44TTS

- Length: 20 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 241 Inches
- Weight: 492 Pounds

PS22U44TTS

- Length: 22 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 265 Inches
- Weight: 508 Inches

PS24U44TTS

- Length: 24 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 289 Inches
- Weight: 524 Pounds

PS26U44TTS

- Length: 26 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 313 Inches
- Weight: 558 Pounds

PS28U44TTS

- Length: 28 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 337 Inches
- Weight: 574 Pounds

PS30U44TTS

- Length: 30 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 361 Inches
- Weight: 590 Pounds

EZ-ACCESS Pathway U-Shaped Ramps | Assembled Ramp Specifications

PS20U55TTS

- Length: 20 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 253 Inches
- Weight: 552 Pounds

PS22U55TTS

- Length: 22 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 277 Inches
- Weight: 568

PS24U55TTS

- Length: 24 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 301 Inches
- Weight: 584 Pounds

PS26U55TTS

- Length: 26 Inches
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 325 Inches
- Weight: 618 Pounds

PS28U55TTS

- Length: 28 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 349 Inches
- Weight: 634 Pounds

PS30U55TTS

- Length: 30 Feet
- Width: 112 Inches
- Height: 34.5 Inches
- Depth: 373 Inches
- Weight: 650 Pounds