

KIMTECH SCIENCE* KIMWIPES* Delicate Task Wipers

PRODUCT NAMES & CODES

KIMTECH SCIENCE* KIMWIPES* Delicate Task Wipers 4103002

PRODUCT DESCRIPTION

KIMTECH SCIENCE* KIMWIPES* Delicate Task Wipers are designed specifically for critical environments. Low level of lint makes them ideal for use on delicate parts and instruments such as electrical and electronic equipment, optical products and laboratory equipment.

TECHNICAL FEATURES

- Key Applications include but not limited to:
 - Wiping sensitive surfaces such as glass, coated plastic lenses and safety glasses
 - Suits most lens cleaning stations
 - Idea for use on precision parts, electronic equipment and laboratory equipment.

TECHNICAL SPECIFICATION

Appearance:

White Tissue Paper

Ingredients:

Cellulose fibre	> 98%
Organic Resins	< 0.5%
Inorganic Salts	Trace only

Not classified as hazardous according to the criteria of NOHSC

HEALTH EFFECTS

When used according to the company's directions and for the purpose for which they were intended, these products will not have any adverse health effects.

FIRST AID

Swallowed:	Remove any material causing obstruction.
Eye:	Remove any foreign matter and flush with water. Seek medical attention if irritation persists.
Skin:	Non irritant
Inhalation:	Remove any material causing obstruction.

PRECAUTIONS FOR USE

- Combustible. Will burn if subject to excessive heat or flame.
- Cellulose dust is considered a "nuisance particulate" TLV at 10mg/m³ of total dust

SAFE HANDLING INFORMATION

- Store away from any source of ignition.
- Store away from oxidising agents eg, pool chlorine.
- Fire fighters should wear personal protection and self-contained breathing apparatus.
- Use fire fighting media suitable for principal cause of fire.
- Do not allow product to enter drains or waterways.
- May react with oxidizing agents.

DISPOSAL INFORMATION

Disposal of uncontaminated wiper to ordinary landfill or properly managed municipal incineration is suitable. Dispose contaminated product in accordance with local regulations.

COMPANY DETAILS

Kimberly-Clark Australia Pty Ltd

52 Alfred Street
Milsons Point, NSW 2061

Customer Enquiries	Australia:	1800 647 994
	New Zealand:	0800 447 814

Kimberly-Clark Thailand Ltd

32nd-33rd Floor, United Center
323 Silom Road, Bangrak
Bangkok 10500

Customer Care: 66 2 230 3003

www.kcprofessional.com