

# Atlas Toilevator®

## Features

- Adds 3.5" to toilet seat height
- Rated for users up to 500lb
- Discreet: Blends in with the toilet base
- Sanitary!
- Easy to install
- Kit contains all parts needed for installation
- Color: Paintable white



## Technical specification

### Atlas Toilevator®

#### Toilet Riser – Fits most round and elongated toilets

- 500lb (227kg) maximum weight capacity
- Adds 3.5" to the height of the base of the toilet
- Discreet – no one will realize it is there
- Fits toilets with the base as close as 7" (18cm) from the wall
- May also be painted to match your existing toilet, we suggest using a primer suitable for ABS plastic
- Will fit toilets where the center of the soil pipe is as close as 7" to the wall

Product Code	ATV300
Dimensions	23.25"x11.44"x3.5"
Addition to the toilet's height	3.5"
Weight Capacity*	white
Available Colors	500lb



### Atlas Toilevator Grande®

#### Toilet Riser – Fits newer, larger base toilets

- 500lb (227kg) maximum weight capacity
- Adds 3.5" to the height of the base of the toilet
- Discreet – no one will realize it is there
- Fits toilets where the center of the soil pipe is as close as 7" to the wall
- Toilevator® can also be painted to match your existing toilet, we suggest using a primer suitable

Product Code	ATV350
Dimensions	28.75"x15.25"x3.5"
Addition to the toilet's height	3.5"
Weight Capacity*	white
Available Colors	500lb

