

OXYGEN CONCENTRATOR (5 L)

5 Liter Concentrator Provides Supplemental Oxygen



WE CARE LIKE FAMILY™

FEATURES AND BENEFITS

- Supplies oxygen directly through a nasal cannula
- Equipped with an audible power failure alarm
- Indicator lights notify users of pressure, flow, and oxygen irregularities
- Flow range: 0.5–5 LPM
- Oxygen purity: 90%–96% at all flow rates
- Quiet operation levels and adjustable screen brightness
- Top handle and 4 rolling casters allow for easy transport
- For individual use indoors only
- **HCPCS: E1390**



ALSO AVAILABLE:
Standard and Sof-Touch
Nasal Cannulas
(Reorder No. 5201-5207)

Reorder No.	Product Name	Quantity
10800	Oxygen Concentrator (5 L)	1/Cs

Compare to:

	Dynarex	Invacare	AirSep	Respironics	Drive	Caire
	Oxygen Concentrator (5L)	Perfecto2™	NewLife® Elite	Everflo™	DeVilbiss® 5 Liter Oxygen Concentrator	Companion 5™
Continuous Oxygen Flow (max LPM)	5	5	5	5	5	5
Oxygen Concentration	90 to 96%	87 to 95.6%	90 to 95.5%	93%	87 to 96%	90%
Oxygen Outlet Pressure (max psi)	6–10	5	7–9	5.5	8.5	4.6
Operating Altitude (feet above sea level)	0 to 7,500 or 13,123 as a setting	0 to 8,000	0 to 10,000	0 to 7,500	0 to 13,123	-1,253 to 9,878
Decibel Level	40	43	48	45	48	40
Weight	37 lb.	40 lb.	54 lb.	31 lb.	36 lb.	36 lb.

*Compare Dynarex Oxygen Concentrator to Perfect O2™ by Invacare, NewLife® Elite by AirSep, Everflo™ by Philips Respironics, DeVilbiss® 5 Liter Oxygen Concentrator by Drive DeVilbiss Healthcare, or Companion 5™ by Caire.



Find Us On:



dynarex is your one-stop medical products brand since 1967

For more information, quotes or sales contact a Vitality Medical representative:
Wiley Dangerfield / 800-397-5899 X 378 / wdangerfield@vitalitymedical.com
Keith Radley / 800-397-5899 X 339 / kradley@vitalitymedical.com