

# multiLift™



A flanged pool lift, with left or right side mounting, and optional folding seat version.

- Third-party tested & verified ADA compliant
- Integrated armrests
- State of California compliant
- 350 lb/159kg lifting capacity
- Retrofit anchor jig is standard
- Optional folding seat assembly
- LiftOperator® Intelligent Controller
- Powder-coated stainless steel and aluminum construction



## Lift Color



GRAY MIST

Due to printing technology actual color may differ.



New Construction Jig with Anchors  
500-5000A



Wheel-A-Way mobility option provides flexibility to transport the lift if needed



Optional folding seat assembly

Model No.	Description	Shipping – (Class 150)		
		Weight	Length	Width Height
575-0000	multiLift	200 lbs/91kg		
575-0000N	multiLift, no anchor	200 lbs/91kg	50"	28" 27"
575-1000	multiLift with activation key	200 lbs/91kg	127cm	71cm 69cm
575-1000N	multiLift with activation key, no anchor	200 lbs/91kg		

### Note:

- Includes battery, charger, battery console cover, water-resistant hand control, footrest, seat belt assembly, armrests and retrofit anchor jig.

### Parts & Accessories

- 1001495 Battery
- 500-5200 Cover
- 500-5500 Wheel-A-Way
- 900-1000 Seat Belt
- 300-6700A Anchors, set of 4
- 300-6900 Retrofit Anchor Jig
- 300-6800A Anchor Bolts, set of 4
- 970-5000 Seat Saver Cover
- 900-6000 Stability Strap
- 170-3000A Armrest Replacement (pair)

## multiLift with Folding Seat

Model No.	Description	Shipping – (Class 77.5)		
		Weight	Length	Width Height
575-0100	multiLift with folding seat*	225 lbs/103kg		
575-0100N	multiLift with folding seat, no anchor*	225 lbs/103kg		
575-0105	multiLift with armrests and folding seat	230 lbs/104kg	50"	24" 28"
575-0105N	multiLift with armrests and folding seat, no anchor	230 lbs/104kg	127cm	61cm 71cm
575-1105	multiLift with armrests, folding seat, and activation key	230 lbs/104kg		
575-1105N	multiLift with armrests, folding seat, and activation key, no anchor	230 lbs/104kg		

### Note:

- Includes battery, charger, battery console cover, water-resistant hand control, footrest, seat belt assembly and retrofit anchor jig.

\* Models without armrests are not compliant with State of California requirements.

### Parts & Accessories

- 1001495 Battery
- 500-5100FC Folding Seat Cover
- 160-6000 Folding Seat
- 900-2000 Stability Vest
- 900-4000 Seat Pad
- 970-0000 Seat Saver Cover
- 170-2320 Armrest Assembly, gray, left & right

## New Construction Guidelines

Scenario	Pool Lift	Anchor Jig
Order pool lift and new construction jig at same time	575-3000	500-5000 (no anchors)
Order new construction jig ahead of pool lift	575-3000N	500-5000A (comes with anchors)

A completed Deck Profile form is required with your pool lift order.

