



# Air*Life*™ Tri-Flo suction catheters and tracheostomy care kits

AirLife<sup>™</sup> brand Tri-Flo suction catheters and tracheostomy care kits help you treat the compromised airway. An array of options allow you to choose the optimal product configurations for your patient needs.



# Tri-Flo suction catheters

Tri-Flo catheters are designed with two triangular eyes and a beveled tip to increase suctioning efficiency while reducing patient trauma. The Tri-Flo catheter portfolio offers catheters in 5/6, 8, 10, 12, 14 and 18 Fr sizes to cover a range of patient needs. Clear connectors are standard on all 5/6 and 8 Fr sizes to help you easily visualize patient secretions.





#### Single catheters



Straight with control port

Tri-Flo catheters are available in looped and straight packaging configurations in 5/6, 8, 10, 12, 14 and 18 Fr sizes. Most single catheters are available with or without a control port, so you can choose continuous or controlled suctioning.

Looped with control port

Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
T263C	5/6	100	T63C	5/6	50
T264C	8	100	T64C	8	50
T261C	10	100	T61C	10	50
T268C	12	100	T68C	12	50
T260C	14	100	T60C	14	50
T262C	18	100	T62C	18	50
Straight wit	hout control p	ort	Looped without control port		
T261	10	100	T64	8	50
T260	14	100	T61	10	50
T262	18	100	T60	14	50
			T62	18	50

#### No-touch catheters



Tri-Flo catheters are packaged in a plastic sleeve to add a layer of safety between the clinician and the patient. No-touch catheters are available in 8, 10, 12, 14 and 18 Fr sizes and include a pop-up basin for your solution.

#### With pop-up basin

Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.	
T164C	8	100	T160C	14	100	
T161C	10	100	T162C	18	100	
T168C	12	100				

#### Sterile solutions



For your suction procedure needs, each peel-back foil-lid cup contains 4 oz (120 mL) of sterile water or 0.9% normal saline solution.

Cat. no.	Description	Qty.
4872B	Water	48
4873A	Normal saline	48

### Cath-N-Glove® economy catheter kits



Tri-Flo catheters are packaged with one or two powder-free vinyl (latex-free) exam gloves and a pop-up basin. All economy kits are packaged in a plastic peel pouch to ease use and speed preparation.

One powder-free vinyl glove			Two powder-free vinyl gloves		
Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
4693T	5/6	100	4893T	5/6	100
4697T	8	100	4897T	8	100
4695T	10	100	4895T	10	100
4698T	12	100	4898T	12	100
4694T	14	100	4894T	14	100
4696T	18	100	4896T	18	100

#### Cath-N-Glove wallet catheter kits



Tri-Flo catheters are packaged with two nitrile (*latex-free*) or latex exam gloves and a pop-up basin wrapped together in an inner wallet. The kit is packaged in a plastic peel pouch to ease use and speed preparation.

Two latex gloves			Two nitrile gloves		
Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
4660A	5/6	100	4863T	5/6	100
4661A	8	100	4867T	8	100
4662A	10	100	4865T	10	100
4668A	12	100	4868T	12	100
4664A	14	100	4864T	14	100

#### Rigid basin kits



Tri-Flo catheters are packaged in a rigid basin with two walleted powder-free vinyl (latex-free) exam gloves. The rigid basin can hold your solution while the inner wallet can provide a sterile field during set-up to support aseptic technique.

Two powde	r-free vinyl glo	ves				
Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.	
44-56	5/6	100	44-12	12	100	
44-08	8	100	44-14	14	100	
44-10	10	100	44-18	18	100	

#### Rigid basin kits with saline



Tri-Flo catheters are packaged in a rigid basin with two walleted powder-free nitrile (*latex-free*) exam gloves and a 4 oz container of sterile saline. The inner wallet can be used as a sterile field during set-up to support aseptic technique.

#### Two powder-free nitrile gloves

Cat. no.	Size (Fr)	Qty.	Cat. no.	Size (Fr)	Qty.
41-56	5/6	12	41-12	12	12
41-08	8	12	41-14	14	12
41-10	10	12			

# Tracheostomy care kits

# Complete tracheostomy cleaning trays



We offer a line of sterile, singleuse trays, with or without gloves. These trays include:

- · Three rigid basins for cleaning, soaking and rinsing
- Pre-sewn trach dressing
- Three gauze sponges
- Three pipe cleaners
- Cleaning brush
- Three cotton-tipped applicators
- Twill tape, 32"
- Plastic forceps
- Water-resistant drape
- Hospital wrap

Cat. no.	Additional components	Qty.
4681A	2 vinyl (latex-free) gloves	20
4682A	Without gloves	20

## Tracheostomy care kits



Our kits are supplied in a tray with a coated paper lid. All kits include:

- Trach gauze dressing
- Twill tape, 34"
- Trach brush
- Four gauze sponges, 4" x 4"

Cat. no.	Additional components		Qty.
3T3030A	2 cotton-tipped applicators	3 pipe cleaners	30
3T4691	2 cotton-tipped applicators 2 pipe cleaners	2 vinyl <i>(latex-free)</i> gloves 1 water-resistant towel	20
3T4691A	2 cotton-tipped applicators 3 pipe cleaners	2 vinyl <i>(latex-free)</i> gloves 1 water-resistant towel	30
3T4692	3 cotton-tipped applicators 3 pipe cleaners 2 vinyl ( <i>latex-free</i> ) gloves	1 water-resistant towel 2 packets hydrogen peroxide, 30 mL 1 trach pad	30

Caution—U.S. Federal Law restricts this device to sale by or on the order of a physician.

Manufactured for:

CareFusion 800.231.2466 toll free 22745 Savi Ranch Parkway 714.283.2228 tel Yorba Linda, CA 92887 714.283.8493 fax

CareFusion Yorba Linda, CA

carefusion.com

