



Home About Contact Products Find a Rep Articles & News

HOME » PRODUCTS » BRAND: CLEANIS® » CAREBAG COMMODOE LINER PACK OF 20



CAREBAG Commode Liner Pack of 20 Carebag Commode Liner with Super Absorbent Pad box 20 liners

Item Number: C7831738 by **Carebag**

Find a Dealer

The Carebag® Commode Line with Super Absorbent Pad (20 counts) by Cleanis will make your life easier by allowing you to use a bedside commode without having to clean the pail. Carebag®, is the gift of comfort and dignity for you, your loved one or their caretaker. Very hygienic, significantly reduces the risk of contamination.

Line the commode pail or pot with the bag, The absorbent pad will lie in the bottom to absorb up to 16oz of fluids, turning the collection into a gel. Odors are also trapped in the gel. When collection is complete, simply pull the ties to gather the bag closed, tie off and dispose like a diaper.

Choose the Carebag to match your collection vessel size.

Product instructions:

1. Open box and unroll first liner by tearing along perforated line
2. Position the liner over the commode with the pad laying flat on the bottom
3. After patient uses bedside commode, bodily fluids turn into a gel within seconds
4. To close: Pull straps tightly and tie a knot. Dispose of the sealed bag in household waster like for a diaper. DO NOT FLUSH

Warranty Info

Related Products

Downloads »

brochure

Cleanis brochure

Cleanis order form

Product Manual

Features

Specifications

Special Usage

Length	20.5 in.	Packaging Length	7.75 in.
Width	15.75 in.	Packaging Width	4.5 in.
Height	0.25 in.	Packaging Height	4.5 in.
Product Weight	1.7 lbs.	Package Volume	0.09 cu. in.
Outside Packaging	7.75x4.5x4.5 in.	Shipping Weight	2 lbs.
Latex	n	Product Assembly	no
Primary Material	PE	Primary Color	white
U of M	pack	Unit Quantity	20
Capacity	1 L		

Related Products



OCEAN



TOILEVATOR GRANDE



SAFETY GLO™



QUATTROPLUS

Warranty Information

Manufacturer provides a 30 day warranty against manufacturing defects. Shelf life is 5 years, expiration date is printed on each carton.

Warranty covers cleanis.

SITE DIRECTORY:

Home
About
Contact
Find a Representative
News & Awareness

PRODUCTS FOR:

Capital Equipment
Caregivers
Pediatrics
Rehab Seating and Positioning
Seniors
Therapist
Veterans

PRODUCT CATEGORIES:

ADL - Staying Independent at Home
Bath Safety
Emergency Medical Services
Mobility
Therapy Gym

PENNSYLVANIA:

Clarke Health Care Products
7830 Steubenville Pike
Oakdale, PA 15071

Toll Free: (888) 347-4537

Phone: (724) 695-2122

Fax: (724) 695-2922

info@clarkehealthcare.com

CALIFORNIA:

Clarke Health Care West
4081 Calle Tesoro Unit H
Camarillo, CA 93010

Phone: (805) 388-0612

Fax: (805) 388-6468