



QUALITY, VALUE AND INNOVATION IN WOUND CARE.

You're challenged daily to care for wounds, helping patients heal and avoid infection. You need high quality, high performance products. We can help.

Our extensive wound care portfolio offers an array of options to meet your needs—from gauze bandages and specialty dressings to tapes and adhesives.

We're always looking for ways to improve our products through innovation. Whether you're looking to reduce cross-contamination, save space or improve efficiency, we have just what you need. Here are a few of our latest innovations:

- » Gauze pads with rulers right on the packaging to save clinicians' time
- » ABD pads are folded in thirds rather than halves to save space
- » Our Bulkee® gauze roll is ultra soft with more crimps per inch for better stretching and greater patient comfort

Our products deliver top level performance and quality in the industry—and they're smartly priced.

Table of Contents

Absorbency	4-16
Preserve	17-18
Protection	19-28
Securement	29-41
Cross Reference	42-55

A tiered product approach for every need.

Medline's Premium products are the most feature-rich items that we carry. Medline Essentials are of equal performance, but designed to be the most cost effective.



ABSORBENCY

Contains drainage

Woven Gauze Sponges

100% cotton.

- » Ideal for wound dressings, wound packing and general wound care
- » Strictly controlled manufacturing assures clean, debris-free packing with folded edges to prevent unraveling
- » Sterile 2's packed in coated paper envelopes to reduce fiber debris
- » Easy-open envelopes carefully sealed to prevent the intrusion of dust and contaminants
- » Not made with natural rubber latex

PREMIUM

- » "H fold" design offers a premium approach to minimizing loose threads and lint"

Item No.	Description	Pkg.
Sterile 2's USP Type VII		
NON21420	2" x 2" (5.1 x 5.1 cm), 8-ply	50 pk/bx, 30 bx/cs, 1500 pk/cs
NON21422	3" x 3" (7.6 x 7.6 cm), 12-ply	40 pk/bx, 30 bx/cs, 1200 pk/cs
NON21448	4" x 4" (10.2 x 10.2 cm), 8-ply	50 pk/bx, 12 bx/cs, 600 pk/cs
NON21424	4" x 4" (10.2 x 10.2 cm), 12-ply	25 pk/bx, 24 bx/cs, 600 pk/cs
NON21842	8" x 4" (20.3 x 10.2 cm), 12-ply	25 pk/bx, 20 bx/cs, 500 pk/cs

Sterile 10's (Trays) USP Type VII		
NON21425	4" x 4" (10.2 x 10.2 cm), 8-ply	10/tray, 128 trays/cs, 1280/cs
NON21426	4" x 4" (10.2 x 10.2 cm), 12-ply	10/tray, 128 trays/cs, 1280/cs
NON21428	4" x 4" (10.2 x 10.2 cm), 16-ply	10/tray, 128 trays/cs, 1280/cs

Non-Sterile (Bulk)		
NON25208	2" x 2" (5.1 x 5.1 cm), 8-ply	200/sleeve, 25 sleeves/cs, 5000/cs
NON25212	2" x 2" (5.1 x 5.1 cm), 12-ply	200/sleeve, 40 sleeves/cs, 8000/cs
NON25312	3" x 3" (7.6 x 7.6 cm), 12-ply	200/sleeve, 20 sleeves/cs, 4000/cs
NON25408	4" x 4" (10.2 x 10.2 cm), 8-ply	200/sleeve, 20 sleeves/cs, 4000/cs
NON25412	4" x 4" (10.2 x 10.2 cm), 12-ply	200/sleeve, 10 sleeves/cs, 2000/cs
NON25416	4" x 4" (10.2 x 10.2 cm), 16-ply	200/sleeve, 10 sleeves/cs, 2000/cs
NON25812	8" x 4" (20.3 x 10.2 cm), 12-ply	200/sleeve, 10 sleeves/cs, 2000/cs

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
Sterile 2's, USP Type VII		
PRM21420	2" x 2" (5.1 cm x 5.1 cm), 8-ply	50 pk/bx, 30bx/cs, 1500 pk/cs
PRM21424	2" x 2" (10.2 x 10.2 cm), 12-ply	25 pk/bx, 24 bx/cs, 600 pk/cs
PRM21448	4" x 4" (10.2 x 10.2 cm), 8-ply	25 pk/bx, 24 bx/cs, 600 pk/cs
Sterile 10's (Soft Pouch), USP Type VII		
PRM21426	4" x 4" (10.2 x 10.2 cm), 12-ply	10/pk, 128 pk/cs, 1280/cs
PRM21428	4" x 4" (10.2 x 10.2 cm), 16-ply	10/pk, 128 pk/cs, 1280/cs



NON21842

NON21424

NON21420



LIGHTER WEIGHT

PRM21424

Woven Gauze Sponges (cont.)

MEDLINE ESSENTIALS

» "C fold" design minimizes loose threads and lint

Item No.	Description	Pkg.
Sterile 2's		
PRM2208	2" x 2" (5.1 x 5.1 cm), 8-ply	50 pk/bx, 30 bx/cs, 1500 pk/cs
PRM3312	3" x 3" (7.6 x 7.6 cm), 12-ply	40 pk/bx, 30 bx/cs, 1200 pk/cs
PRM4408	4" x 4" (10.2 x 10.2 cm), 8-ply	50 pk/bx, 12 bx/cs, 600 pk/cs
PRM4412	4" x 4" (10.2 x 10.2 cm), 12-ply	25 pk/bx, 24 bx/cs, 600 pk/cs
Sterile 1's		
PRM21419	2" x 2" (5.1 x 5.1 cm), 12-ply	100 pk/bx, 24 bx/cs, 2400 pk/cs
PRM21423	4" x 4" (10.2 x 10.2 cm), 12-ply	100 pk/bx, 12 bx/cs, 1200 pk/cs
Non-Sterile		
PRM21208C	2" x 2" (5.1 x 5.1 cm), 8-ply	200/sleeve, 25 pk/cs, 5000/cs
PRM21212C	2" x 2" (5.1 x 5.1 cm), 12-ply	200/sleeve, 40 pk/cs, 8000/cs
PRM21312C	3" x 3" (7.6 x 7.6 cm), 12-ply	200/sleeve, 20 pk/cs, 4000/cs
PRM21408C	4" x 4" (10.2 x 10.2 cm), 8-ply	200/sleeve, 20 pk/cs, 4000/cs
PRM21412C	4" x 4" (10.2 x 10.2 cm), 12-ply	200/sleeve, 10 pk/cs, 2000/cs
PRM21416C	4" x 4" (10.2 x 10.2 cm), 16-ply	200/sleeve, 10 pk/cs, 2000/cs
Cotton Filled Dental Sponges		
PRM22DS	2" x 2" (5.1 x 5.1 cm)	200/sleeve, 25 pk/cs, 5000/cs



X-Ray Detectable Woven Gauze Sponges

100% cotton.

- » Cutting-edge manufacturing technology minimizes lint
- » Highly absorbent
- » Dual X-ray detectable barium sulfate strands layered into the sponges in a unique figure 8 design for easy detection
- » Available in multiple sizes and plies
- » Not made with natural rubber latex
- » Skillfully manufactured to ensure gauze can be removed cleanly from tray

↳ Exceptionally manufactured using cutting-edge technology to minimize lint



NON21430LF

PREMIUM

Item No.	Description	Pkg.
Sterile (Trays) USP Type VII		
NON21430LF	4" x 4" (10.2 x 10.2 cm), 16-ply	10/tray, 128 trays/cs, 1280/cs
NON21432LF	8" x 4" (20.3 x 10.2 cm), 12-ply	10/tray, 80 trays/cs, 800/cs
NON21433	4" x 4" (10.2 x 10.2 cm), 32-ply	10/tray, 48 trays/cs, 480/cs
NON21438	8" x 4" (20.3 x 10.2 cm), 24-ply	10/tray, 20 trays/cs, 200/cs
Non-Sterile (Bulk)		
NON25430	4" x 4" (10.2 x 10.2 cm), 16-ply	10/band, 100/sleeve, 2000/cs

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
USP Type VII X-Ray Detectable, Sterile, Soft Pouch		
PRM21430LF	4" x 4" (10.2 x 10.2 cm), 16-ply	10/pk, 135 pk/cs, 1350/cs

100% Cotton Gauze Bandages

Ideal as a primary or secondary dressing.

- » Open weave design provides fast wicking, absorbency and aeration
- » Tightly finished edge minimizes lint
- » Increased crimps for improved cushioning and patient comfort
- » 6-ply bandage stretches and conforms for a variety of patient needs, especially burn care
- » Not made with natural rubber latex

PREMIUM

Bulkee II

Compare to Kerlix®*

Item No.	Description	Pkg.
Sterile (Singles), Stretched		
NON25865	4.5" x 4.1-yd. (11.4 cm x 3.8 m) 6-ply	100 rolls/cs
NON25861	3.4" x 3.6-yd. (8.6 cm x 3.3 m) 6-ply	96 rolls/cs
NON25850	2.25" x 3-yd. (5.7 cm x 2.7 m) 6-ply	12/bx, 96 rolls/cs

Non-Sterile, Stretched

NON25855	4.5" x 4.1-yd. (11.4 cm x 3.8 m) 6-ply	100 rolls/cs
NON25855P	4.5" x 4.1-yd. (11.4 cm x 3.8 m) 6-ply	5 rolls/pk

* This product is not manufactured by Medtronic, makers of Kerlix® dressings, bandages and sponges.

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
PRM25865	4.5" x 4.1-yd. (11.4 cm x 3.8 m), Sterile 6-ply	100 rolls/cs
PRM25865P	4.5" x 4.1-yd. (11.4 cm x 3.8 m), Sterile 6-ply	5 rolls/pk
PRM25865H	4.5" x 4.1-yd. (11.4 cm x 3.8 m), Sterile 6-ply	1/ea
PRM25855	4.5" x 4.1-yd. (11.4 cm x 3.8 m), Non-Sterile 6-ply	100 rolls/cs



NON25865



PRM25865

30%
LIGHTER IN
WEIGHT

Bulkee Lite® 100% Cotton Conforming Bandages

Compare to Kerlix® Lite*

Recommended as a secondary dressing to secure IVs, dressings, splints.

- » Cushioning and a light stretch helps this bandage conform to body contours
- » 3-ply woven bandage provides mild compression
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile (Singles), Stretched		
NON27496	2" x 3.5-yd. (5.1 cm x 3.2 m) 3-ply	12/bx, 8 bx/cs, 96/cs
NON27497	3" x 4.1-yd. (7.6 cm x 3.8 m) 3-ply	12/bx, 8 bx/cs, 96/cs
NON27498	4" x 4.1-yd. (10.2 cm x 3.8 m) 3-ply	12/bx, 8 bx/cs, 96/cs
NON27499	6" x 4.5-yd. (15.2 cm x 4.1 m) 3-ply	6/bx, 8 bx/cs, 48/cs

Non-Sterile (Bulk), Stretched		
NON27492	2" x 3.5-yd. (5.1 cm x 3.2 m) 3-ply	12/bg, 8 bg/cs, 96/cs
NON27493	3" x 4.1-yd. (7.6 cm x 3.8 m) 3-ply	12/bg, 8 bg/cs, 96/cs
NON27494	4" x 4.1-yd. (10.2 cm x 3.8 m) 3-ply	12/bg, 8 bg/cs, 96/cs
NON27495	6" x 4.5-yd. (15.2 cm x 4.1 m) 3-ply	6/bg, 8 bg/cs, 48/cs

* This product is not manufactured by Medtronic, makers of Kerlix® dressings, bandages and sponges.



NON27496



NON27499

Conforming Bandages

- » One-ply, rayon and polyester blend
- » Gently secures to body contours without restricting movement
- » Accommodates an increase or decrease in swelling
- » Uniquely constructed weave and finished edges help eliminate unraveling and lint
- » Provides stretch memory when applied
- » Ideal for secondary dressing applications
- » Not made with natural rubber latex

PREMIUM

Compare to Kendall Conform®*

Item No.	Description	Pkg.
Sterile (Singles), Relaxed		
NON254955	1" x 75" (2.5 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
NON25496	2" x 75" (5.1 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
NON25497	3" x 75" (7.6 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
NON25498	4" x 75" (10.2 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
NON25499	6" x 80" (15.2 x 203.2 cm), 1-ply	6/bx, 8 bx/cs, 48/cs

Non-Sterile (Bulk), Relaxed

NON25491	1" x 75" (2.5 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
NON25492	2" x 75" (5.1 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
NON25493	3" x 75" (7.6 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
NON25494	4" x 75" (10.2 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
NON25495	6" x 80" (15.2 x 203.2 cm), 1-ply	6/bg, 8 bg/cs, 48/cs

* This product is not manufactured or distributed by Medtronic, makers of Conform® gauze bandages.

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
Sterile – Unstretched		
PRM25496	2" x 75" (5.1 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
PRM25497	3" x 75" (7.6 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
PRM25498	4" x 75" (10.2 x 190.5 cm), 1-ply	12/bx, 8 bx/cs, 96/cs
PRM25499	6" x 80" (15.2 x 203.2 cm), 1-ply	6/bx, 8 bx/cs, 48/cs

Non-Sterile – Unstretched

PRM25492	2" x 75" (5.1 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
PRM25493	3" x 75" (7.6 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
PRM25494	4" x 75" (10.2 x 190.5 cm), 1-ply	12/bg, 8 bg/cs, 96/cs
PRM25495	6" x 80" (15.2 x 203.2 cm), 1-ply	6/bg, 8 bg/cs, 48/cs



NON25498



PRM25497

LIGHTER WEIGHT

Deluxe Non-Woven Gauze Sponges

Unique rayon and polyester blend
Our most absorbent line

- » Unique rayon and polyester blend provides strength and absorbency
- » Increased material weight helps provide a more efficient sponge than our Medline Essentials line by reducing the number of sponges that have to be used for each procedure
- » Not made with natural rubber latex
- » Unique benefits compared to woven cotton gauze:
 - Low linting
 - Increased absorbency
 - Fast wicking
 - Economical
 - Less traumatic to granulating tissue



NON26443

DELUXE

Item No.	Description	Pkg.
Sterile 2's		
NON264442	4" x 4" (10.2 x 10.2 cm), 4-ply	50 pk/bx, 12 bx/cs, 600 pks/cs
Non-Sterile (Bulk)		
NON26223	2" x 2" (5.1 x 5.1 cm), 3-ply	200/sleeve, 40 sleeves/cs, 8000/cs
NON26224	2" x 2" (5.1 x 5.1 cm), 4-ply	200/sleeve, 40 sleeves/cs, 8000/cs
NON26334	3" x 3" (7.6 x 7.6 cm), 4-ply	200/sleeve, 20 sleeves/cs, 4000/cs
NON26443	4" x 4" (10.2 x 10.2 cm), 3-ply	200/sleeve, 10 sleeves/cs, 2000/cs
NON26444	4" x 4" (10.2 x 10.2 cm), 4-ply	200/sleeve, 10 sleeves/cs, 2000/cs

Non-Woven Gauze Sponges

Rayon and polyester blend

- » Produces less lint and is less traumatic to fragile granulating wound tissue than traditional woven cotton gauze
- » Virtually non-adherent
- » Significantly faster wicking ability promotes a healthy healing environment
- » Discover the perfect balance between strength and absorbency
- » Not made with natural rubber latex



NON21444

PREMIUM

Item No.	Description	Pkg.
Sterile 2's		
NON21224	2" x 2" (5.1 x 5.1 cm), 4-ply	50 pk/bx, 30 bx/cs, 1500 pks/cs
NON21334	3" x 3" (7.6 x 7.6 cm), 4-ply	40 pk/bx, 30 bx/cs, 1200 pks/cs
NON21444	4" x 4" (10.2 x 10.2 cm), 4-ply	50 pk/bx, 12 bx/cs, 600 pks/cs
NON21446	4" x 4" (10.2 x 10.2 cm), 6-ply	25 pk/bx, 12 bx/cs, 300 pks/cs
Sterile 10's (Trays)		
NON21449	4" x 4" (10.2 x 10.2 cm), 4-ply	10/tray, 128 trays/cs, 1280/cs
Non-Sterile (Bulk)		
NON25223	2" x 2" (5.1 x 5.1 cm), 3-ply	200/sleeve, 25 sleeves/cs, 5000/cs
NON25224	2" x 2" (5.1 x 5.1 cm), 4-ply	200/sleeve, 40 sleeves/cs, 8000/cs
NON25334	3" x 3" (7.6 x 7.6 cm), 4-ply	200/sleeve, 20 sleeves/cs, 4000/cs
NON25443	4" x 4" (10.2 x 10.2 cm), 3-ply	200/sleeve, 20 sleeves/cs, 4000/cs
NON25444	4" x 4" (10.2 x 10.2 cm), 4-ply	200/sleeve, 10 sleeves/cs, 2000/cs
NON25449	4" x 4" (10.2 x 10.2 cm), 4-ply	25/sleeve, 100 sleeves/cs, 2500/cs

Non-Woven Gauze Sponges (cont.)

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
Sterile 2's		
PRM21224	2" x 2" (5.1 x 5.1 cm), 4-ply	50 pk/bx, 30 bx/cs, 1500 pk/cs
PRM21444	4" x 4" (10.2 x 10.2 cm), 4-ply	50 pk/bx, 12 bx/cs, 600 pk/cs
Non-Sterile Bulk		
PRM25224	2" x 2" (5.1 x 5.1 cm), 4-ply	200/sleeve, 40 sleeves/cs, 8000/cs
PRM25444	4" x 4" (10.2 x 10.2 cm), 4-ply	200/sleeve, 10 sleeves/cs, 2000/cs



Bulkee II® Super Fluff Sponges

Compare to Kerlix® Fluff Sponges*

All the great features of Bulkee II in a sponge.

- » An excellent primary dressing with ideal loft, bulk and absorbency
- » Effectively manages exudate in and around the wound
- » Ideal for pre-op prepping and cleansing, debriding and packing
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile		
NON25852	6" x 6.75" (15.2 x 17.1 cm), diagonally	2/pk, 20 pk/bx, 12 bx/cs, 480 sponges/cs
NON25853	6" x 6.75" (15.2 x 17.1 cm), diagonally	5/tray, 120 trays/cs, 600 sponges/cs
NON25854	6" x 6.75" (15.2 x 17.1 cm), diagonally	10/tray, 48 trays/cs, 480 sponges/cs
NON25856	7.75" x 8.75" (19.7 x 22.2 cm), diagonally	2/pk, 25 pk/bx, 8 bx/cs, 400 sponges/cs
NON25844	4" x 4" (10.2 x 10.2 cm), 12-ply	2/pk, 25 pk/bx, 24 bx/cs, 1200 sponges/cs
NON25884	4" x 4" (10.2 x 10.2 cm), 12-ply	10/tray, 128 trays/cs, 1280 sponges/cs
Non-Sterile		
NON1167B	6" x 6¾" (15.2 x 17.1 cm)	100/sleeve, 12 sleeves/cs, 1200/cs

* This product is not manufactured by Medtronic, makers of Kerlix dressings, bandages and sponges.



Drain Sponges

Rayon and polyester blend.

- » Precisely fenestrated for easy and consistent application
- » Unique fenestration allows sponge to lie flat when wrapped around a drainage tube
- » Pre-cut sponge eliminates the need for sterile scissors and conserves nursing time
- » 6-ply, highly absorbent, low-linting, fast wicking
- » All the benefits of non-woven
- » Not made with natural rubber latex

PREMIUM

Item No.	Description	Pkg.
NON256000	4" x 4" (10.2 x 10.2 cm), Sterile 2's	25/bx, 12 bx/cs 300 pk/cs
NON256001	2" x 2" (5.1 x 5.1 cm), Sterile 2's	35/bx, 20 bx/cs 700 pk/cs

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
PRM256001	2" x 2" (5.1 x 5.1 cm)	35/bx, 20 bx/cs, 700 pk/cs
PRM256000	4" x 4" (10.2 x 10.2 cm)	25/bx, 12 bx/cs, 300 pk/cs
PRM256000H	4" x 4" (10.2 x 10.2 cm)	25/bx

100% Cotton Drain Sponges

- » Natural breathability
- » Tightly woven into fenestrated shape for your convenience
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON255000	4" x 4" (10.2 x 10.2 cm) 16-ply, Sterile	1/pk, 25/bx, 12 bx/cs

Post-Op Gauze Sponges

More absorbent than cotton gauze sponges.

- » Filled with extra-absorbent cellulose fibers for situations that require increased absorbency
- » Outer layer is constructed of absorbent non-woven rayon and polyester fibers
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile 2's		
NON21442	4" x 4" (10.2 x 10.2 cm)	25/bx, 24 bx/cs, 600 pk/cs
NON21441	4" x 3" (10.2 x 7.6 cm)	25/bx, 24 bx/cs, 600 pk/cs
Non-Sterile (Bulk)		
NON25442	4" x 4" (10.2 x 10.2 cm)	100/sleeve, 20 sleeves/cs, 2000/cs
NON25441	4" x 3" (10.2 x 7.6 cm)	100/sleeve, 20 sleeves/cs, 2000/cs



Drain sold separately



30% LIGHTER WEIGHT

PRM256000



NON255000



NON21441

Abdominal (ABD) Pads

Features extra absorbency and wicking for heavily draining wounds.

- » Soft non-woven outer layer quickly wicks fluid to a cellulose center
- » Cellulose quickly absorbs and disperses fluids laterally to prevent pooling
- » Hydrophobic back impedes strike-through
- » Sealed edges prevent linting
- » Not made with natural rubber latex

PREMIUM

Item No.	Description	Pkg.
Sterile 1's		
NON21450	5" x 9" (12.7 x 22.9 cm)	25/bx, 16 bx/cs, 400/cs
NON21453	8" x 7.5" (20.3 x 19.1 cm)	20/bx, 12 bx/cs, 240/cs
NON21454	8" x 10" (20.3 x 25.4 cm)	18/bx, 20 bx/cs, 360/cs
Non-Sterile (Bulk)		
NON21451	5" x 9" (12.7 x 22.9 cm)	576/cs
NON21452	8" x 7.5" (20.3 x 19.1 cm)	576/cs
NON21456	8" x 10" (20.3 x 25.4 cm)	432/cs
NON21457	12" x 16" (30.5 x 40.6 cm)	144/cs

MEDLINE ESSENTIALS

Item No.	Description	Pkg.
Sterile 1's		
PRM21450	5" x 9" (12.7 x 22.9 cm)	25/bx, 16 bx/cs, 400/cs
PRM21453	8" x 7.5" (20.3 x 19.1 cm)	20/bx, 12 bx/cs, 240/cs
PRM21454	8" x 10" (20.3 x 25.4 cm)	16/bx, 20 bx/cs, 320/cs
Non-Sterile Bulk		
PRM21451	5" x 9" (12.7 x 22.9 cm)	576/cs
PRM21456	8" x 10" (20.3 x 25.4 cm)	432/cs



NON21450



LIGHTER WEIGHT

PRM21450

Non-Adherent Pads

Sterile, cotton and polyester blend.

- » Protects granulating tissue by not sticking to the lightly draining wound
- » Soft, perforated mylar film is porous, allowing air to circulate
- » Absorbent pad allows fluid to collect
- » Not made with natural rubber latex

PREMIUM

CURAD®

Compare to Telfa®*

Item No.	Description	Pkg.
Sterile 1's		
NON25700	2" x 3" (5.1 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25710	3" x 4" (7.6 x 10.2 cm)	100/bx, 12 bx/cs, 1200/cs
NON25720	3" x 8" (7.6 x 20.3 cm)	50/bx, 12 bx/cs, 600/cs

Sterile 1's with Adhesive Tabs

NON25700A	2" x 3" (5.1 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25710A	3" x 4" (7.6 x 10.2 cm)	100/bx, 12 bx/cs, 1200/cs

*This product is not manufactured by Medtronic, makers of Telfa® wound dressings.



MEDLINE ESSENTIALS

PRM25700	2" x 3" (5.1 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
PRM25710	3" x 4" (7.6 x 10.2 cm)	100/bx, 12 bx/cs, 1200/cs

Multi-Trauma Sterile Dressing

Item No.	Description	Pkg.
NON21459	10" x 30" (25.4 x 76.2 cm), ABD Pad	50/cs
DYNDRAINPACK	Set includes: (1) ABD Pad, 10" x 12" (25.4 x 30.5 cm) (6) Bulkee Fluff Sponges, 4"x 4" 12-ply, (10.2 x 10.2 cm)	48/cs



DYNDRAINPACK

Burn Dressings

100% cotton.

- » Provides loft and extra cushioning
- » Wide mesh dressings absorb wound drainage
- » Designed for use on a wide range of burns
- » Available in a variety of sizes, ply, sewn or unsewn edges
- » Not made with natural rubber latex



Item No.	Description	Pkg.
Sewn		
NON7911	18" x 18" (45.7 x 45.7 cm), 10-ply	1/pk, 80 pk/cs
NON7912	18" x 18" (45.7 x 45.7 cm), 2-ply	2/pk, 100 pk/cs
NON7913	18" x 18" (45.7 x 45.7 cm), 4-ply	5/pk, 40 pk/cs
NON216645	6" x 6", (15.2 x 15.2 cm) 45-ply	6/tray, 20 trays/cs

Unsewn

NON7915	18" x 18" (45.7 x 45.7 cm), 50-ply	1/pk, 20 pk/cs
NON7916	36" x 36" (91.4 x 91.4 cm), 18-ply	1/pk, 20 pk/cs



NON7911



SORBEX™ Sterile Absorbent Dressings
Compare to Exu-Dry®*

A highly absorbent primary or secondary wound dressing.

- » Does not stick or cause trauma to wound sites during dressing changes
- » Softly conforms to body contours
- » Cushions and protects the wound
- » Recommended for the management of superficial to partial thickness wounds
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NONEXDISC	3" (7.6 cm) Disc	100/cs
NONEX0023F	2" x 3" (5.1 x 7.6 cm) Slit	25/bx, 2 bx/cs, 50/cs
NONEX0034	3" x 4" (7.6 x 10.2 cm) Full	50/bx, 2 bx/cs, 100/cs
NONEX0046	4" x 6" (10.2 x 15.2 cm) Full	10/bx, 12 bx/cs, 120/cs
NONEX0069	6" x 9" (15.2 x 22.8 cm) Full	12/bx, 4 bx/cs, 48/cs
NONEX0915	9" x 15" (22.8 x 38.1 cm) Full	30/cs
NONEX1524	15" x 24" (38.1 x 60.9 cm) Full	30/cs

* This product is not manufactured by Smith and Nephew, maker of the Exu-Dry® wound dressings

Cotton Rolls

One pound of 100% cotton wrapped and neatly rolled in paper.

- » Highly absorbent surgical grade cotton
- » Ideal for wound cleansing, application of medications, padding of wound area
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON6028	1 lb. Roll (0.5 kg), Sterile	10 rolls/cs
NON6027	1 lb. Roll (0.5 kg), Non-Sterile	25 rolls/cs



NON6028

Cotton Balls and Rayon Balls

Soft and absorbent.

- » Designed for many uses around any facility
- » Not made with natural rubber latex

Item No.	Description	Pkg.
----------	-------------	------

Cotton Balls, Non-Sterile

MDS21460	1" diameter (2.5 cm)	2000/bg, 2 bg/cs, 4000/cs
MDS21461	1" diameter (2.5 cm)	200/bg, 20 bg/cs, 4000/cs
MDS21462	1¼" diameter (3.2 cm)	1000/bg, 2 bg/cs, 2000/cs
MDS21463	1¼" diameter (3.2 cm)	50/bg, 80 bg/cs, 4000/cs

Rayon Balls, Non-Sterile

MDS21470	1" diameter (2.5 cm)	2000/bg, 2 bg/cs, 4000/cs
MDS21472	1.25" diameter (3.2 cm)	1000/bg, 2 bg/cs, 2000/cs



MDS21460

Eye Pads

Sterile.

- » Made with a soft non-woven cover over an absorbent cotton filler
- » Compressed edges help reduce linting
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON21600	1.625" x 2.625" (4.1 x 6.7 cm)	1/pk, 50 pk/bx, 12 bx/cs, 600/cs
NON21601	2.125" x 2.625" (5.4 x 6.7 cm)	1/pk, 50 pk/bx, 12 bx/cs, 600/cs



NON21600

Nursing Pad

Soft and highly absorbent.

- » Helps wick fluids away from the skin
- » Adhesive on the back for easy securement to undergarments
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NURSEPAD1	Disposable nursing pad	12/bx, 288/cs



NURSEPAD1



CURAD® Packing Strips

Use for packing or as drainage conduits in wound care management applications.

- » Tightly woven, 100% cotton gauze
- » Ravel resistant; won't fray or lint
- » Available plain or impregnated with iodoform, an antiseptic agent
- » Packaged as one 5-yard strip per bottle
- » Recommended for: Nasal or sinus packing, boils, abscesses, fistulas and other draining or tunneling wounds
- » Not made with natural rubber latex



NON256015

NON255015

Item No.	Description	Pkg.
Sterile, Plain		
NON255015	1" x 5-yd. (2.5 cm x 4.6 m)	12/cs
NON255025	2" x 5-yd. (5.1 cm x 4.6 m)	12/cs
NON255125	.5" x 5-yd. (1.3 cm x 4.6 m)	12/cs
NON255145	.25" x 5-yd. (0.6 cm x 4.6 m)	12/cs
Sterile, Iodoform		
NON256015	1" x 5-yd. (2.5 cm x 4.6 m)	12/cs
NON256025	2" x 5-yd. (5.1 cm x 4.6 m)	12/cs
NON256125	.5" x 5-yd. (1.3 cm x 4.6 m)	12/cs
NON256145	.25" x 5-yd. (0.6 cm x 4.6 m)	12/cs

PRESERVE

Safeguard granulating tissue



CUR250383



CUR253180

CURAD® Oil Emulsion Dressing

Compare to Adaptic® Non-Adhering Dressing*
Knitted, high porosity cellulose acetate is evenly impregnated with U.S.P. white petrolatum.

- » Emulsion blend permits the flow of exudates without adhering to granulating tissue
- » Minimizes skin maceration
- » Highly conformable
- » Use directly on wounds beneath an absorbent secondary dressing
- » Non-occlusive dressing
- » Not made with natural rubber latex
- » Ideal for most draining wounds including:
 - Minor burns
 - Suture lines
 - Lacerations
 - Skin ulcers
 - Episiotomy incisions
 - Abrasions

Item No.	Description	Pkg.
CUR250383	3" x 8" (7.6 x 20.3 cm)	3/pk, 108/bx, 648/cs
CUR250330	3" x 3" (7.6 x 7.62 cm)	50/bx, 600/cs
CUR250381	3" x 8" (7.6 x 20.3 cm)	1/pk, 24/bx, 144/cs
CUR250317	3" x 16" (7.6 x 40.6 cm)	36/bx, 216/cs

*This product is not manufactured or distributed by Johnson & Johnson Corporation, distributors of Adaptic® non-adherent dressings.

CURAD® Xeroform Gauze Dressing

Petrolatum gauze with the addition of 3% Xeroform (Bismuth Tribromophenate), a deodorizing agent.

- » Contraindicated for heavily draining wounds
- » Occlusive dressing
- » Not made with natural rubber latex
- » Ideal for non-draining to lightly draining wounds including:
 - Donor sites
 - Lacerations
 - Minor or partial thickness burns
 - Abrasions
 - Skin graft sites
 - Suture wounds

Item No.	Description	Pkg.
Overwrap Gauze Dressing, Sterile		
CUR253590W	5" x 9" (12.7 x 22.9 cm)	50/bx, 200/cs
CUR253180W	1" x 8" (2.5 x 20.3 cm)	50/bx, 200/cs
Gauze Dressing, Sterile		
CUR253180	1" x 8" (2.5 x 20.3 cm)	50/bx, 200/cs
CUR253590	5" x 9" (12.7 x 22.9 cm)	50/bx, 200/cs
CUR253220	2" x 2" (5.1 x 5.1 cm)	25/bx, 150/cs
CUR253440	4" x 4" (10.2 x 10.2 cm)	25/bx, 150/cs

CURAD® Petrolatum Gauze Dressing

An effective, non-adherent primary dressing that maintains a moist wound bed for autolytic healing.

- » U.S.P. cotton gauze is evenly impregnated with U.S.P. white petrolatum
- » Soothing and conforming
- » Occlusive dressing
- » Won't stick to the wound site
- » Not made with natural rubber latex
- » Ideal for use on non-draining wounds including:
 - Tube sites
 - Umbilical bandage
 - Stomas
 - Non-infected wounds
 - Skin grafts
 - Donor sites
 - Circumcision
 - Lacerations
 - Minor burns
 - Abrasions

Item No.	Description	Pkg.
Overwrap Gauze Dressing, Sterile		
CUR251136W	1" x 36" (2.5 x 91.4 cm)	12/bx, 72/cs
CUR251272W	½" x 72" (1.3 x 182.9 cm)	12/bx, 72/cs
CUR251390W	3" x 9" (7.6 x 22.9 cm)	12/bx, 72/cs
CUR251318W	3" x 18" (7.6 x 45.7 cm)	12/bx, 72/cs

Gauze Dressing, Sterile		
CUR251136	1" x 36" (2.5 x 91.4 cm)	12/bx, 72/cs
CUR251272	½" x 72" (1.3 x 182.9 cm)	12/bx, 72/cs
CUR251336	3" x 36" (7.6 x 91.4 cm)	12/bx, 72/cs
CUR251390	3" x 9" (7.6 x 22.9 cm)	50/bx, 200/cs
CUR251180	1" x 8" (2.5 x 20.3 cm)	50/bx, 200/cs
CUR251318	3" x 18" (7.6 x 45.7 cm)	12/bx, 72/cs
CUR251636	6" x 36" (15.2 x 91.4 cm)	12/bx, 72/cs

Sterile Disposable Drapes

Create a sterile field with a sterile paper drape.

- » Repellent plastic liner helps prevent strike through
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON21001	18" x 26" (45.7 x 66 cm), Paper drape	50/bx, 300/cs
NON21002	18" x 26" (45.7 x 66 cm), Paper Drape with 2¾" Fenestration	50/bx, 300/cs



CUR251136



NON21001

PROTECTION

Guards from external contaminants

CURAD® Plastic Adhesive Bandages

Rounded edges for great fit and contour.

- » Perforated for breathability
- » Even adhesion
- » Absorbent non-stick island pad
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
NON25500	¾" x 3" (1.9 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25600	1" x 3" (2.5 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25501	Spot ⅞" (2.2 cm)	100/bx, 12 bx/cs, 1200/cs
NON25504	X-Lg 2" x 4" (5.1 x 10.2 cm)	50/bx, 12 bx/cs, 600/cs
NON25509	Jr. ⅝" x 1½" (0.9 x 3.8 cm)	100/bx, 36 bx/cs, 3600/cs
NON25515	Patch, 1½" x 1½" (3.8 x 3.8 cm)	100/bx, 12 bx/cs, 1200/cs



NON25600

MEDLINE ESSENTIALS

Plastic Adhesive Bandages

Ideal for bandaging minor wounds and scrapes.

- » Tan plastic backing is water resistant
- » Long-lasting adhesive with minimal residue
- » Absorbent non-adherent continuous pad
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
PRM25500	¾" x 3" (1.9 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
PRM25500H	¾" x 3" (1.9 x 7.6 cm)	100/bx
PRM25600	1" x 3" (2.5 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
PRM25600H	1" x 3" (2.5 x 7.6 cm)	100/bx



PRM25500

RETAIL FRIENDLY

CURAD® Plastic Bandages

Soft and strong—ideal for everyday use.

- » Made from flexible plastic reinforced with a layer of non-woven material
- » Skin-friendly adhesive won't leave a sticky residue on the skin
- » Ensures maximum comfort as you heal
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR02278RB	Plastic, ¾" x 2½" (1.9 x 7.3 cm)	80/bx, 24 bx/cs
CUR45157RB	Plastic Bandages, Assorted:	80/bx, 24 bx/cs
	20 – Spot ⅞" (2.2 cm)	
	20 – ¾" x 2½" (1.9 x 5.4 cm)	
	30 – ¾" x 2½" (1.9 x 7.3 cm)	
CUR0655RB	10 – 1½" x 2½" (2.9 x 5.4 cm)	8/bx, 24bx/cs
	Plastic, Knee and Elbow, 2¾" x 2¾" (69 x 69 mm)	
CUR43021RB	Waterproof, 1" x 3¼" (2.5 x 8.3 cm)	20/bx, 24 bx/cs
CUR44010RB	Clear, ¾" x 3" (1.9 x 7.6 cm)	30/bx, 24 bx/cs



CUR45157RB

CURAD® Flex-Fabric™ Adhesive Bandages

Material naturally breathes.

- » Stretches and conforms to wounds
- » Even adhesion
- » Absorbent non-stick island pad
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
NON25650	¾" x 3" (1.9 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25660	1" x 3" (2.5 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25502	Spot ⅞" (2.2 cm)	100/bx, 12 bx/cs, 1200/cs
NON25510	Knuckle, 1½" x 3" (3.8 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25513	Fingertip, 1¼" x 2" (4.4 x 5.1 cm)	100/bx, 12 bx/cs, 1200/cs
NON25524	X-Large 2" x 4" (5.1 x 10.2 cm)	50/bx, 12 bx/cs, 600/cs
NON25650B	¾" x 3" (1.9 x 7.6 cm)	2025/bx, 4 bx/cs, 8100/cs
NON25660B	1" x 3" (2.5 x 7.6 cm)	1800/bx, 4 bx/cs, 7200/cs



NON25510

RETAIL FRIENDLY

CURAD® Flex-Fabric™ Bandages

Designed for active lifestyles.

- » Stays in place even when wet
- » Flex-Fabric stretches and breathes
- » Water-repellent coating keeps the bandages dry and secure
- » Skin-friendly adhesive won't leave a sticky, itchy residue on the skin
- » Maximum comfort as you heal
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR0700RB	Flex-Fabric, Assorted	100/bx, 24 bx/cs
	6 – 2" x 4" (5.1 x 10.2 cm)	
	10 – 1" x 3" (2.5 x 7.6 cm)	
	44 – ¾" x 3" (1.9 x 7.6 cm)	
CUR45246RB	40 – ⅝" x 2¼" (1.6 x 5.7 cm)	20/bx, 24 bx/cs
	10 – Fingertip, 1¼" x 2" (4.4 x 5.1 cm)	
	10 – Knuckle, 1½" x 3" (3.8 x 7.6 cm)	



CUR0700RB

CURAD® Comfort Fabric Adhesive Bandages

Constructed of a soft, non-woven gentle material.

- » Breathable, comfortable and durable
- » Even adhesion stretches with your skin
- » Absorbent non-stick island pad
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
NON25750	¾" x 3" (1.9 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25760	1" x 3" (2.5 x 7.6 cm)	100/bx, 12 bx/cs, 1200/cs
NON25744	X-Large 2" x 4" (5.1 x 10.2 cm)	50/bx, 12 bx/cs, 600/cs



NON25750

CURAD® Sensitive Skin Bandages

Ideal for sensitive skin.

- » Specially designed to be gentle on the skin
- » Ultra-light fabric is breathable, comfortable, hypoallergenic
- » Absorbent non-stick pad is touchless to remove
- » Innovative adhesive stays on longer
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR45230NRB	50 – Spot 7/8" (2.2 cm)	50/bx, 24 bx/cs
CUR45227NRB	30 – 3/4" x 3" (1.9 x 7.6 cm)	30/bx, 24 bx/cs



CUR45227NRB



CUR45230NRB

CURAD® Extreme Hold Bandages

Designed to protect hardworking hands.

- » Heavy-duty fabric provides maximum protection
- » Sure stick adhesive ensures the bandage stays in place
- » Special shapes for knuckles and fingertips
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR14924RB	Extreme Hold, Assorted 4 – 2" x 3 3/4" (5.1 x 9.5 cm) 4 – Knuckle, 1 1/2" x 3" (3.8 x 7.6 cm) 4 – Fingertip, 1 1/4" x 2" (4.4 x 5.1 cm) 18 – 7/8" x 3" (2.2 x 7.6 cm)	30/bx, 24 bx/cs
CUR14926RB	3 3/4" x 2" (9.5 x 5.1 cm)	10/bx, 24 bx/cs
CUR14925RB	1" x 3" (2.5 x 7.6 cm)	20/bx, 24 bx/cs
CUR01101RB	3/4" x 4 3/4" (1.9 x 12.1 cm), Extreme Length	20/bx, 24 bx/cs



CUR14924RB

CURAD® Sheer Bandages

Designed to blend with your skin.

- » Nearly invisible protection while you heal
- » Skin-friendly adhesive won't leave a sticky, itchy residue on the skin
- » Maximum comfort as you heal
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR02279RB	3/4" x 2 7/8" (1.9 x 7.3 cm)	80/bx, 24 bx/cs
CUR02277RB	2" x 3 3/4" (5.0 x 9.5 cm)	10/bx, 24 bx/cs
CUR01726RB	3" x 4" (7.6 x 10.1 cm)	10/bx, 24 bx/cs
CUR02276RB	Sheer, Assorted 10 – 1" x 3" (2.5 x 7.6 cm) 20 – 3/4" x 3" (1.9 x 7.6 cm) 30 – 5/8" x 2 1/4" (1.6 x 5.7 cm)	80/bx, 24 bx/cs
CUR45243RB	Sheer, Assorted 20 – Spot 7/8" (2.2 cm) 20 – 3/4" x 2 1/8" (1.9 x 5.4 cm) 30 – 3/4" x 2 7/8" (1.9 x 7.3 cm) 10 – 1 1/8" x 2 7/8" (2.9 x 7.3 cm)	80/bx, 24 bx/cs



CUR02279RB

CURAD® Pressure Adhesive Bandages

Once in contact with fluid, the absorbent cellulose island pad expands to over 7 times its normal height.

- » Helps provide gentle pressure to stop bleeding
- » Waterproof PE backing prevents strike-through
- » Promotes comfort for dialysis, chemotherapy and cardiac cath patients
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
NON85100	Large, 1" x 2¾" (2.5 x 7 cm)	100/bx, 20 bx/cs, 2000/cs



CURAD® Detectable Adhesive Bandages

Highly visible color meets the needs of the food service industry.

- » Innovative, long-lasting adhesive is gentle on the skin
- » Contains detectable metal to easily trigger metal detectors used by food packagers
- » Highly identifiable blue color stands out against a variety of foods
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
Woven Blue Cloth, Sterile		
NON25650BL	¾" x 3" (1.9 x 7.6 cm)	100/bx, 12 bx/cs
NON25660BL	1" x 3" (2.5 x 7.6 cm)	100/bx, 12 bx/cs
NON25513BL	Fingertip, 1¾" x 2" (4.4 x 5.1 cm)	100/bx, 12 bx/cs
NON25510BL	Knuckle, 1½" x 3" (3.8 x 7.6 cm)	100/bx, 12 bx/cs



Blue Plastic, Sterile

NON25500BL	¾" x 3" (1.9 x 7.6 cm)	100/bx, 12 bx/cs
NON25600BL	1" x 3" (2.5 x 7.6 cm)	100/bx, 12 bx/cs

CURAD® Neon Adhesive Bandages

Assortment of five neon colors.

- » Each box includes 10 plastic strips each of pink, blue, green, purple and yellow
- » Non-stick island pad protects granulating tissue
- » Adhesive provides long lasting protection.
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
NON256131	¾" x 3" (1.9 x 7.6 cm)	50/bx, 24 bx/cs, 1200/cs



NON256131


RETAIL FRIENDLY

CURAD® Owl Bandages

Provides maximum comfort for kids.

- » These fun, kid friendly bandages feature colorful owl designs
- » Non-stick pad with long lasting adhesion
- » Assorted sizes
- » Sterile; Not made with natural rubber latex

Item No.	Description	Pkg.
CUR00004RB	Owl Bandages, Assorted 12 – 1" x 1" (2.5 x 2.5 cm) 8 – 3/4" x 3" (1.9 x 7.6 cm)	20/bx, 24bx/cs



CUR00004RB

CURAD® Busy Bugs Bandages

Ideal for everyday use for kids.

- » Curad Busy Bugs bandages are colorful and kid-friendly
- » Assorted sizes
- » Sterile; Not made with natural rubber latex

Item No.	Description	Pkg.
CUR00001RB	Busy Bugs, Assorted 12 – 1" x 1" (2.5 x 2.5 cm) 8 – 3/4" x 3" (1.9 x 7.6 cm)	20/bx, 24 bx/cs

Sesame Street and associated characters, trademarks and design elements are owned and licensed by Sesame Workshop. ©2008 Sesame Workshop. All Rights Reserved.



CUR00001RB

CURAD® Medtoons® Adhesive Bandages

Children will adore the smiling bumble bees, soccer-playing dogs and snorkeling fish.

- » Each box includes 50 strips
- » Constructed from PVC-free vinyl
- » Perforated for increased moisture vapor transmission
- » Acrylic adhesive is skin friendly and provides long-lasting protection
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
NON256130	¾" x 3" (1.9 x 7.6 cm)	50/bx, 24 bx/cs, 1200/cs



NON256130

CURAD® Cupcake Cover Bandages

The perfect bandage for children of all ages.

- » Curad Cupcake Cover bandages are colorful and kid-friendly
- » Assorted sizes
- » Sterile; Not made with natural rubber latex

Item No.	Description	Pkg.
CUR00002RB	Cupcakes, Assorted 12 – 1" x 1" (2.5 x 2.5 cm) 8 – 3/4" x 3" (1.9 x 7.6 cm)	20/bx, 24bx/cs



CUR00002RB

CURAD® Camp Camo Adhesive Bandages

Green, pink and blue camo patterns.

- » Woven fabric stretches with movement
- » Strong adhesive offers durable protection
- » Absorbent non-stick island pad is painless to remove
- » Ideal for knees and elbows
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR45701RB	Green, 3/4" x 3" (1.9 x 7.6 cm)	25/bx, 24 bx/cs
CUR45702RB	13 – Pink, 3/4" x 3" (1.9 x 7.6 cm) 12 – Blue, 3/4" x 3" (1.9 x 7.6 cm)	25/bx, 24 bx/cs



CUR45701RB

CUR45702RB

CURAD® Blister Heal™

Great for shallow wounds or preventing blisters.

- » Hydrocolloid
- » Smooth backing
- » Tapered edge
- » Longer wear time
- » Sterile; not made with natural rubber latex

Item No.	Description	Pkg.
CUR76122RB	3 – 2 3/4" x 1.7" (7 x 4.2 cm) 3 – 2.4" x 1.5" (6 x 3.8 cm)	6/bx, 24 bx/cs



CUR76122RB

CURAD® Adhesive Moleskin

Tan flannel with a napped cotton backing.

- » Zinc oxide-based adhesive
- » Cushions the skin to help manage blisters, corns, calluses and chafing
- » Cut to desired size

Item No.	Description	Pkg.
NON5148	9" x 4-yd. (22.9 cm x 3.7 m) Not made with natural rubber latex	1 roll/tube, 12 tubes/cs
NON5147	9" x 4-yd. (22.9 cm x 3.7 m) Contains natural rubber latex	1 roll/tube, 12 tubes/cs

NON5148



CURAD® Medi-Strips® Reinforced Wound Closure Strips

Compare to 3M™ Steri-Strips®*

Spun-bond fiber backing with reinforced filaments for high tensile strength.

- » Offers excellent adhesion while being gentle to the skin
- » Use in addition to or in place of sutures
- » Color-coded for quick and easy size identification
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON250318	1/8" x 3" (0.3 x 7.6 cm)	5/pk, 50 pk/bx, 4 bx/cs
NON250114	1/4" x 1 1/2" (0.6 x 3.8 cm)	6/pk, 50 pk/bx, 4 bx/cs
NON250412	1/2" x 4" (1.3 x 10.2 cm)	6/pk, 50 pk/bx, 4 bx/cs
NON250314	3/4" x 3" (0.6 x 7.6 cm)	3/pk, 50 pk/bx, 4 bx/cs
NON250414	1/4" x 4" (0.6 x 10.2 cm)	5 strips/card, 2 cards/pk, 50 pk/bx, 4 bx/cs
NON250501	1" x 5" (2.5 x 12.7 cm)	4/pk, 25 pk/bx, 4 bx/cs

NON250318



*This product is not manufactured or distributed by 3M Healthcare, makers of Steri-Strips®.

CURAD® Nasal Tube Fastener

One-piece, flexible non-woven fastening strip.

- » Hypoallergenic adhesive minimizes skin irritation
- » Remains securely in place until removed

Item No.	Description	Pkg.
CURNASTUB50	Adult, One-piece, One Size Fits Most	10/pk, 5 pk/bx, 50/bx
CURNSO4430	Nasal Strip	18/cs

CURNASTUB50



CURAD® Mediplast®

Effective for removing corns, calluses and warts.

- » Adhesive cut-to-fit pads
- » Soft, convenient and easy to use
- » Active ingredient Salicylic Acid (40%), highly effective and widely recommended by physicians and pharmacists
- » Adhesive contains natural rubber latex

Item No.	Description	Pkg.
CUR01496	2" x 3" pads (5.1 x 7.6 cm)	25/bx, 6 bx/cs

CUR01496





MDS066004

Cotton Undercast Padding

100% cotton.

- » Provides a soft environment under a hard cast
- » Low-linting with more loft
- » Provides greater patient comfort
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile (Singles)		
DYNJ066002	2" x 4-yd. (5.1 cm x 3.7 m)	20 rolls/cs
DYNJ066003	3" x 4-yd. (7.6 cm x 3.7 m)	20 rolls/cs
DYNJ066004	4" x 4-yd. (10.2 cm x 3.7 m)	20 rolls/cs
DYNJ066006	6" x 4-yd. (15.2 cm x 3.7 m)	20 rolls/cs
Non-Sterile		
MDS066002	2" x 4-yd. (5.1 cm x 3.7 m)	24/bg, 3 bg/cs, 72/cs
MDS066003	3" x 4-yd. (7.6 cm x 3.7 m)	12/bg, 6 bg/cs, 72/cs
MDS066004	4" x 4-yd. (10.2 cm x 3.7 m)	12/bg, 6 bg/cs, 72/cs
MDS066006	6" x 4-yd. (15.2 cm x 3.7 m)	6/bg, 6 bg/cs, 36/cs



MDS067002

Synthetic Undercast Padding

- » 100% polyester needle punched construction
- » Ideal for use under all forms of synthetic hard casting material
- » Absorbs excess moisture and dries quickly
- » Provides a comfortable fit
- » Not made with natural rubber latex

Item No.	Description	Pkg.
MDS067002	2" x 4-yd. (5.1 cm x 3.7 m)	12/bg, 6 bg/cs, 72/cs
MDS067003	3" x 4-yd. (7.6 cm x 3.7 m)	12/bg, 6 bg/cs, 72/cs
MDS067004	4" x 4-yd. (10.2 cm x 3.7 m)	12/bg, 6 bg/cs, 72/cs
MDS067006	6" x 4-yd. (15.2 cm x 3.7 m)	6/bg, 6 bg/cs, 36/cs
MDS067008	8" x 4-yd. (20.3 cm x 3.7 m)	6/bg, 6 bg/cs, 36/cs



NONFAK100

Bloodborne Pathogen Kit

30 pieces, personal protection use

Guards caregivers of ill or injured patients and protects during biohazard clean-up. Sturdy plastic case is easily portable or wall-mountable.

Contents:

- (1) Fluid control solidifier, 21 gm pack
- (1) Biohazard scoop
- (1) Disposable gown w/full sleeves
- (2) Disposable shoe covers
- (1) Disposable bonnet
- (2) 24" x 24" Biohazard bags, 10 gallon capacity
- (2) 8" x 12" Clear plastic bags
- (2) Twist ties
- (3) Disposable clean-up towels
- (2) Exam quality vinyl gloves, 1 pair
- (1) Eye shield with ear loop mask
- (3) Antiseptic cleansing wipes (sting free)
- (2) 5" x 8" Germicidal (kills germs) wipe, covers 7.5 sq. ft.
- (1) CPR Pack: 1 CPR one-way valve faceshield, latex free, 2 exam quality vinyl gloves, 3 antiseptic cleansing wipes (sting free)
- (1) Bodily fluid pick-up guide

Kit dimensions: 9" x 8³/₈" x 2¹/₂" (22.9 x 21.3 x 6.4 cm)

Item No.	Description	Pkg.
NONFAK100	Bloodborne Pathogen Kit	1 kit



NONFAK200

General First Aid Kit

106 pieces, 25 person use

Ideal for worksites with up to 25 employees. It carries 17 critical products including a 4 oz. eye wash. The partitioned, water-resistant plastic case is wall mountable.

Contents:

- (3) 1³/₄" x 3" Large fingertip fabric bandages
- (2) 1¹/₂" x 3" Knuckle fabric bandages
- (50) 3/4" x 3" Adhesive plastic bandages
- (1) 5" x 9" Trauma pad
- (6) 2" x 2" Gauze dressing pads, (3) 2-pks
- (1) 2" Conforming gauze bandage roll
- (1) Sterile eye pad
- (1) 1/2" x 5-yd. First aid tape roll
- (1) Triangular sling/bandage
- (12) Antiseptic cleansing wipes (sting free)
- (6) Triple antibiotic ointment packs
- (12) Aspirin tablets, (6) 2-pks
- (1) 4" x 5" Instant cold compress
- (2) Exam quality vinyl gloves, 1 pair
- (1) 4¹/₂" Scissors, nickel plated
- (1) 4" Tweezers, plastic
- (1) Eye wash, 4-oz.
- (3) Insect sting relief pads
- (1) First aid guide

Kit dimensions: 9" x 8³/₈" x 2¹/₂" (22.9 x 21.3 x 6.4 cm)

Item No.	Description	Pkg.
NONFAK200	General First Aid Kit	1 kit



NONFAK300

General Workplace Cabinet

206 pieces, 50 person use

The unique, simple cabinet design is an effective system for recording and restocking, eliminating disorganization and missing supplies. Clearly labeled compartments create a cabinet easy to use. With expansion pocket, you can customize your cabinet to meet your specific needs.

Contents:

- | | |
|---|---|
| (10) 1½" x 1½", Patch plastic bandage | (20) Alcohol cleansing pad |
| (25) 1" x 3", Fabric bandage | (20) Antiseptic cleansing wipe (sting free) |
| (25) ¾" x 3", Plastic bandage | (10) Castile soap towelette |
| (10) Fingertip fabric bandage | (10) First aid/burn cream |
| (10) Knuckle fabric bandage | (10) Antibiotic ointment |
| (1) Elastic bandage wrap | (1) ½" x 10-yd., First aid tape |
| (10) Butterfly wound closure | (1) Scissors |
| (1) 40" x 40" x 56", Triangular sling/bandage | (10) 2", Moleskin (blister prevention) |
| (10) 2" x 2", Gauze dressing pad | (1) First aid guide |
| (10) 3" x 3", Gauze dressing pad | (2) Exam quality vinyl gloves |
| (1) 2" Conforming gauze roll | (1) 3" Stainless steel tweezers |
| (1) 3" Conforming gauze roll | (1) CPR face shield and vinyl gloves |
| (1) 5" x 9", Trauma pad | |

Kit dimensions: 13" x 14" x 4⅞" (33 x 35.6 x 10.5 cm)

Item No.	Description	Pkg.
NONFAK300	General Workplace Cabinet	1 kit

SECUREMENT

Hold dressings in place

Co-Flex® LF2

Compare to 3M Coban®*

A unique quick-stick self-adhering bandage.

- » Soft foam layer
- » Performs better than latex
- » Features patented hand tear technology
- » Will not constrict
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile, Stretched, Tan		
DYNJ089001	1" x 5-yd. (2.5 cm x 4.6 m)	20 rolls/cs
DYNJ089002	2" x 5-yd. (5.1 cm x 4.6m)	20 rolls/cs
DYNJ089003	3" x 5-yd. (7.6 cm x 4.6 m)	20 rolls/cs
DYNJ089004	4" x 5-yd. (10.2 cm x 4.6 m)	20 rolls/cs
DYNJ089006	6" x 5-yd. (15.2 cm x 4.6 m)	20 rolls/cs

Non-Sterile, Stretched, Tan

MDS089001	1" x 5-yd. 2/pkg (2.5 cm x 4.6 m)	30 rolls/cs
MDS089002	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
MDS089003	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
MDS089004	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
MDS089006	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs

Co-Flex® LF2 Colorpack, Non-Sterile, Stretched

Contains rolls of purple, light blue, blue, red, neon green and pink

MDS089001CP	1" x 5-yd. (2.5 cm x 4.6 m)	30 rolls/cs
MDS089002CP	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
MDS089003CP	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
MDS089004CP	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
MDS089006CP	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs

* This product is not manufactured or distributed by 3M Healthcare, makers of Coban®.

Co-Flex® NL

Kid's pack contains rolls of:

- » Neon Pink/Purple Hearts
- » Light Blue/Dark Blue Trucks and Cars
- » Neon Green/Dark Green Dinosaurs
- » Yellow/Black Bumblebee

Item No.	Description	Pkg.
MDS088002KP	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
MDS088003KP	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs



MDS089004

MDS089002



MDS089002CP



MDS088002KP

Co-Flex® Med

Self-adherent, hand-tear compression wrap.
Contains natural rubber latex.

Item No.	Description	Pkg.
Co-Flex® Med, Sterile, Stretched		
DYNJ086002	2" x 5-yd. (5.1 cm x 4.6 m)	20 rolls/cs
DYNJ086003	3" x 5-yd. (7.6 cm x 4.6 m)	20 rolls/cs
DYNJ086004	4" x 5-yd. (10.2 cm x 4.6 m)	20 rolls/cs
DYNJ086006	6" x 5-yd. (15.2 cm x 4.6 m)	20 rolls/cs
Co-Flex® Med, Non-Sterile, Stretched		
MDS086001	1" x 5-yd., 2/pk (2.5 cm x 4.6 m)	30 rolls/cs
MDS086002	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
MDS086003	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
MDS086004	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
MDS086006	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs
Co-Flex® Med, Color Pack, Non-Sterile, Stretched		
Contains rolls of purple, light blue, blue, red, neon green and pink		
MDS086001CP	1" x 5-yd. (2.5 cm x 4.6 m)	30 rolls/cs
MDS086002CP	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
MDS086003CP	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
MDS086004CP	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
MDS086015CP	1.5" x 5-yd. (3.8 cm x 4.6 m)	48 rolls/cs



Co-Flex® AFD

High absorption and retention of exudate.

- » Keeps wound site moist and clean
- » All-in-one dressing/cohesive bandage
- » Soft and comfortable CoFlex LF2 base
- » Waterproof film layer over foam pad
- » Expands to conform to wound cavity
- » Provides controlled compression and easy application
- » Non-adherent to wound bed
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Co-Flex® AFD, Sterile		
DYNJ089202	2" x 2.5-yd. (5.1 cm x 2.3 m)	16 rolls/cs
DYNJ089204	4" x 2.5-yd. (10.2 cm x 2.3 m)	8 rolls/cs
DYNJ089207	4" x 7.5-yd. (10.2 cm x 6.9 m)	6 rolls/cs
Co-Flex® AFD, Non-Sterile		
MDS089202	2" x 2.5-yd. (5.1 cm x 2.3 m)	16 rolls/cs
MDS089204	4" x 2.5-yd. (10.2 cm x 2.3 m)	8 rolls/cs
MDS089207	4" x 7.5-yd. (10.2 cm x 6.9 m)	6 rolls/cs



MDS089204

MEDLINE ESSENTIALS

Self-Adherent Wrap

Sticks to itself but not to skin, hair or dressings. No pins or clips required.

- » Secondary dressing functions like a tape to retain dressings without tension around digits or limbs
- » Promotes pain-free removal that won't damage fragile skin
- » Comfortable, porous and lightweight construction
- » Ideal for wrapping splints, gauze pads, I.V. and dressing securement
- » Serves as a compression wrap after blood draws
- » Each bandage is individually wrapped

Item No.	Description	Pkg.
----------	-------------	------

Tan, Non-sterile; not made with natural rubber latex

PRM088001	1" x 5-yd. (2.5 cm x 4.6 m)	30 rolls/cs
PRM088002	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
PRM088003	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
PRM088004	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
PRM088006	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs

Tan, Non-sterile; contains natural rubber latex

PRM086001	1" x 5-yd. (2.5 cm x 4.6 m)	30 rolls/cs
PRM086002	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
PRM086003	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
PRM086004	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
PRM086006	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs

Tan, Sterile; contains natural rubber latex

DYNJPRM086003	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
DYNJPRM086004	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
DYNJPRM086006	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs

Tan, Sterile; not made with natural rubber latex

DYNJPRM088003	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs
DYNJPRM088004	4" x 5-yd. (10.2 cm x 4.6 m)	18 rolls/cs
DYNJPRM088006	6" x 5-yd. (15.2 cm x 4.6 m)	12 rolls/cs

Color Pack, Non-sterile; not made with natural rubber latex

PRM088001CP	1" x 5-yd. (2.5 cm x 4.6 m)	30 rolls/cs
PRM088002CP	2" x 5-yd. (5.1 cm x 4.6 m)	36 rolls/cs
PRM088003CP	3" x 5-yd. (7.6 cm x 4.6 m)	24 rolls/cs



PRM088001



PRM088002CP



MDS087004LF

Matrix® Elastic Bandage

Our finest elastic bandage. Cotton and polyester blend.

- » Ideal for medium compression
- » Excellent elasticity and recoverability
- » Hook and loop closure at both ends
- » Color-coded for easy size identification
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile (Singles), Stretched		
DYNJ05152LF	2" x 5-yd. (5.1 cm x 4.6 m)	20/cs
DYNJ05153LF	3" x 5-yd. (7.6 cm x 4.6 m)	20/cs
DYNJ05154LF	4" x 5-yd. (10.2 cm x 4.6 m)	20/cs
DYNJ05155LF	4" x 10-yd. (10.2 cm x 9.1 m)	20/cs
DYNJ05156LF	6" x 5-yd. (15.2 cm x 4.6 m)	20/cs
DYNJ05157LF	6" x 10-yd. (15.2 cm x 9.1 m)	20/cs
DYNJ05158LF	6" x 15-yd. (15.2 cm x 13.7 m)	20/cs
Non-Sterile, Stretched		
MDS087002LF	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS087003LF	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS087004LF	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS087006LF	6" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS087104LF	4" x 10-yd. (10.2 cm x 9.1 m)	4/bx, 5 bx/cs, 20/cs
MDS087106LF	6" x 10-yd. (15.2 cm x 9.1 m)	4/bx, 5 bx/cs, 20/cs
MDS087156LF	6" x 15-yd. (15.2 cm x 13.7 m)	4/bx, 5 bx/cs, 20/cs



MDS077003

Swift-Wrap® Elastic Bandage

A high quality elastic bandage.

Cotton and polyester blend.

- » Convenient hook and loop closure
- » No need for clips or tape
- » Ideal for high compression
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Sterile (Singles), Stretched		
DYNJ05144	2" x 5-yd. (5.1 cm x 4.6 m)	20/cs
DYNJ05145	3" x 5-yd. (7.6 cm x 4.6 m)	20/cs
DYNJ05146	4" x 5-yd. (10.2 cm x 4.6 m)	20/cs
DYNJ05147	6" x 5-yd. (15.2 cm x 4.6 m)	20/cs
Non-Sterile, Stretched		
MDS077002	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS077003	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS077004	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS077006	6" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS055002Z	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 20/cs
MDS055003Z	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 20/cs
MDS055004Z	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 20/cs
MDS055006Z	5" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 20/cs



MDS046004



MDS057004

Soft-Wrap® Elastic Bandage

Knitted and breathable. Cotton and polyester blend.

- » Recommended for support and low compression
- » An ideal secondary bandage
- » Clips included
- » Not made with natural rubber latex

Item No.	Description	Pkg.
----------	-------------	------

Sterile (Singles), Stretched

DYNJ05125	2" x 5-yd. (5.1 cm x 4.6 m)	20/cs
DYNJ05126	3" x 5-yd. (7.6 cm x 4.6 m)	20/cs
DYNJ05131	4" x 5-yd. (10.2 cm x 4.6 m)	20/cs
DYNJ05132	6" x 5-yd. (15.2 cm x 4.6 m)	20/cs

Non-Sterile, Stretched

MDS046002	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS046003	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS046004	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS046006	6" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs

Sure-Wrap® Elastic Bandage

Our most cost-effective elastic bandage.

Cotton and polyester blend.

- » Delivers consistent pressure
- » Appropriate for applications requiring higher compression
- » Clips included
- » Not made with natural rubber latex

Item No.	Description	Pkg.
----------	-------------	------

Non-Sterile, Stretched, White

MDS057002	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS057003	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS057004	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS057006	6" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs

Non-Sterile, Stretched, Beige

MDS055002	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS055003	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS055004	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS055006	6" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 5 bx/cs, 50/cs
MDS077002Z	2" x 5-yd. (5.1 cm x 4.6 m)	10/bx, 20/cs
MDS077003Z	3" x 5-yd. (7.6 cm x 4.6 m)	10/bx, 20/cs
MDS077004Z	4" x 5-yd. (10.2 cm x 4.6 m)	10/bx, 20/cs
MDS077006Z	5" x 5-yd. (15.2 cm x 4.6 m)	10/bx, 20/cs

CURAD® Elastic Net

Compare to Surgilast®*

Rayon/polyester knit. Favored by wound care nurses everywhere. Tubular retainer net holds dressings firmly in place without compression. It is an essential product for fragile, at-risk skin, or for patients with tape allergies.

- » Easy to cut
- » Reduces dressing time
- » Not made with neoprene or natural rubber latex



NONNET04

Item No.	Description	Size	Stretched Width	Relaxed Width	Pkg.
Elastic Net, 25-yd. Stretched Length, Domestic					
NONNET012	Small Fingers, Toes, Wrists, Circumcision	½	4" (10.2 cm)	⅝" (1.6 cm)	1/bx
NONNET01	Large Fingers, Toes, Wrists	1	6" (15.2 cm)	0.6" (1.5 cm)	1/bx
NONNET02	Small Hand, Arm, Leg, Foot	2	7⅞" (19.5 cm)	0.8" (2 cm)	1/bx
NONNET03	Medium Hand, Arm, Leg, Foot, Multiple Fingers	3	9½" (23.7 cm)	1" (2.5 cm)	1/bx
NONNET04	Large Hand, Arm, Leg, Foot	4	11½" (29.6 cm)	1.2" (3 cm)	1/bx
NONNET05	Average Arm, Small Leg/Knee	5	12⅞" (31.4 cm)	1.5" (3.8 cm)	1/bx
NONNET06	Large Leg, Amputation, Arm Board, Small Head	6	17" (61 cm)	1.9" (4.8 cm)	1/bx
NONNET07	Average Head, Face, Ears, Neck	7	28⅜" (72.8 cm)	2.3" (5.8 cm)	1/bx
NONNET08	Small Chest, Back, Perineum, Axilla, Shoulder, Large Head	8	31⅞" (79.6 cm)	2.8" (7.1 cm)	1/bx
NONNET09	Medium Chest, Back, Perineum, Axilla, Shoulder	9	35" (88.9 cm)	3.2" (8.1 cm)	1/bx
NONNET10	Large Chest, Back, Perineum, Axilla, Shoulder	10	41" (93.1 cm)	3.6" (9.1 cm)	1/bx
NONNET11	X-Large Chest, Back, Perineum, Axilla, Shoulder	11	46" (101.6 cm)	4.3" (10.9 cm)	1/bx
NONNET12	XX-Large Chest, Back, Perineum, Axilla, Shoulder	12	56" (142.2 cm)	4.7" (11.9 cm)	1/bx
NONNET22	Bariatric/Custom Size	22	80" (203.2 cm)	8" (20.3 cm)	1/bx


*This product is not manufactured by Western Medical, Inc., makers of Surgilast®.

TIPS!

If you lose your applicator, simply take a 20cc syringe (without the needle), remove the plunger and place the netting over the syringe. Place the digit inside the syringe and slowly move netting on to the digit.

MEDLINE ESSENTIALS

Item No.	Description	Size	Stretched Width	Relaxed Width	Pkg.
Elastic Net, 25-yd. Length, Imported					
NONSTNDNET01	Large Fingers, Toes, Wrists	1	6" (15.2 cm)	0.6" (1.5 cm)	1/bx
NONSTNDNET02	Small Hand, Arm, Leg, Foot	2	7" (17.8 cm)	0.8" (2 cm)	1/bx
NONSTNDNET03	Medium Hand, Arm, Leg, Foot, Multiple Fingers	3	9" (22.9 cm)	1" (2.5 cm)	1/bx
NONSTNDNET04	Large Hand, Arm, Leg, Foot	4	11" (27.9 cm)	1.2" (3 cm)	1/bx
NONSTNDNET05	Average Arm, Small Leg/Knee	5	12" (30.5 cm)	1.5" (3.8 cm)	1/bx
NONSTNDNET06	Large Leg, Amputation, Arm Board, Small Head	6	17" (43.2 cm)	1.9" (4.8 cm)	1/bx
NONSTNDNET07	Average Head, Face, Ears, Neck	7	28" (71.1 cm)	2.3" (5.8 cm)	1/bx
NONSTNDNET08	Small Chest, Back, Perineum, Axilla, Shoulder, Large Head	8	31" (78.7 cm)	2.8" (7.1 cm)	1/bx
NONSTNDNET09	Medium Chest, Back, Perineum, Axilla, Shoulder	9	35" (88.9 cm)	3.2" (8.1 cm)	1/bx
NONSTNDNET10	Large Chest, Back, Perineum, Axilla, Shoulder	10	41" (104.1 cm)	3.6" (9.1 cm)	1/bx
NONSTNDNET11	X-Large Chest, Back, Perineum, Axilla, Shoulder	11	46" (116.8 cm)	4.3" (10.9 cm)	1/bx
NONSTNDNET12	XX-Large Chest, Back, Perineum, Axilla, Shoulder	12	56" (142.2 cm)	4.7" (11.9 cm)	1/bx
NONSTNDNET22	Bariatric/Custom Size	22	80" (203.2 cm)	8" (20.3 cm)	1/bx



CURAD® Pre-Cut Elastic Net

Individually wrapped.

- » Pre-cut sizes
- » Easy to use
- » Ideal for home health settings

Item No.	Description	Size	Pkg.
Individual 24" Lengths (Relaxed)			
NONNET124P	Stretched Width: 6" (15.2 cm)	1	25/cs
NONNET224P	Stretched Width: 7 ³ / ₈ " (19.5 cm)	2	25/cs
NONNET324P	Stretched Width: 9 ¹ / ₂ " (23.7 cm)	3	25/cs
NONNET424P	Stretched Width: 11 ³ / ₈ " (29.6 cm)	4	25/cs
NONNET524P	Stretched Width: 12 ³ / ₈ " (31.4 cm)	5	25/cs
NONNET624P	Stretched Width: 17" (43.2 cm)	6	25/cs
NONNET724P	Stretched Width: 28 ³ / ₈ " (72.8 cm)	7	25/cs
NONNET824P	Stretched Width: 31 ¹ / ₂ " (79.6 cm)	8	25/cs
NONNET924P	Stretched Width: 35" (88.9 cm)	9	25/cs
NONNET1024P	Stretched Width: 41" (104.1 cm)	10	25/cs
NONNET1124P	Stretched Width: 46" (116.8 cm)	11	25/cs
NONNET1224P	Stretched Width: 56" (142.2 cm)	12	25/cs
NONNET2224P	Stretched Width: 80" (203.2 cm)	22	25/cs
Individual 36" Lengths (Relaxed)			
NONNET636P	Stretched Width: 17" (43.2 cm)	6	50/cs
NONNET736P	Stretched Width: 28 ³ / ₈ " (72.8 cm)	7	50/cs



NONNET924P

Paper Tape

A gentle and breathable tape for sensitive skin.

- » Recommended for pediatric and geriatric patients
- » Flexible and conformable
- » Adheres well to damp skin
- » Single use rolls reduce waste and cross contamination
- » Not made with natural rubber latex

PREMIUM

CURAD®

Compare to Micropore®*

Item No.	Description	Pkg.
NON270012	½" x 10-yd. (1.3 cm x 9.1 m)	24 rolls/bx, 10 bx/cs, 240 rolls/cs
NON270001	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 10 bx/cs, 120 rolls/cs
NON270002	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 10 bx/cs, 60 rolls/cs
NON270003	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 10 bx/cs, 40 rolls/cs
NON270001S	1" x 1.5-yd. (2.5 cm x 1.4 m)	100 rolls/bx, 5 bx/cs, 500 rolls/cs

* This product is not manufactured or distributed by 3M Healthcare, makers of Micropore®.



NON270001

MEDLINE ESSENTIALS

Compare to DUKAL Paper Tape**

- » Specification curated to meet lighter adhesive needs

Item No.	Description	Pkg.
PRM260001	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260001H	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx
PRM260002	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260002H	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx
PRM260003	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs
PRM260003H	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx

** This product is not manufactured or distributed by Dukal Corporation.



PRM260001

CURAD® Sensitive Skin Paper Tape

Revolutionary gentle adhesion.

- » Adhesive is present only on every alternate wave
- » Unique wave pattern provides plenty of surface area for the tape to stick, yet the smaller amount of adhesive is gentler to the skin when the tape comes off
- » Gaps in the adhesive also allow for breathability
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON260002E	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
NON260003E	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs



NON260002E

Cloth "Silk" Tape

Strong woven construction. Made from silk-like cloth.

- » Ideal for affixing dressings.
- » Water repellent
- » Easy to tear
- » Not made with natural rubber latex

PREMIUM

CURAD®

Compare to Durapore®*

Item No.	Description	Pkg.
NON270112	½" x 10-yd. (1.3 cm x 9.1 m)	24 rolls/bx, 10 bx/cs, 240 rolls/cs
NON270101	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 10 bx/cs, 120 rolls/cs
NON270102	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 10 bx/cs, 60 rolls/cs
NON270103	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 10 bx/cs, 40 rolls/cs

* This product is not manufactured or distributed by 3M Healthcare, makers of Durapore®.



NON270101

MEDLINE ESSENTIALS

Compare to DUKAL Cloth**

- » Specification curated to meet lighter adhesive needs

Item No.	Description	Pkg.
PRM260101	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260101H	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx
PRM260102	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260102H	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx
PRM260103	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs

** This product is not manufactured or distributed by Dukal Corporation.



PRM260101

CURAD® Waterproof Tape

Easy to write on. Cotton backing with polyethylene coating.

- » Waterproof and oil-resistant
- » Provides high tensile strength
- » Won't permit moisture to penetrate the dressing
- » Strong, durable adhesive holds securely
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON260501	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
NON260502	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
NON260503	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs
CUR47440	½" x 5-yd. (1.3 cm x 4.6 m)	24 rolls/cs
CUR47441	½" x 10-yd. (1.3 cm x 9.1 m)	24 rolls/cs



NON260501

Transparent Tape

A breathable perforated plastic tape.

- » Permits skin examination without tape removal
- » “Water repellent” and provides a clean and easy bilateral tear for custom widths
- » Provides a clean and easy bilateral tear for custom widths
- » Excellent for holding facial dressings, IV sites, tubing retention
- » Not made with natural rubber latex

PREMIUM

CURAD®

Compare to Transpore®*

Item No.	Description	Pkg.
NON270212	½" x 10-yd. (1.3 cm x 9.1 m)	24 rolls/bx, 10 bx/cs, 240 rolls/cs
NON270201	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 10 bx/cs, 120 rolls/cs
NON270202	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 10 bx/cs, 60 rolls/cs
NON270203	3" x 10-yd. (7.6 cm x 9.1 m)	4 rolls/bx, 10 bx/cs, 40 rolls/cs
NON270201S	1" x 1.5-yd. (2.5 cm x 1.4 m)	100 rolls/bx, 5 bx/cs, 500 rolls/cs
NON270202S	2" x 1.5-yd. (5.1 cm x 1.4 m)	50 rolls/bx, 5 bx/cs, 250 rolls/cs

* This product is not manufactured or distributed by 3M Healthcare, makers of Transpore®.



NON270201

MEDLINE ESSENTIALS

Compare to DUKAL Transparent Tape*

- » Specification curated to meet lighter adhesive needs

Item No.	Description	Pkg.
PRM260201	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
PRM260201H	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx
PRM260202	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
PRM260202H	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx

* This product is not manufactured or distributed by Dukal Corporation.



PRM260201

CURAD® Ortho-Porous Sports Tape

Strong woven cotton backing provides high tensile strength and support to highly active patients.

- » Powerful rubber-based adhesive for long lasting adhesion with minimal residue
- » Highly porous construction allows the skin to breathe
- » Permits moisture vapor transmission
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON260301	1" x 10-yd. (2.5 cm x 9.1 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
NON260315	1½" x 10-yd. (3.8 cm x 9.1 m)	8 rolls/bx, 12 bx/cs, 96 rolls/cs
NON260302	2" x 10-yd. (5.1 cm x 9.1 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
NON260303	3" x 10-yd. (7.5 cm x 9.1 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs



NON260301

CURAD® Elastic Adhesive Tape

100% cotton elastic cloth tape backing.

- » Provides high tensile strength
- » Strong rubber-based adhesive permits body movement without loosening of tape
- » Porous construction maintains skin integrity
- » Permits moisture vapor transmission
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON260401	1" x 5-yd. (2.5 cm x 4.6 m)	12 rolls/bx, 12 bx/cs, 144 rolls/cs
NON260402	2" x 5-yd. (5.1 cm x 4.6 m)	6 rolls/bx, 12 bx/cs, 72 rolls/cs
NON260403	3" x 5-yd. (7.6 cm x 4.6 m)	4 rolls/bx, 12 bx/cs, 48 rolls/cs
NON260404	4" x 5-yd. (10.2 cm x 4.6 m)	6 rolls/bx, 6 bx/cs, 36 rolls/cs



NON260401

CURAD® Elastic Foam Tape

Highly conformable, multi-directional stretch.

- » Promotes patient comfort, accommodates swelling
- » Moderate adhesion secures dressings over challenging areas
- » Provides moderate compression
- » Smooth, cushioned elastic foam backing protects wound site and dressings
- » Water resistant for greater durability
- » Adhesive and backing leave minimal skin residue
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON260601	1" x 5.5-yd. (2.5 cm x 5 m)	12 rolls/bx, 6 bx/cs, 72 rolls/cs
NON260602	2" x 5.5-yd. (5.1 cm x 5 m)	6 rolls/bx, 6 bx/cs, 36 rolls/cs
NON260603	3" x 5.5-yd. (7.6 cm x 5 m)	4 rolls/bx, 6 bx/cs, 24 rolls/cs
NON260604	4" x 5.5-yd. (10.2 cm x 5 m)	3 rolls/bx, 6 bx/cs, 18 rolls/cs



NON260601

Montgomery Straps

Ideal for frequent dressing changes.

- » No need to remove and reapply tape
- » Breathable, strong, comfortable
- » Helps reduce the risk of patient skin trauma associated with frequent taping
- » Cut to size
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON5129	7½" x 11½" (18.4 x 28.3 cm)	24/bx (12 pairs), 6 bx/cs, 144 ea/cs



NON5129

Adhesive Surgical Dressing

Helps reduce nursing time by combining a sterile sponge, a non-adhesive pad, and a cloth "silk" tape into one dressing.

- » Non-adherent contact layer permits free flow of drainage
- » Won't stick to granulating tissue
- » Ideal for light to medium draining wounds
- » Tape backing holds dressing in place on two sides
- » Easily cut to any size
- » Ideal for long surgical incisions, lacerations, abrasions, invasive lines
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NON4311	4" x 6" (4" x 3" pad) (10.2 cm x 15.2 cm)	25/bx, 8 bx/cs, 200 ea/cs
NON4313	8" x 6" (8" x 3" pad) (20.3 cm x 15.2 cm)	25/bx, 4 bx/cs, 100 ea/cs



NON4311

Protective Sleeves

93% Cotton, 7% Lycra® Spandex.

- » Ideal for mild compression
- » Provides warmth
- » Acts as a protective shield for fragile skin
- » Available in white and beige
- » Not made with natural rubber latex

Item No.	Description	Pkg.
Arm Sleeve w/ Thumb Loop		
NONSLEEVE	18" (45.7 cm), White	1 pair/pk
NONSLEEVEVETSM	15.5" (39.3 cm), Beige	1 pair/pk
NONSLEEVEVET	16.5" (41.9 cm), Beige	1 pair/pk
NONSLEEVEVETLG	18" (45.7 cm), Beige	1 pair/pk
NONSLEEVEVETXL	20" (50.8 cm), Beige	1 pair/pk

Leg Sleeve

NONSLEEVEL	15" (38.1 cm), White	1 pair/pk
NONSLEEVELT	15" (38.1 cm), Beige	1 pair/pk



NONSLEEVE



NONTUB112



NONTRI200

Tubular Gauze Bandages

100% knitted cotton.

- » Available in a range of sizes for contoured parts of the body
- » Useful for dressing retention, simple fractures, splinting
- » Can be used with finger applicator NONTUB12
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NONTUB058	Small Finger, Toe, ¾" x 50-yd. (1.6 cm x 45.7 m)	1/bx
NONTUB100	Large Finger, Toe, 1" x 50-yd. (2.5 cm x 45.7 m)	1/bx
NONTUB112	Small Arm, Hand, Foot, 1½" x 50-yd. (3.8 cm x 45.7 m)	1/bx
NONTUB258	Arm, Lower Leg, Small Thigh, 2¾" x 50-yd. (7.9 cm x 45.7 m)	1/bx
NONTUB358	Leg, Thigh, Head, Shoulder, 3¾" x 50-yd. (11.1 cm x 45.7 m)	1/bx
NONTUB12	Finger Applicator	6/bx

Triangular Bandage

Non-sterile, individually packaged.

- » Includes 2 safety pins
- » Conveniently unfolds for arm sling
- » Ideal for EMS and First Aid Kits
- » Not made with natural rubber latex

Item No.	Description	Pkg.
NONTRI200	36" x 36" x 51" (91.4 x 91.4 x 129.5 cm)	12/bg, 240/cs

Woven Gauze Sponges	Medline	J&J	Covidien (Kendall)	Cardinal Health (Allegiance)	Medline Essentials	McKesson	Dukal	Covidien Dermacea
2" x 2", 8-ply, Sterile 2's 50 pk/bx, 30 bx/cs, 3000/cs	NON21420	2318	1806	GZ2208-2	PRM21420	16-42228		441204
3" x 3", 12-ply, Sterile 2's 40 pk/bx, 30 bx/cs, 2400/cs	NON21422	1903	N/A	GZ3312-2	PRM21422			441208
4" x 4", 8-ply, Sterile 2's 50 pk/bx, 12 bx/cs, 1200/cs	NON21448	2315	2187	C-SG448ZS	PRM21448	16-42448		442213
4" x 4", 12-ply, Sterile 2's 25 pk/bx, 24 bx/cs, 1200/cs	NON21424	2317	3033	C-SG44121S	PRM21424	16-42442		442212
4" x 4", 12-ply, Sterile 10's 10/tray, 128 pk/cs, 1280/cs	NON21426	2322 960/cs	6939	GZ4412-10	PRM21426			442214
4" x 4", 16-ply, Sterile 10's 10/tray, 128 pk/cs, 1280/cs	NON21428	2329 720/cs	7605	GZ4416-10	PRM21428			442215
4" x 4", 8-ply, Sterile 10's 10/tray, 128 pk/cs, 1280/cs	NON21425	2321	6318					
Gauze 4" x 8", 12-ply, Sterile 2's 25 pk/bx, 20 bx/cs, 1000/cs	NON21842	2259					6842	

Medline Essentials Woven Gauze Sponges	Medline	Henry Schein Inc.	Dukal	McKesson	Covidien Dermacea
2" x 2", 8-ply 50 pk/bx, 30 bx/cs, 1500 pk/cs	PRM2208		6208	16-4228	
3" x 3", 12-ply 40 pk/bx, 30 bx/cs, 1200 pk/cs	PRM3312		6312		
4" x 4", 12-ply 25 pk/bx, 24 bx/cs, 600 pk/cs	PRM4412		6412	16-4242	441000
4" x 4", 8-ply 50 pk/bx, 12 bx/cs, 600 pk/cs	PRM4408		6408	16-4248	441001
4" x 4", 12-ply 100 pk/bx, 12 bx/cs, 1200 pk/cs	PRM21423			16-42441	
2" x 2" Cotton Filled Dental Sponges 200/sleeve, 25 pk/cs, 5000/cs	PRM22DS	90-8221-I / 90-8223 / 90-8225-I			

Non-Sterile Woven Gauze Sponges	Medline	J&J	Covidien	Cardinal Health	Medline Essentials	McKesson	Dukal	Covidien Dermacea
2" x 2", 8-ply, N/S 200/pk, 25 pk/cs, 5000/cs	NON25208	7612	2146	GZ2208	PRM21208C		2085	442308
2" x 2", 12-ply, N/S 200/pk, 40 pk/cs, 8000/cs	NON25212		2252	GZ2212	PRM21212C		2128	441218
3" x 3", 12-ply, N/S 200/pk, 20 pk/cs, 4000/cs	NON25312		2346	GZ3312	PRM21312C	33122000	6312	441219
4" x 4", 8-ply, N/S 200/pk, 20 pk/cs, 4000/cs	NON25408	7610	2556	C-NSG448E	PRM21408C	44082000	4084	441217
4" x 4", 12-ply, N/S 200/pk, 10 pk/cs, 2000/cs	NON25412	7623	2634	GZ4412	PRM21412C	44122000	4122	441215
4" x 4", 16-ply, N/S 200/pk, 10 pk/cs, 2000/cs	NON25416		2733	GZ4416	PRM21416C	44162000		441216
8" x 4", 12-ply, N/S 200/pk, 10 pk/cs, 2000/cs	NON25812		2835	GZ8412				

Cotton Drain Sponges	Medline							
4" x 4", 16-ply, Sterile 1/pk, 25/bx, 12 bx/cs	NON255000							

X-Ray Detectable Gauze Sponges	Medline	J&J	Covidien	Cardinal Health	Medline Essentials	Covidien Dermacea
4" x 4", 16-ply, Sterile 10's 10/tray, 128 trays/cs, 1280/cs	NON21430LF	2515 720/cs	7317 1350/cs	GT4416-10	PRM21430LF	441002
8" x 4", 12-ply, Sterile 10's 10/tray, 80 trays/cs, 800/cs	NON21432LF	2520 480/cs	7318	GT8412-10 48 trays/cs, 480/cs	PRM21432LF	441602
4" x 4", 32-ply, Sterile 10's 10/tray, 48 trays/cs, 480/cs	NON21433		7320			
8" x 4", 24-ply, Sterile 10's 10/tray, 20 trays/cs, 200/cs	NON21438		7321	C-XR4824S	PRM21438LF	441606

Non-Sterile X-Ray Detectable Gauze Sponges	Medline	J&J	Covidien		Medline Essentials	Covidien Dermacea
4" x 4" 16 ply, N/S 10 sponges/band, 100/sleeve, 2000/cs	NON25430		7114		PRM25430	441601
8" x 4" 12 ply, N/S 10 sponges/band, 100/sleeve, 2000/cs	NON25432		7462		PRM25432	441603

Deluxe Non-Woven Gauze Sponges	Medline	J&J "Nu-Gauze™"				
4" x 4", 4-ply, Sterile 2's 50 pk/bx, 12 bx/cs, 1200/cs	NON264442	2337				
2" x 2", 3-ply, N/S 200/pk, 40 pk/cs, 8000/cs	NON26223	7635				
2" x 2", 4-ply, N/S 200/pk, 40 pk/cs, 8000/cs	NON26224	7632				
3" x 3", 4-ply, N/S 200/pk, 20 pk/cs, 4000/cs	NON26334	7633				
4" x 4", 3-ply, N/S 200/pk, 10 pk/cs, 2000/cs	NON26443	7636				
4" x 4", 4-ply, N/S 200/pk, 10 pk/cs, 2000/cs	NON26444	7634				

Non-Woven Gauze Sponges	Medline	Dukal	Covidien "Versalon®"	Cardinal Health	McKesson	
2" x 2", 4-ply, Sterile 2's 50 pk/bx, 30 bx/cs, 3000/cs	NON21224	7224P	8042	NW2204-2		
3" x 3", 4-ply, Sterile 2's 40 pk/bx, 30 bx/cs, 2400/cs	NON21334	6123/7334P	8043	NW3304-2 1200/cs		
4" x 4", 4-ply, Sterile 2's 50 pk/bx, 12 bx/cs, 1200/cs	NON21444	7444P	8044	NW4404-2		
4" x 4", 6-ply, Sterile 2's 25 pk/bx, 12 bx/cs, 600/cs	NON21446		7084	NW4406-2		
4" x 4", 4 ply, Sterile 10's 10/tray 128 trays/cs	NON21449	6125	8047	N/A		
2" x 2", 4-ply, N/S 200/pk, 40 pk/cs, 8000/cs	NON25224	6112	9022	NW2204		
3" x 3", 4-ply, N/S 200/pk, 20 pk/cs, 4000/cs	NON25334	6113	9023	C-NWS334		
4" x 4", 4-ply, N/S 200/pk, 10 pk/cs, 2000/cs	NON25444	6114	9024	NW4404	94442000	

Medline Essential Non-Woven Gauze Sponges	Medline	Dukal		Cardinal Health		Covidien Dermacea
2" x 2", 4-ply, Sterile 2's 50 pk/bx, 30 bx/cs, 1500 pk/cs	PRM21224	6122				441403
4" x 4", 4-ply, Sterile 2's 50 pk/bx, 12 bx/cs, 600 pk/cs	PRM21444	6124				441412
2" x 2", 4-ply, N/S 200/sleeve, 40 sleeves/cs, 8000/cs	PRM25224	2128		C-NWS224		441400

Bulkee Super Fluff Sponges	Medline	J&J "Kling Fluff®"	Covidien "Kerlix Super™"	Cardinal Health			
6" x 6¾" Diagonal, Sterile 2's 20/bx, 12 bx/cs, 480/cs	NON25852	2335 20/bx, 9 bx/cs	2585	C-SUS662S			
6" x 6¾" Diagonal, Sterile 5's 5/tray, 120 trays/cs, 600/cs	NON25853	2336 5/pk, 90 pk/cs	3085	FLFM-5			
6" x 6¾" Diagonal, Sterile 10's 10/tray, 48 trays/cs, 480/cs	NON25854	2334	7310	C-SUS66ZS			
Large, 7¾" x 8¾" Diagonal Sterile 2's 25 pk/bx, 16 bx/cs, 400/cs	NON25856		4308				
4" x 4", 12 ply Sponge, Sterile 2's 25 pk/bx, 24 bx/cs, 1200/cs	NON25844		5072				
4" x 4", 12 ply Sponge, Sterile 10's 10/tray, 128 trays/cs, 1280/cs	NON25884		6120				

Drain Sponges	Medline	J&J	Covidien	Cardinal Health	McKesson		
4" x 4", 6-ply, Sterile 2's 25 pk/bx, 12 bx/cs, 600/cs	NON256000	2391	7086	C-NWDS226S	16-42426		
2" x 2", 6-ply, Sterile 2's 35 pk/bx, 20 bx/cs, 1400/cs	NON256001	2392	7087	VS2206-2	16-4226		

ABD Pads	Medline	J&J "Surgipad"	Covidien "Tendersorb/Curity"	Cardinal Health	McKesson		
5" x 9", Sterile 1's 25/bx, 16 bx/cs, 400/cs	NON21450	2145	9190/7196 256/cs/400/ cs	AB0509-1	16-4250		
8" x 7½", Sterile 1's 20/bx, 12 bx/cs, 240/cs	NON21453	2144 200/cs 240/cs	9192/7197 500/cs	AB0875-1			
8" x 10", Sterile 1's 18/bx, 20 bx/cs, 360/cs	NON21454	2148 320/cs	9194/7198 144/cs/360/ cs	AB0810-1	16-4254		
5" x 9", N/S Bulk (576/cs)	NON21451		8190/6196	AB0509			
8" x 7½", N/S Bulk (576/cs)	NON21452		8192/6197	C-ABP758			
8" x 10", N/S Bulk (432/cs)	NON21456		8194/6198	C-ABP810			

CURAD Non-Adherent Pads	Medline	J&J "Release"	Covidien (Kendall) "Telfa®"	Cardinal Health (Allegiance)	McKesson	Medline Essentials	Dukal	Covidien LP
2" x 3", Sterile 1's 100/bx, 12 bx/cs, 1200/cs	NON25700	2052	1961 2400/cs	NA0302-1		PRM25700	123	
2" x 3", Sterile 1's with Adhesive Tab 100/bx, 12 bx/cs, 1200/cs	NON25700A						125	
3" x 4", Sterile 1's 100/bx, 12 bx/cs, 1200/cs	NON25710	2054/4731	2132/1050 2400/cs	NA0304-1	16-4293	PRM25710	134	7643
3" x 4", Sterile 1's with Adhesive Tab 100/bx, 12 bx/cs, 1200/cs	NON25710A						124	
3" x 8", Sterile 1's 50/bx, 12 bx/cs, 600/cs	NON25720	2058	1238	NA0308-1		PRM25720	138	

Burn Dressing	Medline	Covidien "Curity Burn"	Medical Action "Burn Dressings"	De Royal "Burn Dressings"
18" x 18", 10-ply, Sewn 1/pk, 80 pk/cs	NON7911	7911	1291815	10-5107
18" x 18", 2-ply, Sewn 2/pk, 100 pk/cs	NON7912		1291830	
18" x 18", 4-ply, Sewn 5/pk, 40 pk/cs	NON7913		1292705	10-5175
18" x 18", 50-ply, Unsewn 1/pk, 20 pk/cs	NON7915	7960	1291828	10-5106
36" x 36", 18-ply, Unsewn 1/pk, 20 pk/cs	NON7916		1291872	10-5147
6" x 6", 45-ply, Sewn 6/tray, 20 tray/cs	NON216645	2963		

Medline SORBEX™ Absorbent Dressing, Sterile Not made with natural rubber latex	Medline	Smith and Nephew "Exu-Dry"		
3" (7.62 cm) Disc 100/cs	NONEXDISC	5999003		
2 x 3" (5.08 x 7.62 cm) Slit, 25/bx, 2 bx/cs, 50/cs	NONEX0023F	5999PTM		
3 x 4" (7.62 x 10.2 cm) Full, 50/bx, 2 bx/cs, 100/cs	NONEX0034	5999034		
4 x 6" (10.2 x 15.2 cm) Full, 10/bx, 12 bx/cs, 120/cs	NONEX0046	5999004120		
6 x 9" (15.2 x 22.8 cm) Full, 12/bx, 4 bx/cs, 48/cs	NONEX0069	5999006		
9 x 15" (22.8 x 38.1 cm) Full 30/cs	NONEX0915	5999009		
15 x 24" (38.1 x 60.9 cm) Full 30/cs	NONEX1524	5999024		

Cotton Balls	Medline			Cardinal Health
1" diameter 2000/bg, 2 bg/cs, 4000/cs	MDS21460			C-COTBM
1¼" diameter 1000/bg, 2 bg/cs, 2000/cs	MDS21462			C-COTBL

Eye Pads	Medline	J&J	Covidien	Cardinal Health
1½" x 2½" Pad, Sterile 1's 50/bx, 12 bx/cs, 600/cs	NON21600	8771	2841	EP2581
2½" x 2½" Pad, Sterile 1's 50/bx, 12 bx/cs, 600/cs	NON21601	8773		EP2582

Post-Op Sponges	Medline	J&J "Topper"	Covidien "Cover"	Cardinal Health
4" x 4", Sterile 2's, 25 pk/bx, 24 bx/cs, 1200/cs	NON21442	2436	2913	C-SG44S
4" x 3", Sterile 2's, 25 pk/bx, 24 bx/cs, 1200/cs	NON21441	2435/2345	3157	CV0403-2

CURAD Packing Strips	Medline	Covidien "Curity"	Cardinal Health	McKesson	Dynarex	
Plain, 1" x 5-yd.	NON255015	7633	59321		3403	
Plain, 2" x 5-yd.	NON255025	7634	59421		3404	
Plain, ½" x 5-yd.	NON255125	7632	C-PG125P	61-59220	3402	
Plain, ¼" x 5-yd.	NON255145	7631	59121	61-59120	3401	
Iodoform, 1" x 5-yd.	NON256015	7833	59346		3413	
Iodoform, 2" x 5-yd.	NON256025	7834	59446		3414	
Iodoform, ½" x 5-yd.	NON256125	7832	59246		3412	
Iodoform, ¼" x 5-yd.	NON256145	7831	59146	61-59145	3411	

CURAD Oil Emulsion Gauze	Medline	Covidien	Cardinal Health	McKesson		J&J "Adaptic"
3" x 8" 3/pk, 108/bx, 648/cs	CUR250383	6115				2013
3" x 3" 50/bx, 600/cs	CUR250330	6112		61-77041		2012
3" x 8" 1/pk, 24/bx, 144/cs	CUR250381	6113	PETROCL38	61-77042		2015
3" x 16" 36/bx, 216/cs	CUR250317	6114				2014

CURAD Xeroform Gauze	Medline	Covidien	Cardinal Health	McKesson		
Overwrap, 5" x 9" 50/bx, 200/cs	CUR253590W	8884431605	C-WMD18P			
Overwrap, 1" x 8" 50/bx, 200/cs	CUR253180W	8884431302				
5" x 9" 50/bx, 200/cs	CUR253590	8884433605	C-WMD59			
2" x 2" 25/bx, 150/cs	CUR253220	8884433400	C-WMD22			
4" x 4" 25/bx, 150/cs	CUR253440	8884433500				
1" x 8" 50/bx, 200/cs	CUR253180	8884433301		61-77034		

CURAD Petrolatum Gauze	Medline	Covidien	Cardinal Health			
Overwrap, 1" x 36" (12/bx, 72/cs)	CUR251136W	8884422600				
Overwrap, ½" x 72" (12/bx, 72/cs)	CUR251272W	8884421600				
Overwrap, 3" x 36" (12/bx, 72/cs)	CUR251336W	8884425600				
Overwrap, 3" x 9" (12/bx, 72/cs)	CUR251390W	8884423600				
Overwrap, 3" x 18" (12/bx, 72/cs)	CUR251318W	8884424600				
Overwrap, 1" x 8" (50/bx, 200/cs)	CUR251180W	8884427601				
1" x 36" (12/bx, 72/cs)	CUR251136	8884412600				
½" x 72" (12/bx, 72/cs)	CUR251272	8884411600				
3" x 36" (12/bx, 72/cs)	CUR251336	8884415600				
3" x 9" (50/bx, 200/cs)	CUR251390	8884413605	C-WNM39			
1" x 8" (50/bx, 200/cs)	CUR251180	8884417601				
3" x 18" (12/bx, 72/cs)	CUR251318	8884414600				
6" x 36" (12/bx, 72/cs)	CUR251636	8884416600				

CURAD Plastic Adhesive Bandages	Medline	J&J "Band-Aid®"	Covidien (Kendall) "Curity"	Cardinal Health	McKesson
¾" x 3" 100/bx, 12 bx/cs, 1200/cs	NON25500	4634	3065 3600/cs	BS3403 9600/cs	
1" x 3" 100/bx, 12 bx/cs, 1200/cs	NON25600		3270 3600/cs	BS0103 9600/cs	
X-Lg, 2" x 4" 50/bx, 12 bx/cs, 600/cs	NON25504	5716 2" x 4½" 900/cs	7539 1200/cs	BS0245	16-4825
¾" x 1½", Junior 100/bx, 36 bx/cs, 3600/cs	NON25509	4633			
⅞" Spot 100/bx, 12 bx/cs, 1200/cs	NON25501	4930 3600/cs	6403 3600/cs	C-BDPSP	

Medline Essentials Plastic Adhesive Bandages	Medline				McKesson
¾" x 3" 100/bx, 12 bx/cs, 1200/cs	PRM25500				16-4823
1" x 3" 100/bx, 12 bx/cs, 1200/cs	PRM25600				16-4821

CURAD Flex-Fabric Adhesive Bandages	Medline	J&J "Band-Aid"	Covidien	Cardinal Health	McKesson
¾" x 3", Comfort Cloth 100/bx, 12 bx/cs, 1200/cs	NON25650	4434	5024	BF3403 7680/cs	16-4813
1" x 3", Comfort Cloth 100/bx, 12 bx/cs, 1200/cs	NON25660	4444	6135	BF0103 7680/cs	
Knuckle, Comfort Cloth 100/bx, 12 bx/cs, 1200/cs	NON25510	4438	6411	C-BDFKN	16-4814
Fingertip, Comfort Cloth 100/bx, 12 bx/cs, 1200/cs	NON25513	4436	6409	C-BDFFT	
2" x 4", Flex-Fabric 50/bx, 12 bx/cs, 600/cs	NON25524		6136		16-4816/ 16-4817
Spot ⅞" 100/bx, 12 bx/cs, 1200/cs	NON25502				16-4812

CURAD Adhesive Bandages – Sesame Street	Medline	J&J "Band-Aid"	Covidien	Cardinal Health	
¾" x 3" 50/bx, 24 bx/cs, 1200/cs	NON47070	AB67	9519 91300	1076538 1076137	

CURAD Adhesive Bandages – Medtoons	Medline	J&J "Band-Aid"	Covidien	Cardinal Health	
¾" x 3" 50/bx, 24 bx/cs, 1200/cs	NON256130	AB67	9519 91300	1075737 1073737	

CURAD Medi-Strips Wound Closure Strips	Medline	J&J "Proxi-Strip®"	Covidien "Curi-Strip®"	Dukal	3M "Steri-Strip"
¼" x 3" 5/pk, 50 pk/bx, 4 bx/cs, 200 pk/cs	NON250318	100K	9890	5150	R1540
¼" x 3" 3/pk, 50 pk/bx, 4 bx/cs, 200 pk/cs	NON250314	101K	9892	5151	R1541
¼" x 1½" 6/pk, 50 pk/bx, 4 bx/cs, 200 pk/cs	NON250114	102K	9891	5152	R1542
¼" x 4" 10/pk, 50 pk/bx, 4 bx/cs, 200 pk/cs	NON250414	103K	9893	5156	R1546
½" x 4" 6/pk, 50 pk/bx, 4 bx/cs, 200 pk/cs	NON250412	104K	9894	5157	R1547
1" x 5", 4/pk, 25 pk/bx, 4 bx/cs, 100 pk/cs	NON250501				R1548

CURAD Nasal Tube Fastener, Non-Sterile	Medline				Derma Sciences
Adult, One Size Fits Most 10/pk, 5 pk/bx, 50/bx	CURNASTUB50				NG-50

Bulkee II Gauze Bandages	Medline	J&J "Kling Fluff"	Covidien (Kendall) "Kerlix"	Cardinal Health (Allegiance) "Hi-Bulk"	
4.5" x 4.1-yd., Sterile 1's, (100/cs)	NON25865	6930	6715	FR4541-1	
3.4" x 3.6-yd., Sterile 1's, (96/cs)	NON25861	6964	6725	C-FR366S	
2¼" x 3-yd., Sterile 1's, 12/bx, 96 rolls/cs	NON25850			C-FR236S	
4.5" x 4.1-yd., N/S, (100/cs)	NON25855	6937	1892	C-FR446	

Medline Essentials Cotton Gauze Bandages	Medline			Cardinal Health (Allegiance) "Hi-Bulk"	McKesson
4.5" x 4.1-yd., Sterile 1's, (100/cs)	PRM25865			C-FR438S	16-4264

Bulkee Lite® Conforming Bandages	Medline	J&J "Kling"	Covidien "Kerlix Lite™"	Cardinal Health	
2" x 3.5-yd., N/S 12/bx, 8 bx/cs, 96/cs	NON27492	6902	8052	PR0402	
3" x 4.1-yd., N/S 12/bx, 8 bx/cs, 96/cs	NON27493	6903	8053	PR0403	
4" x 4.1-yd., N/S 12/bx, 8 bx/cs, 96/cs	NON27494	6904	8054	PR0404	
6" x 4.5-yd., N/S 6/bx, 8 bx/cs, 48/cs	NON27495	6906	8056	PR0406	
2" x 3.5-yd., Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON27496	6922	8072	PR0402-1	
3" x 4.1-yd., Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON27497	6923	8073	PR0403-1	
4" x 4.1-yd., Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON27498	6924	8074	PR0404-1	
6" x 4.5-yd., Sterile 1's 6/bx, 8 bx/cs, 48/cs	NON27499	6926	8076	PR0406-1	

Sterile Conforming Bandages	Medline	J&J "Sof-Kling"	Covidien (Kendall) "Conform"	Cardinal Health (Allegiance)	
1" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON254955				
2" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON25496	6992	2231	GR0402-1	
3" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON25497	6993 72/cs	2232	GR0403-1	
4" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	NON25498	6994 72/cs	2236	GR0404-1	
6" x 80", Sterile 1's 6/bx, 8 bx/cs, 48/cs	NON25499	6996 48/cs	2238	GR0406-1	

Medline Essential Sterile Conforming Bandages	Medline	Dukal			
2" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	PRM25496	8513			
3" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	PRM25497	8514			
4" x 75", Sterile 1's 12/bx, 8 bx/cs, 96/cs	PRM25498	8515			
6" x 80", Sterile 1's 6/bx, 8 bx/cs, 48/cs	PRM25499	8516			

Non-Sterile Conforming Bandages	Medline	J&J "Sof-Kling"	Covidien "Conform"	Cardinal Health	
1" x 75", N/S 12/bx, 8 bx/cs, 96/cs	NON25491	6981	2239	GR0401	
2" x 75", N/S 12/bx, 8 bx/cs, 96/cs	NON25492	6982	2242	GR0402	
3" x 75", N/S 12/bx, 8 bx/cs, 96/cs	NON25493	6983	2244	GR0403	
4" x 75", N/S 12/bx, 8 bx/cs, 96/cs	NON25494	6984	2247	GR0404	
6" x 80", N/S 6/bx, 8 bx/cs, 48/cs	NON25495	6986	2249	GR0406	

Medline Essential Non-Sterile Conforming Bandages	Medline	Dukal			
2" x 75", N/S 12/bx, 8 bx/cs, 96/cs	PRM25492	8517			
3" x 75", N/S 12/bx, 8 bx/cs, 96/cs	PRM25493	8518			
4" x 75", N/S 12/bx, 8 bx/cs, 96/cs	PRM25494	8519			
6" x 80", N/S 6/bx, 8 bx/cs, 48/cs	PRM25495	8520			

Tubular Gauze	Medline	Medical Action		Derma Sciences	
5/8" x 50-yd.	NONTUB058	58200 61033		GL-241	
1" x 50-yd.	NONTUB100	58201		GL-220	
1½" x 50-yd.	NONTUB112	58202 61024		GL-243	
2½" x 50-yd.	NONTUB258	58203		GL-222	
3½" x 50-yd.	NONTUB358	58204		GL-223	

CURAD Elastic Net, Latex-Free, 25-yds.	Medline	Glenwood/ Western Medical "Surgilast"	Medi-Tech "Spandage" TM	DeRoyal 1% Latex "Stretch Net"	Medical Action "SePro Net"
Fingers, Toes, Wrist, Size 1	NONNET01	GLLF2501	MT01	10-7421	53110LF
Small Hand, Arm, Leg, Foot, Size 2	NONNET02	GLLF2502	MT02	10-7422	53120LF
Medium Hand, Arm, Leg, Size 3	NONNET03	GLLF2503	MT03	10-7433	53130LF
Large Hand, Arm, Leg, Foot, Size 4	NONNET04	GLLF2504	MT04	10-7424	53140LF
Small Head, Shoulder, Thigh, Size 5	NONNET05	GLLF2505	MT05	10-7425	53150LF
Large Head, Shoulder, Thigh, Size 6	NONNET06	GLLF2506	MT06	10-7427	53160LF
Small Chest, Back, Perineum, Axilla, Size 7	NONNET07	GLLF2507	MT07	10-7428	53170LF
Medium Chest, Back, Perineum, Axilla, Size 8	NONNET08	GLLF2508	MT08	10-7429	53180LF
Large Chest, Back, Axilla, Perineum, Size 9	NONNET09	GLLF2509	MT09	10-7430	53190LF
X-Large chest, Back, Axilla, Perineum, Size 10	NONNET10	GLLF25010	MT010	10-7431	53210LF
XX-Large Chest, Back, Axilla, Perineum, Size 11	NONNET11	GLLF25011	MT011	10-7432	53220LF
XXX-Large Chest, Back, Axilla, Perineum, Size 12	NONNET12	GLLF25012	MT012	10-7433	
Custom Size	NONNET22		MT022		

CURAD Pre-Cut Elastic Net, Individual Sizes	Medline		Medi-Tech "Spandage" TM		
Finger, Toe, Wrist, Size 1, 24" Long	NONNET124P		MT1X24		
Small Hand, Arm, Leg, Foot, Size 2, 24" Long	NONNET224P		MT2X24		
Medium Hand, Arm, Leg, Foot, Size 3, 24" Long	NONNET324P		MT3X24		
Large Hand, Arm, Leg, Foot, Size 4, 24" Long	NONNET424P		MT4X24		
Small Head, Shoulder, Thigh, Size 5, 24" Long	NONNET524P		MT5X24		
Large Head, Shoulder, Thigh, Size 6, 24" Long	NONNET624P		MT6X24		
Large Head, Shoulder, Thigh, Size 6, 36" Long	NONNET636P		MT6X36		
Small Chest, Back, Perineum, Axilla, Size 7, 24" Long	NONNET724P		MT7X24		
Small Chest, Back, Perineum, Axilla, Size 7, 36" Long	NONNET736P		MT7X36		
Medium Chest, Back, Perineum, Axilla, Size 8, 24" Long	NONNET824P		MT8X24		
Large Chest, Back, Perineum, Axilla, Size 9, 24" Long	NONNET924P		MT9X24		
X-Large Chest, Back, Perineum, Axilla, Size 10, 24" Long	NONNET1024P		MT10X24		
XX-Large Chest, Back, Perineum, Axilla, Size 11, 24" Long	NONNET1124P		MT11X24		
XXX-Large Chest, Back, Perineum, Axilla, Size 12, 24" Long	NONNET1224P		MT12X24		
Custom Size, Size 22, 24" Long	NONNET2224P		MT22X24		

Matrix Elastic Bandage, Sterile Not made with natural rubber latex	Medline	Hartmann		Cardinal Health (Allegiance)	Avcor	Tetra
2" x 5-yd. (20/cs)	DYNJ05152LF	59720000		23593-12LF	12LF	66112S
3" x 5-yd. (20/cs)	DYNJ05153LF	59730000		23593-13LF	13LF	66113S
4" x 5-yd. (20/cs)	DYNJ05154LF	59740000		23593-14LF	14LF	66114S
6" x 5-yd. (20/cs)	DYNJ05156LF	59760000		23593-16LF	16LF	66116S
6" x 10-yd. (20/cs)	DYNJ05157LF	59790000				S61126
6" x 15-yd. (20/cs)	DYNJ05158LF			23593-S960	96LF	S61116

Matrix Elastic Bandage, Non-Sterile Not made with natural rubber latex	Medline	Hartmann		Cardinal Health	Avcor	Tetra
2" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS087002LF	59120000		23593-02LF	02LF	661120
3" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS087003LF	59130000		23593-03LF	03LF	661130
4" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS087004LF	59140000		23593-04LF	04LF	661140
6" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS087006LF	59160000		23593-06LF	06LF	661160
4" x 10-yd. 4/bx, 5 bx/cs, 20/cs	MDS087104LF	59180000			940	661124
6" x 10-yd. 4/bx, 5 bx/cs, 20/cs	MDS087106LF	59190000				661126
6" x 15-yd. 4/bx, 5 bx/cs, 20/cs	MDS087156LF			23593-960	960	661116

Soft-Wrap, Sterile Not made with natural rubber latex	Medline	Hartmann		Cardinal Health	Avcor	Tetra
2" x 5-yd. (20/cs)	DYNJ05125	27200000			23599-12LF	01312S
3" x 5-yd. (20/cs)	DYNJ05126	27300000		23700-63LF	23599-13LF	01313S
4" x 5-yd. (20/cs)	DYNJ05131	27400000			23599-14LF	01314S
6" x 5-yd. (20/cs)	DYNJ05132	27600000		23700-66LF	23599-16LF	01316S

Soft-Wrap, Non-Sterile Not made with natural rubber latex	Medline	Hartmann	Covidien	Cardinal Health	Avcor	Tetra
2" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS046002	38200000	3302	23606-020	23599-02LF	15120
3" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS046003	39300000	3550	23606-030	23599-03LF	15130
4" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS046004	38400000	3616	23606-040	23599-04LF	15140
6" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS046006	39600000	3930	23606-060	23599-06LF	15160

Swift-Wrap, Sterile	Medline	Hartmann			Avcor	Tetra	
2" x 5-yd. (20/cs)	DYNJ05144	59920000			23599-120	S620LF	
3" x 5-yd. (20/cs)	DYNJ05145	59930000			23599-130	S630LF	
4" x 5-yd. (20/cs)	DYNJ05146	59940000				S640LF	
6" x 5-yd. (20/cs)	DYNJ05147	59960000				S660LF	

Swift -Wrap, Non-Sterile	Medline	Hartmann	Covidien	McKesson	Avcor	Tetra	
2" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS077002	59520000	4202		23599-32LF	6620LF	
3" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS077003	59530000	4203	16-1033-3	23599-33LF	6630LF	
4" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS077004	59540000	4204	16-1033-4	23599-34LF	6640LF	
6" x 5-yd. 10/bx, 5 bx/cs, 50/cs	MDS077006	59560000	4206	16-1033-6	23599-36LF	6660LF	

Sure-Wrap, Non-Sterile Not made with natural rubber latex	Medline	Hartmann	Covidien	McKesson	Avcor	Tetra	Cardinal Health
2" x 5-yd., Beige 10/bx, 5 bx/cs, 50/cs	MDS055002	54200000	4422		23599-021	13120	23602-02LF
3" x 5-yd., Beige 10/bx, 5 bx/cs, 50/cs	MDS055003	54300000	4423		23599-031	13130	23602-03LF
4" x 5-yd., Beige 10/bx, 5 bx/cs, 50/cs	MDS055004	54400000	4424	13-254	23599-041	13140	23602-04LF
6" x 5-yd., Beige 10/bx, 5 bx/cs, 50/cs	MDS055006	54600000	4426		23599-061	13160	23602-06LF
2" x 5-yd., White 10/bx, 5 bx/cs, 50/cs	MDS057002	16200000				51020	23602-021
3" x 5-yd., White 10/bx, 5 bx/cs, 50/cs	MDS057003	16300000				51030	23602-031
4" x 5-yd., White 10/bx, 5 bx/cs, 50/cs	MDS057004	16400000				51040	23602-041
6" x 5-yd., White 10/bx, 5 bx/cs, 50/cs	MDS057006	16600000				51060	23602-061

Wytex Undercast Padding, Sterile	Medline	J&J	Covidien				
2" x 4-yd. 20 rolls/cs	DYNJ066002	9082S	2283				
3" x 4-yd. 20 rolls/cs	DYNJ066003	9083S	2394				
4" x 4-yd. 20 rolls/cs	DYNJ066004	9084S	2502				
6" x 4-yd. 20 rolls/cs	DYNJ066006	9086S	2554				

Wytex Undercast Padding, Non-Sterile	Medline	J&J	Covidien				
2" x 4-yd. 24/bg, 3 bg/cs, 72/cs	MDS066002	9082	1418				
3" x 4-yd. 12/bg, 6 bg/cs, 72/cs	MDS066003	9083	2059				
4" x 4-yd. 12/bg, 6 bg/cs, 72/cs	MDS066004	9084	3175				
6" x 4-yd. 6/bg, 6 bg/cs, 36/cs	MDS066006	9086	3489				

Syn-Tex Undercast Padding, Synthetic, Non-Sterile	Medline	J&J		Covidien	3M
2" x 4-yd. 12/bg, 6 bg/cs, 72/cs	MDS067002	9052		6242	CMW02
3" x 4-yd. 12/bg, 6 bg/cs, 72/cs	MDS067003	9033		6243	CMW03
4" x 4-yd. 12/bg, 6 bg/cs, 72/cs	MDS067004	9034		6244	CMW04
6" x 4-yd. 6/bg, 6 bg/cs, 36/cs	MDS067006	9036		6246	CMW06
8" x 4-yd. 6/bg, 6 bg/cs, 36/cs	MDS067008				

Self-Adherent Bandages, Non-Sterile, Latex-Free	Medline		Andover "Co-Flex LF2"		3M "Coban"
Co-Flex LF2, Stretched, Tan, 1" x 5-yd. 30 rolls/cs	MDS089001		9100		2081
Co-Flex LF2, Stretched, Tan, 2" x 5-yd. 36 rolls/cs	MDS089002		9200		2082
Co-Flex LF2, Stretched, Tan, 3" x 5-yd. 24 rolls/cs	MDS089003		9300		2083
Co-Flex LF2, Stretched, Tan, 4" x 5-yd. 18 rolls/cs	MDS089004		9400		2084
Co-Flex LF2, Stretched, Tan, 6" x 5-yd. 12 rolls/cs	MDS089006		9600		2086

Self-Adherent Bandages, Sterile, Latex-Free	Medline		Andover "Co-Flex LF2"		3M "Coban"
Co-Flex LF2, Stretched, Tan, 2" x 5-yd. 20 rolls/cs	DYNJ089002		9200S		
Co-Flex LF2, Stretched, Tan, 3" x 5-yd. 20 rolls/cs	DYNJ089003		9300S		2083S
Co-Flex LF2, Stretched, Tan, 4" x 5-yd. 20 rolls/cs	DYNJ089004		9400S		2084S
Co-Flex LF2, Stretched, Tan, 6" x 5-yd. 18 rolls/cs	DYNJ089006		9600S		2086S

Self-Adherent Bandages, Sterile	Medline			Covidien (Kendall) "Flex-Wrap"	3M "Coban"	
Co-Flex Med, Stretched, Tan, 2" x 5-yds. 20 rolls/cs	DYNJ086002					
Co-Flex Med, Stretched, Tan, 3" x 5-yds. 20 rolls/cs	DYNJ086003			4583S	1583S	
Co-Flex Med, Stretched, Tan, 4" x 5-yd. 20 rolls/cs	DYNJ086004			4584S, 18/pk	1584S, 18/pk	
Co-Flex Med, Stretched, Tan, 6" x 5-yd. 20 rolls/cs	DYNJ086006			4586S, 12/pk	1586S, 12/pk	

Self-Adherent Bandages, Non-Sterile	Medline		Andover "Co-Flex Med"	Covidien "Flex-Wrap"	3M "Coban"	
Co-Flex Med, Stretched, Tan, 1" x 5-yd. 30 rolls/cs	MDS086001		7100TN	4581	1581	
Co-Flex Med, Stretched, Tan, 2" x 5-yd. 36 rolls/cs	MDS086002		7200TN	4582	1582	
Co-Flex Med, Stretched, Tan, 3" x 5-yd. 24 rolls/cs	MDS086003		7300TN	4583	1583	
Co-Flex Med, Stretched, Tan, 4" x 5-yd. 18 rolls/cs	MDS086004		7400TN	4584	1584	
Co-Flex Med, Stretched, Tan, 6" x 5-yd. 12 rolls/cs	MDS086006		7600TN	4586	1586	

CURAD Paper Tape	Medline	J&J "Dermiform"	Covidien "Tenderskin®"	Cardinal Health	3M "Micropore"	McKesson
½" x 10-yd. 24/bx, 10 bx/cs, 240/cs	NON270012	5180	1596	2TPPR00/ 2TPPR00T	1530-0	16-47305
1" x 10-yd. 12/bx, 10 bx/cs, 120/cs	NON270001	5181	1914	2TPPR01T	1530-1 120/cs	16-47310
2" x 10-yd. 6/bx, 10 bx/cs, 60/cs	NON270002	5182	2419	32262	1530-2 60/cs	16-47320
3" x 10-yd., 4/bx, 10 bx/cs, 40/cs	NON270003	5183	3394	2TPPR03	1530-3	
1" x 1.5-yd., 100/bx, 5 bx/cs, 500/cs	NON270001S				1530S-1	
2" x 1.5-yd., 50/bx, 5 bx/cs, 250/cs	NON270002S				1530S-2	

Medline Essential Paper Tape	Medline					McKesson
1" x 10-yd. 12 rolls/bx, 12 bx/cs, 144 rolls/cs	PRM260001					16-47310

CURAD Cloth Tape	Medline	J&J "Dermicel"	Covidien "Curasilk®"	Cardinal Health	3M "Durapore"	McKesson
½" x 10-yd. 24/bx, 10 bx/cs, 240/cs	NON270112	5143	7137	2TRCL00	1538-0 240/cs	16-47205
1" x 10-yd. 12/bx, 10 bx/cs, 120/cs	NON270101	5144	7138	32266 120/cs	1538-1 120/cs	
2" x 10-yd. 6/bx, 10 bx/cs, 60/cs	NON270102	5145	7139	2TRCL02	1538-2 60/cs	16-47120
3" x 10-yd. 4/bx, 10 bx/cs, 40/cs	NON270103	5146 40/cs	7140	2TRCL03	1538-3	
Cloth "Silk" Tape, 1" x 1.5-yd. 100/bx, 5 bx/cs, 500/cs	NON260101S				1538S-1	

CURAD Cloth Tape	Medline					McKesson
1" x 10-yd. 12 rolls/bx, 12 bx/cs, 144 rolls/cs	PRM260101					16-47110

CURAD Transparent Tape	Medline	J&J "Dermiview"	Covidien "Curity®"	Cardinal Health	3M "Transpore"	McKesson
½" x 10-yd. 24/bx, 10 bx/cs, 240/cs	NON270212		8533	2TPPB00	15270-0 240/cs	16-47320
1" x 10-yd. 12/bx, 10 bx/cs, 120/cs	NON270201	5161 120/cs	8534	32271 120/cs	1527-1 120/cs	
2" x 10-yd. 6/bx, 10 bx/cs, 60/cs	NON270202	5162 60/cs	8535	3TPLT02	1527-2 60/cs	
3" x 10-yd. 4/bx, 10 bx/cs, 40/cs	NON270203	5163 40/cs	8536	2TPPB03/ 3TPLT03	1527-3 40/cs	
1" x 1.5-yd. 100/bx, 5 bx/cs, 500/cs	NON270201S			3TPLT00	1527S-1	

Medline Essential Transparent Tape	Medline			Cardinal Health		McKesson
1" x 10-yd. 12 rolls/bx, 12 bx/cs, 144 rolls/cs	PRM260201			2TPPB01		16-47210
2" x 10-yd. 6 rolls/bx, 12 bx/cs, 72 rolls/cs	PRM260202			2TPPB02		

CURAD Elastic Foam Tape	Medline			Cardinal Health		
1" x 5.5-yd. 12 rolls/bx, 6 bx/cs, 72 rolls/cs	NON260601			2TFM01		
2" x 5.5-yd. 6 rolls/bx, 6 bx/cs, 36 rolls/cs	NON260602			2TFM02		
3" x 5.5-yd. 4 rolls/bx, 6 bx/cs, 24 rolls/cs	NON260603			2TFM03		
4" x 5.5-yd. 3 rolls/bx, 6 bx/cs, 18 rolls/cs	NON260604			2TFM04		

CURAD Ortho-Porous Sports Tape	Medline			Cardinal Health		
1" x 10-yd. 12 rolls/bx, 12 bx/cs, 144 rolls/cs	NON260301			2TSCL01		
2" x 10-yd. 6 rolls/bx, 12 bx/cs, 72 rolls/cs	NON260302			2TSCL02		
3" x 10-yd. 4 rolls/bx, 12 bx/cs, 48 rolls/cs	NON260303			2TSCL03		

ADVANCING THE HEALTH OF HEALTHCARE.™

In a complex healthcare world, Medline's mission is to help our customers achieve clinical and financial success. We do that through a personalized approach to listen and better understand our customers' needs in an environment that tirelessly demands lower costs and better outcomes.

We are both a manufacturer and distributor with 1,200 direct sales reps who are dedicated points of contact for customers across the continuum of care. Because of our capabilities and engagement with other clinical and industry thought leaders, we deliver the most robust product portfolio and patient-care solutions to help healthcare providers perform at their very best.

Our responsiveness and commitment to advancing the health of healthcare is evident in our actions every day. This is who we are. This is why customers choose us as their trusted business partner.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US 