




# NEW Products From Cure Medical®



## MORE OPPORTUNITIES TO SUPPORT RESEARCH

Since Cure Medical **donates 10% of net income to medical research** in pursuit of a cure for spinal cord injuries (SCI) and central nervous system disorders (CNS/D), use of the new, extra long (**M14XL**) **Cure Medical® Pocket Catheter** and any **hydrophilic Cure Catheter®** serves as another opportunity to advance research toward a cure.

### CURE MEDICAL® POCKET CATHETER (M14XL)



NEW - Extra Long, 25" M14XL discreetly fits in a pant pocket

Now, the family of U-shaped Cure Medical® Pocket Catheters includes an extra long, 25" option – the M14XL. This unique catheter facilitates easy draining by eliminating the need to transfer from a wheelchair or use additional extension tubing. The package includes lubricant as an added benefit.

### HYDROPHILIC CURE CATHETER®



Features polished eyelets, tip options, gripper sleeve, and purified water pouch



The popular Cure Catheter® line now includes hydrophilic options in 12FR - 16FR with straight or Coude tip designs. The hydrophilic Cure Catheter® enables easy, clean catheter hydration while in the package. Additionally, a textured gripper sleeve allows for controlled, clean use.

## EXCEPTIONAL BENEFITS

All Cure Medical® products reflect innovations and benefits designed to improve the user experience:

- Smooth, polished eyelets on every catheter offers increased comfort
- High quality materials – not made with DEHP\*, BPA, or natural rubber latex – offer peace of mind
- Support of research in pursuit of a cure for SCI and CNS/D possibly offers the best benefit of all!

## FREE SAMPLES

To request free samples of the new (M14XL) Cure Medical® Pocket Catheter, hydrophilic Cure Catheter® or any Cure Medical® product, contact your local supplier of quality healthcare products or visit [www.curemedical.com](http://www.curemedical.com).



## Cure Catheter®

### NOW INCLUDES HYDROPHILIC OPTIONS

Offered as a case of 300 units (10 boxes of 30 units in each box)

Part No.	Description	FR Size
M8	Male, 16", Straight Tip	8
M10	Male, 16", Straight Tip	10
M10C	Male, 16", Coude Tip	10
M12	Male, 16", Straight Tip	12
M12C	Male, 16", Coude Tip	12
M14	Male, 16", Straight Tip	14
M14C	Male, 16", Coude Tip	14
M16	Male, 16", Straight Tip	16
M16C	Male, 16", Coude Tip	16
M18	Male, 16", Straight Tip	18
M18C	Male, 16", Coude Tip	18
HM12	Male, 16", Straight Tip, Hydrophilic	12
HM12C	Male, 16", Coude Tip, Hydrophilic	12
HM14	Male, 16", Straight Tip, Hydrophilic	14
HM14C	Male, 16", Coude Tip, Hydrophilic	14
HM16	Male, 16", Straight Tip, Hydrophilic	16
HM16C	Male, 16", Coude Tip, Hydrophilic	16
F8	Female, 6", Straight Tip	8
F10	Female, 6", Straight Tip	10
F12	Female, 6", Straight Tip	12
F14	Female, 6", Straight Tip	14
F14NC	Female, 6", Straight Tip, No Connector	14
F16	Female, 6", Straight Tip	16
F18	Female, 6", Straight Tip	18
P8	Pediatric, 10", Straight Tip	8
P10	Pediatric, 10", Straight Tip	10
P12	Pediatric, 10", Straight Tip	12
P14	Pediatric, 10", Straight Tip	14

## Cure Catheter® Closed System

Offered as a case of 100 units

Part No.	Description	FR Size
CB8	Closed System Single, Straight Tip	8
CB10	Closed System Single, Straight Tip	10
CB12	Closed System Single, Straight Tip	12
CB14	Closed System Single, Straight Tip	14
CB16	Closed System Single, Straight Tip	16
CS8	Closed System Kit, Straight Tip	8
CS10	Closed System Kit, Straight Tip	10
CS12	Closed System Kit, Straight Tip	12
CS14	Closed System Kit, Straight Tip	14
CS14C	Closed System Kit, Coude Tip	14
CS16	Closed System Kit, Straight Tip	16

## Cure Medical® Catheter Insertion Kit

Offered as a case of 100 units

Part No.	Description
K1	Catheter Insertion Kit
K2-100	Catheter Insertion Kit with universal connector on the bag

Offered as a case of 90 units

Part No.	Description
K2-90	Catheter Insertion Kit with universal connector on the bag

## Cure Medical® Pocket Catheter

### NOW INCLUDES A NEW, EXTRA LONG STYLE!

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
M12U	Male, 16", Straight Tip	12
M14U	Male, 16", Straight Tip	14
M14XL	Male, 25" (Extra Long), Straight Tip with Lubricant and Gripper Sleeve	14
M14UL	Male, 16", with Lubricant	14
M16U	Male, 16", Straight Tip	16

Offered as a case of 90 units (3 boxes of 30)

M14UK	Male, 16", Straight Tip, Kit	14
-------	------------------------------	----

## Cure Twist®

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
T8	Female, 6", Straight Tip	8
T10	Female, 6", Straight Tip	10
T12	Female, 6", Straight Tip	12
T14	Female, 6", Straight Tip	14
T16	Female, 6", Straight Tip	16

## Cure Ultra®

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
ULTRA8	Female, 6", Straight Tip	8
ULTRA10	Female, 6", Straight Tip	10
ULTRA12	Female, 6", Straight Tip	12
ULTRA14	Female, 6", Straight Tip	14
ULTRA16	Female, 6", Straight Tip	16

