



# Cure Ultra<sup>®</sup> FOR MEN AND WOMEN

Pre-lubricated coude tip with a blue control stripe

Since Cure Medical **donates 10% of net income to medical research** in pursuit of a cure for spinal cord injuries (SCI) and central nervous system disorders (CNS/D), use of the new **Cure Ultra<sup>®</sup>** serves as another opportunity to advance research toward a cure.

The sterile, single use Cure Ultra is a ready-to-use catheter featuring many user-driven benefits:

- Available in straight for women, plus straight and coude for men
- The small, flexible package features an easy tear top and enables disposal of minimal material
- Smooth polished eyelets on the pre-lubricated catheter tip increase comfort
- Even, complete distribution of lubricant offers ‘No Drip/No Mess’ comfortable cathing
- A unique ‘No Roll’ connector/funnel end helps increase confidence when draining
- High quality materials – not made with DEHP\*, BPA, or natural rubber latex offer peace of mind
- Support of research in pursuit of a cure for SCI and CNS/D possibly offers the best benefit of all!

The Cure Ultra is offered in 6” length with a straight tip for women and in a 16” length with the option for straight or coude tip for men.

## NEW! ESPECIALLY FOR MEN

The Cure Ultra for men features a proprietary CoverAll™ application process for even distribution of lubricant over the catheter as it is removed from the package – making it ready-to-use with ultra-low resistance as you self-cath. Additionally, it includes a unique gripper for increased ease of use.



## FREE SAMPLES

To request free samples of the new Cure Ultra or any Cure Medical catheter, contact your local distributor of quality healthcare products or visit [www.curemedical.com](http://www.curemedical.com). Thank you for choosing a Cure!



[Click here to view detailed product descriptions and prices.](#) **Call 800-397-5899 to order today!**

## Cure Catheter®

Offered as a case of 300 units (10 boxes of 30 units in each box)

Part No.	Description	FR Size
HM12	Hydrophilic, Male, 16", Straight Tip	12
HM12C	Hydrophilic, Male, 16", Coude Tip	12
HM14	Hydrophilic, Male, 16", Straight Tip	14
HM14C	Hydrophilic, Male, 16", Coude Tip	14
HM16	Hydrophilic, Male, 16", Straight Tip	16
HM16C	Hydrophilic, Male, 16", Coude Tip	16

M8	Male, 16", Straight Tip	8
M10	Male, 16", Straight Tip	10
M10C	Male, 16", Coude Tip	10
M12	Male, 16", Straight Tip	12
M12C	Male, 16", Coude Tip	12
M14	Male, 16", Straight Tip	14
M14C	Male, 16", Coude Tip	14
M16	Male, 16", Straight Tip	16
M16C	Male, 16", Coude Tip	16
M18	Male, 16", Straight Tip	18
M18C	Male, 16", Coude Tip	18

F8	Female, 6", Straight Tip	8
F10	Female, 6", Straight Tip	10
F12	Female, 6", Straight Tip	12
F14	Female, 6", Straight Tip	14
F14NC	Female, 6", Straight Tip, No Connector	14
F16	Female, 6", Straight Tip	16
F18	Female, 6", Straight Tip	18

HP8	Hydrophilic, Pediatric, 10", Straight Tip	8
HP10	Hydrophilic, Pediatric, 10", Straight Tip	10
HP12	Hydrophilic, Pediatric, 10", Straight Tip	12
HP14	Hydrophilic, Pediatric, 10", Straight Tip	14

P8	Pediatric, 10", Straight Tip	8
P10	Pediatric, 10", Straight Tip	10
P12	Pediatric, 10", Straight Tip	12
P14	Pediatric, 10", Straight Tip	14

## Cure Catheter® Closed System

Offered as a case of 100 units

Part No.	Description	FR Size
CB8	Closed System Single, Straight Tip	8
CB10	Closed System Single, Straight Tip	10
CB12	Closed System Single, Straight Tip	12
CB14	Closed System Single, Straight Tip	14
CB16	Closed System Single, Straight Tip	16

CS8	Closed System Kit, Straight Tip	8
CS10	Closed System Kit, Straight Tip	10
CS12	Closed System Kit, Straight Tip	12
CS14	Closed System Kit, Straight Tip	14
CS14C	Closed System Kit, Coude Tip	14
CS16	Closed System Kit, Straight Tip	16

## Cure Medical® Catheter Insertion Kit

Offered as a case of 100 units

Part No.	Description
K1	Insertion Kit
K2-100	Insertion Kit w/universal connector

Offered as a case of 90 units

Part No.	Description
K2-90	Insertion Kit w/universal connector

## Cure Medical® Pocket Catheter

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
M12U	Male, 16", Straight Tip	12
M12ULC	Male, 16", Coude Tip w/Lubricant	12
M14U	Male, 16", Straight Tip	14
M14XL	Male, 25" (Extra Long), Straight Tip	14
M14UL	Male, 16", Straight Tip w/Lubricant	14
M14ULC	Male, 16", Coude Tip w/Lubricant	14
M16U	Male, 16", Straight Tip	16
M16ULC	Male, 16", Coude Tip w/Lubricant	16
M14UK	Male, 16", Straight Tip, Kit	14

## Cure Twist®

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
T8	Female, 6", Straight Tip	8
T10	Female, 6", Straight Tip	10
T12	Female, 6", Straight Tip	12
T14	Female, 6", Straight Tip	14
T16	Female, 6", Straight Tip	16

Offered as a case of 90 units (3 boxes of 30)

T8K	Female, 6", Straight Tip w/Insert. Kit	8
T10K	Female, 6", Straight Tip w/Insert. Kit	10
T12K	Female, 6", Straight Tip w/Insert. Kit	12
T14K	Female, 6", Straight Tip w/Insert. Kit	14
T16K	Female, 6", Straight Tip w/Insert. Kit	16

## Cure Ultra®

Offered as a case of 300 units (10 boxes of 30)

Part No.	Description	FR Size
ULTRA8	Female, 6", Straight Tip	8
ULTRA10	Female, 6", Straight Tip	10
ULTRA12	Female, 6", Straight Tip	12
ULTRA14	Female, 6", Straight Tip	14
ULTRA16	Female, 6", Straight Tip	16
ULTRA M8	Male, 16", Straight Tip	8
ULTRA M10	Male, 16", Straight Tip	10
ULTRA M12	Male, 16", Straight Tip	12
ULTRA M12C	Male, 16", Coude Tip	12
ULTRA M14	Male, 16", Straight Tip	14
ULTRA M14C	Male, 16", Coude Tip	14
ULTRA M16	Male, 16", Straight Tip	16
ULTRA M16C	Male, 16", Coude Tip	16
ULTRA M18	Male, 16", Straight Tip	18
ULTRA M18C	Male, 16", Coude Tip	18