



## DH OFFLOADING WALKER™

- Shear force reduction by stabilizing foot and ankle at 90°
- Wide low-profile rocker bottom provides a smooth roll-through and increased stability for altered gaits
- Designed around Össur's unique Active Offloading Insole
- The softgood's patented inner bootie prevents the foot from migrating on the insole
- Provides a safe environment for the simultaneous use of wound care material
- Create an Instant TCC using a single roll of 2" cast tape for patient compliance

### INDICATIONS

- Designed for use in the treatment of plantar ulcers

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE
DH0200BLK	XSmall	2-4	3.5-5.5
DH0400BLK	Small	4.5-7	6-8
DH0600BLK	Medium	7.5-10.5	8.5-11.5
DH0800BLK	Large	10.5-12.5	11.5-13.5
DH1000BLK	XLarge	12.5+	13.5+

### Replacement Parts

XSMALL	SMALL	MEDIUM	LARGE	XLARGE	DESC.
Call	10958	10959	10960	Call	Softgood
11742C	PP04C	PP06C	11740C	11741C	Insole

Protected by U.S. Patent #: 5,329,705; 5,761,834; D473,654  
Other U.S. and foreign patent pending



## DH OFFLOADING POST-OP STYLE SHOE

- Designed to accommodate either the left or the right foot
- Designed around Össur's patented Active Offloading Insole
- Designed to accommodate bulky dressings
- Super-soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe piece aides in wound accessibility for daily wound site inspection
- Semi-rigid rocker bottom outsole

### INDICATIONS

Designed for use in the treatment of plantar wounds, this shoe can be an excellent alternative to post-op style shoes with static pads (wool, felt, foam, etc.). It accommodates the ambulatory and compliant patient for short-term treatment (1 to 4 months).

PART#	SIZE	MEN'S SHOE	WOMEN'S SHOE
10340	Small	4-6	6.5-8.5
10341	Medium	6.5-8.5	9-11
10342	Large	9-11	12+
10343	XLarge	11.5-14	13+

### Replacement Insoles

XSMALL	SMALL	MEDIUM	LARGE	XLARGE
Call	12365	12366	12367	12368

Patent #: 5,329,705; 5,778,565; 6,792,699

