

DreamStation helps support long-term patient use while creating **new efficiencies for your business**.

Connectivity options offer more monitoring choices than any leading CPAP brand with integrated Bluetooth® on every device, and offering optional Wi-Fi and cellular modems.*

Performance Check simplifies in-home device evaluation and can help reduce the number of normally functioning devices that are returned for service.

Daily Progress Feedback displays a simple trend of patients' nightly use to help track progress on the path to adherence.

DreamStation features

- · Mask type resistance control
- · Mask fit check
- · Oximetry capable
- Encore patient management system with optional Adherence Profiler
- DreamMapper patient self-management system compatibility

Humidification

- · Fixed, adaptive, heated tube
- · Pre-heat humidifier option
- \cdot Modular design with one-piece, dishwasher-safe water chamber

Protection programs

- · Flex Promise program
- · Bi-level Rescue program
- "One-Touch" warranty programs

CPAP features

- C-Flex
- SmartRamp
- · Basic AHI/Leak reporting
- Standard integrated Bluetooth

DreamStation CPAP

	Comfort features		Encore data capture and reporting		
Modes	Features	Flex Comfort	Advanced Event Detection	Flow Waveforms	
Fixed CPAP	N/A	C-Flex	N/A	N/A	

Specifications

DreamStation CPAP		
Pressure range	4 to 20 cm H2O	
Flex pressure relief	0 to 3	
Ramp time	0 to 45 min (5-minute increments)	
Starting ramp pressure	4 cm H2O to CPAP Pressure (CPAP)	
Humidification	Heated humidification: Fixed, Adaptive, Heated Tube modes; Heated Tube (optional)	
Data storage capacity (minimum)	SD card: >1 year On-board: <6 months	
Altitude compensation	Automatic	
Filters	Reusable pollen; disposable ultra-fine	
Device controls	LCD, control dial/push button	
Warranty	2 years (US)	
Dimensions	15.7x19.3x8.4cm (w/o humidifier) 29.7x19.3x8.4cm (w/ humidifier)	
Weight	1.33kg/2.94lbs (w/o humidifier) 1.98kg/4.37lbs (w/humidifier) Incl. power supply	
Electrical requirements	100 – 240 VAC, 50/60 Hz, 2.0-1.0 A	

Ordering information

DreamStation CPAP	Part number
DreamStation CPAP	DSX200S11
DreamStation CPAP with humidifier	DSX200H11
DreamStation CPAP with humidifier and heated tube	DSX200T11
Power supply	Part number
DreamStation 80-watt power supply	1118499
Power cord	1005894

Ordering information

Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubes	Part number
15MM Std Tube – DreamStation - RP	PR15
15MM Heated Tube	HT15
Filters	Part number
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519
Accessories	Part number
Link module	1120293
Nonin SpO2, assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Modems	Part number
Cellular modem, oximetry capable	100600C
Cellular modem, non-oximetry	100605C
Wi-Fi accessory, oximetry capable	100700W
Wi-Fi accessory, non-oximetry	100705W

© 2015 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.



Caution: US federal law restricts these devices to sale by or on the order of a physician.

edoc EB 11/12/15 MCI 4107064 1010 Murry Ridge Lane, Murrysville, PA 15668 800-345-6443 • 724-387-4000