

extriCARE® Dressings

Negative Pressure Wound Therapy System

Dressings Overview

Our anatomically fitted dressings come in a variety of sizes to fit most needs. They feature a non-adhesive contact layer which can be applied directly to the wound bed. This layer is centered on an adhesive film which creates a tight seal, and may promote healing. A drainage tube connects each dressing to a canister connected to our pump to remove wound material. Our foam dressings come in two sizes in a kit that includes adhesive film, connector tubing, and a suction bell.

Features and Benefits:

- 
Versatile: Available in many sizes, can be used with all extriCARE® pumps.
- 
Simple: Easy to use, quick application
- 
Effective: Creates a tight seal, important for creating negative pressure for treatment of wounds.



Extra Large Oval
EC2400-IM



Small Oval
EC2400-SO



Mini Oval
EC2400-MN



Large Oval
EC2400-LO



Sacral
EC2400-S



Foot
EC2400-F



Large Foam Dressing
EC-FOAM-L



Small Foam Dressing
EC-FOAM-S

LB30.0004 Rev H 20160115



ALLEVA MEDICAL CO.

Suite M-Q, 12th Floor, Kings Wing Plaza 2, 1 On Kwan St. Shek Mun, Shatin, N.T., Hong Kong
Tel: +852 3166 7200 Fax: +852 2355 7663 www.allevamedical.com



© 2018