



Our exceptionally easy-to-use **EZ-SHAMPOO® Hair Washing Tray** is contoured to fit the neck comfortably and rest on the client's shoulders for convenient shampooing in a seated position.

- Compact and lightweight design
- Reduces awkward lifting and bending for caregiver
- Raised edge prevents water from spilling over
- 35" stabilizing strap included
- Designed for use with chair or wheelchair
- Made of durable plastic and designed for years of trouble-free use
- One-year limited warranty

The **EZ-SHAMPOO® Hair Washing Basin** allows for a refreshing, in-bed shampoo. Comfortable and convenient, the inflatable basin provides secure, form fitting support that cushions the head and shoulders and includes an adjustable, built-in headrest. A caregiver can quickly shampoo and rinse a patient's hair without awkward lifting and bending. An attached one-way drain hose easily and conveniently removes suds and water.

- Inflatable and lightweight for portability and convenient storage
- Features a double-tube design to help prevent splashes and spills
- Flexible drain hose with custom stopper
- 24" wide x 20" long x 8" deep, inflated size
- 6½" neck opening (one size fits most)
- Made of heavy-duty vinyl
- Latex free
- One-year limited warranty



Designed to be used where a water source is unavailable, the **EZ-SHOWER** is an overhead bedside shower and is the ideal accessory to the EZ-SHAMPOO® Hair Washing Basin.

- Holds up to 2½ gallons of warm tap water – ample enough for a complete shampoo and rinse
- Flexible hose directs a gentle stream of water exactly where it's needed
- Compact and easy to store when not in use
- Includes a flexible 2' hose with an on/off switch to control water flow
- Built-in handle to carry to the bedside
- Hangs on a bedpost or an I.V. pole
- Made of heavy-duty vinyl
- Latex free
- One-year limited warranty

customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com