CONSULT® diagnostics

MSKESSON

STREP A TESTS

McKesson Consult Strep A tests are designed to aid healthcare professionals in the early detection of group A streptococcal antigen, commonly termed strep throat. The external positive and negative controls help produce accurate results. Consult Strep A test kits quickly identify the streptococcus organism and allow medical professionals to form an accurate patient diagnosis.

FEATURES AND BENEFITS

MFR# 4999 and 5003

- Results available in five minutes for quick patient diagnosis
- 97% sensitivity and 95% specificity to help ensure accuracy of results
- Available in both 25 and 50 count boxes, for caregiver convenience and economy
- Individually packaged and designed for single patient use
- CLIA waived for simple administration
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on Care SM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2019 McKesson Medical-Surgical Inc.

MSKESSON



FEATURES AND BENEFITS

MFR# 181-34125

- Color-coded control and test lines for easy-to-read results
- 96.2% sensitivity and 98.7% specificity to help ensure accuracy of results
- 24-month shelf life
- External positive and negative quality controls included
- Individually packaged and designed for single patient use
- CLIA waived for simple administration
- Not made with natural rubber latex



MFR#	DESCRIPTION	QTY
4999	Strep A Test, Dipstick	50 ea/bx
5003	Strep A Test, Dipstick	25 ea/bx
181-34125	Strep A Test, Cassette	25 ea/bx

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908