

FRIO® Facts

- ~Keeps insulin cool and safe
- ~Ice packs or refrigeration NOT needed - just activate with water!
- ~Light and compact
- ~Perfect for power outages
- ~Reusable (a very **green** product)
- ~Perfect for travel. TSA Friendly!
- ~Comes in six colors
- ~Sizes for vials, pens or pumps
- ~Great in **HOT** weather!!! FRIOS® have been used in temperatures in excess of 120°F with complete success.
- ~Tests prove FRIOS® keep contents at around 78°F (25°C) for a minimum of two days in temps of 100°F (38°C). When over 100°F, soak more frequently.
- ~FRIOs® provide insulation, helping to protect insulin from freezing at cold temps.
- ~Low Cost!!!
- ~FRIO® offers convenience, freedom and, most importantly, peace of mind!



Scan to view our Website!

**Activates
with
Water!**



Pump



With FRIO® you have
**Freedom
to Travel the World!**

www.FRIOCase.com



Duo

Resellers Welcome

FRIO®

distributed by

ReadyCare, LLC

1301 Ramsay Circle

Walnut Creek, CA 94597

(925) 937-5288

(888) 909-7737

**FRIO® Insulin
Cooling Cases**

*No Ice Packs...
Ever Again!*



**Convenience,
Freedom &
Peace of Mind**

www.FRIOCase.com

FRIO® Cooling Case

The FRIO® Cooling Case will improve the life of anyone dependent on insulin or many other temperature-sensitive medications.

The FRIO® is a reusable evaporative cooler whose cooling properties do not come from an icepack - or anything that needs refrigeration. Its cooling properties come from the evaporation of water. When activated, it keeps its low temperature for a minimum of two days, even in temperatures of 100°F.

For a FRIO® to be activated, simply soak it in water for 5-10 minutes. It retains its cooling properties for a minimum of two full days. After that, re-soak it in water and it is good for another two days - and so on, and so on...

While a traditional medication cooling case gives you a six to eight hour "chain" to refrigeration, FRIO® cuts that chain. The FRIO® offers the security of not being dependent on refrigeration for the life-critical cooling of your insulin.* It gives you freedom and flexibility to enjoy life as you choose.

Furthermore, because you do not need an icepack or the surrounding insulating materials, the FRIO® wallet is lighter and smaller than traditional medication-cooling cases. There are various sizes of FRIO® available - some small enough to fit in your jean pocket!



Large

PERFECT FOR
*Emergency Kits
Camping & Bike Trips
Working Outdoors
Picnics & Theme Parks
Travel & Daily Use*

These make great gifts!

These are some examples of how the FRIO® will change your life...

Travel~ When you travel, you need to find a way to re-freeze a traditional cooling case throughout the day or overnight. Although many hotels have in-room refrigerators/freezers, many do not. With the FRIO®, you need not worry if there is a refrigerator. All you need is water! The FRIO® is TSA friendly.

Camping~ With a traditional cooling case, camping or backpacking becomes difficult or even impossible. After soaking your FRIO® in water, you can go on a full two-day trip. As long as there is water, the trip can be as long as you want it to be - with no worries!

Confidentiality & Convenience~ It is difficult to maintain confidentiality when you need to take out a large carrying case. You can carry a smaller-sized FRIO® in your pocket. The confidentiality and convenience are both a welcome difference.

Safety ~ After the hurricanes in New Orleans, there were people with Diabetes that died - because of the lack of electricity. Without electricity, there was no refrigeration. Traditional cases could not be kept cool - and insulin went bad. FRIO® users had no such issues. With a FRIO®, you are no longer dependent on ice packs or refrigeration.*

The FRIO® is a basic necessity for the emergency preparedness of anyone using temperature-sensitive medication such as Insulin, Byetta or Victoza.



Individual

FRIO® Products

Mini ~ 3" x 4"
1 vial

Individual ~ 2.5" x 7"
1 pen

Extra Small/ Eye Drop ~ 5" x 5"
3 vials or eye drop bottles

Small ~ 5.5" x 6"
4 vials

Duo ~ 3.25" x 7"
2 pens or 2 vials

Large ~ 5.5" x 7.5"
Up to 4 pens/vials or combo

Extra Large ~ 6" x 8.25"
Up to 8 pens/vials or combo


Pump ~ 3.75" x 4.75"
Comes w/liner. Holds most pumps.

Other FRIO® Products Available:

- Head/Neck Coolers
- Wrist Bands
- Ankle Bands
- Migra Head Soothers
- Drink Coolers

Available Colors

- Black
- Blue
- Burgundy
- Green
- Purple
- Red



(Head/Wrist/Ankle Bands & Pump not available in Red)

**Insulin to be stored beyond 30 days should be refrigerated
Belts, Pens & Vials not included*

ACTIVATES WITH WATER!

REUSABLE!

www.FRIOCase.com