

## Everyday Hip Chair

The Lumex Everyday Hip Chair is ideal for people who have arthritis or for patients who have undergone hip or knee surgery. The 27" high seat helps to reduce stress and fatigue on joints while sitting and when entering and exiting the chair.

- High density foam in the upholstered seat and back makes sitting more comfortable for longer periods of time
- New adjustable height footrest provides greater comfort while sitting. Footrest height is easily and securely adjusted by two spring-loaded adjustment knobs located on either side of frame
- Padded full-length arms make sitting and standing easier
- The anti-slip strip on the adjustable footrest provides added safety and is an assist during repositioning
- Sturdy metal frame and padded seat and back are easy to clean and maintain
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED

- GF4405427** Blue Ridge 1 ea  
**GF4405857** Jade 1 ea  
**GF4405863** Rosewood 1 ea

**Specifications**

|                              |                           |
|------------------------------|---------------------------|
| Overall Height.....          | 40 ¾"                     |
| Overall Width.....           | 25"                       |
| Overall Depth.....           | 24 ½"                     |
| Width Between Arms.....      | 17 ¾"                     |
| Arm Height From Floor.....   | 32"                       |
| Arm Height From Seat.....    | 5"                        |
| Seat Height.....             | 27 ¼"                     |
| Seat Width.....              | 19"                       |
| Seat Depth.....              | 16"                       |
| Shipping Weight.....         | 52 lb                     |
| Maximum Weight Capacity..... | 350 lb EVENLY DISTRIBUTED |

| Footrest position | From top of footrest to floor | From top of seat to top of footrest |
|-------------------|-------------------------------|-------------------------------------|
| High              | 9 ¼"                          | 18"                                 |
| Mid               | 8 ¼"                          | 19"                                 |
| Low               | 7 ¼"                          | 20"                                 |



Footrest can be adjusted to three different heights.



The AEGIS Microbe Shield used on Lumex Long Term Care and Patient Room Seating works by chemically bonding itself to the treated surface and rupturing the membranes of any bacteria that it comes in contact with which inhibits the growth and spread of microbial growth.

