



Welch Allyn  
**Green Series™ 777**  
Integrated Wall System



**WelchAllyn®**  
Advancing Frontline Care™

# Welch Allyn Green Series™ 777 Integrated Wall System

Allows for reliable diagnostic patient data in a standardized and energy efficient solution to improve your bottom line.



## *Welch Allyn*

*Green Series™ 777 Wall Transformer—configured to each facility's needs.*

- Help reduce facility cost with a more energy efficient product than the 767 Wall Transformer
- Help improve workflow and save energy with handles that automatically turn on/off
- Quickly and easily upgrade from the Welch Allyn 767 Wall Transformer to the Green Series 777 Wall Transformer
- No on/off switch reduces contact, to help reduce the risk of cross-contamination

## *Welch Allyn Otoscopies & Ophthalmoscopes*

- Standardize or choose the technology you need for each facility
- Upgrade to the Welch Allyn LED lamp for greater energy efficiency and lower cost of ownership

### **Welch Allyn MacroView™ Otoscope**

- Almost twice the field of view and 30% greater magnification than traditional otoscopes
- Allows for simplified diagnosis of Otitis Media through superior visualization of the tympanic membrane



### **Welch Allyn Pneumatic Otoscope**

- Rotating lens with ultraseal enables easy pneumatic otoscopy
- Easier instrumentation to allow a faster and more comfortable exam for the patient



### **Welch Allyn PanOptic™ Ophthalmoscope**

- Fast, easy entry into undilated pupils with 5X larger view of the fundus for a streamlined exam
- 25° field of view vs. the standard 5° field of view provides optimal visualization



### **Welch Allyn Coaxial Ophthalmoscope**

- 18 unique aperture/filter combinations for greater versatility
- 28 focusing lenses with a range of -25 to +40 diopters







*Welch Allyn industry-leading blood pressure and thermometry devices provide you with fast, accurate and reliable patient data you can trust—to help improve patient care and clinical decisions.*



#### Welch Allyn Wall Aneroid

- Recessed, laser-engraved dial face with high-contrast pointer provides certified accuracy to reduce errors and improve patient care
- Lifetime calibration warranty reduces cost of ownership



#### Welch Allyn Connex® ProBP™ Digital Blood Pressure Device

- Connectivity to most leading EHR systems, to improve patient documentation accuracy and efficiency
- Improve patient comfort and increase throughput with fast, 15-second NIBP readings
- Portable, so you can take BP readings just about anywhere



#### Welch Allyn Spot Vital Signs® Device

- Improve workflow by using a device that includes NIBP, pulse rate and MAP
- Optional SureTemp® thermometry and Masimo® or Nellcor® pulse oximetry



#### Welch Allyn Spot Vital Signs® LXi Device

- Improve patient comfort and increase throughput with 15-second NIBP readings taken as the cuff inflates
- Connectivity to most leading EHR systems, to improve patient vitals documentation accuracy and efficiency
- Improve workflow by using a device that includes NIBP, pulse rate, MAP, thermometry and SpO<sub>2</sub>



#### Welch Allyn SureTemp® Plus Electronic Thermometer

- Provides accurate oral temperatures in 4 to 6 seconds, axillary and rectal temperatures in 10 to 13 seconds and adult axillary temperatures in 12 to 15 seconds
- Reduce risk of cross-contamination with single-use probe covers
- Improve efficiency with convenient storage for 75 probe covers
- Takes approximately 6,000 readings on 3 AA batteries and provides last temperature recall

*Standardize across facilities for greater efficiency and consistent data capture or customize for your facility's needs.*

| Part #  | WALL BOARD | GS 777 WALL TRANSFORMER | COAXIAL OPHTHALMOSCOPES | PANOPTIC™ OPHTHALMOSCOPES | DIAGNOSTIC MACROVIEW™ OTOSCOPES | PNEUMATIC OTOSCOPE | KLEENSPEC® DISPENSERS | WALL ANEROID | CONNEX® PROBP™ DEVICE | SURETEMP® PLUS THERMOMETER |
|---|------------|-------------------------|-------------------------|---------------------------|---------------------------------|--------------------|-----------------------|--------------|-----------------------|----------------------------|
| Part #  | 77710      | 11720                   | 11810                   | 11820                     | 23810                           | 23820              | 52400-PF              | 52100-PF     | 34XFVTFB              | 01690-300                  |
| <b>GS 777 Integrated Wall Systems</b>   |            |                         |                         |                           |                                 |                    |                       |              |                       |                            |
| 77798   | •          | •                       |                         |                           |                                 | •                  |                       | •            |                       |                            |
| 77791-1MPX  | •          |                         | •                       |                           | •                               |                    | •                     | •            |                       | •                          |
| 77791-2MPX  | •          |                         |                         | •                         | •                               |                    | •                     | •            |                       | •                          |
| 77791-MX  | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77792-1MP   | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77792-2MP   | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77792-M   | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77795-2MPX  | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77795-MX  | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77796-M   | •          |                         |                         |                           | •                               |                    | •                     | •            |                       | •                          |
| 77797-2PX   | •          |                         |                         |                           | •                               |                    |                       | •            |                       | •                          |
| 77797-X   | •          |                         |                         |                           |                                 |                    |                       | •            |                       | •                          |
| 77791-MXNOBP  | •          |                         |                         |                           | •                               |                    |                       |              |                       | •                          |
| 77792-MNOBP   | •          |                         |                         |                           | •                               |                    |                       |              |                       | •                          |
| <b>GS 777 Integrated Wall Systems with Connex® ProBP™ Digital Blood Pressure Device (includes wall mounting bracket and device)</b>                                     |            |                         |                         |                           |                                 |                    |                       |              |                       |                            |
| 77791-MXPROBP   | •          |                         |                         |                           | •                               |                    |                       |              | •                     | •                          |
| 77791-2MPXPROBP   | •          |                         |                         |                           | •                               |                    |                       |              | •                     | •                          |
| <b>GS 777 Integrated Wall Systems with Spot Vital Signs® or Spot Vital Signs® LXI Devices (includes wall mounting bracket only, devices must be ordered separately)</b> |            |                         |                         |                           |                                 |                    |                       |              |                       |                            |
| 77794-MSPOT   | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| 77794-MLXI  | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| 77794-2MSPOT  | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| 77794-2MPLXI  | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| <b>GS 777 Wall Transformer Sets</b>   |            |                         |                         |                           |                                 |                    |                       |              |                       |                            |
| 77510   | •          |                         |                         |                           |                                 |                    |                       |              |                       |                            |
| 77910   | •          |                         |                         |                           |                                 |                    |                       |              | •                     |                            |
| 77710-71M   | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| 77710-81M   | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| 77710-82M   | •          |                         |                         |                           | •                               |                    |                       |              |                       |                            |
| 77710-TAKX  | •          |                         |                         |                           |                                 |                    |                       |              | •                     |                            |

The United States Green Building Councils Leadership in Energy and Environmental Design (LEED) is an internationally-recognized green building certification system providing benchmarks for the design, construction and operation of high-performance green buildings.

Welch Alllyn is proud to have earned a Gold rating in the LEED certification program.



To place an order, call our Customer Service Department at: 1.800.535.6663 or contact your local Welch Alllyn representative.

© 2012 Welch Alllyn SM4085 Rev A www.welchallyn.com



# Welch Allyn Green Series™ 777 Wall Transformer and Integrated Wall System

## Cost Savings Clinical Evidence

### The Cost Savings is Clear

Healthcare organizations spend more than \$8.5 billion on energy each year to meet patient needs so increasing energy efficiency of frequently used medical products has a significant potential for cost savings<sup>1</sup>. Welch Allyn is committed to supporting healthcare organizations achieve their sustainability and cost savings goals by developing the Welch Allyn Green Series™ 777 Wall Transformer and Integrated Wall System (IWS) which incorporate the ENERGY STAR certified power supply.

The Welch Allyn Green Series 777 Wall Transformer and IWS have the benefit of 54% greater energy efficiency over the 767 Wall Transformer and IWS helping decrease your organization's energy usage and costs over the life of the product while continuing to provide high quality Welch Allyn diagnostic products.

The elimination of the power switch also helps your facility avoid the time and costs associated with manually powering on and off the Wall Transformer when not in use to decrease energy consumption. The table below reflects the estimated energy cost savings of the Welch Allyn Green Series 777 Wall Transformer and IWS over the Welch Allyn 767 Wall Transformer and IWS.



### Estimated Cost Savings for Welch Allyn Green Series 777 Wall Transformer & Integrated Wall System

| Item   | Number                |
|--|-----------------------|
| Annual energy savings of 777 vs 767 Wall Transformer or IWS    | 73.53 Kilowatt Hours  |
| Energy cost  | \$0.12 /Kilowatt Hour |
| Estimated 10 year Cost Savings per 777 Wall Transformer or IWS | \$122.00              |
| Number of 777 Wall Transformer or IWS Purchased                | 20 Units              |
| <b>Total Estimated 10 Year Cost Savings</b>                    | <b>\$2,440.00</b>     |

*Assumptions: Energy cost is rising consistently at a rate of 7% per year | 8 hour usage per day/261 days/year/10 years*

<sup>1</sup> EIA, CBCECS 2003, adjusted for inflation to 2009 dollars

**WelchAllyn**<sup>®</sup>

Advancing Frontline Care™