

ADDITIONAL OFFERINGS FEATURING A REUSABLE INNER CANNULA

Shiley™ Reusable Inner Cannula Tracheostomy Tubes

An economical alternative to disposable tubes, these tubes offer:

- An integral 15 mm twist-lock connector
- Reusable inner cannula; may be cleaned and reused up to 29 days on a single patient†
- Four configurations for ventilatory support and weaning: cuffed, cuffless, fenestrated, or nonfenestrated
- Swivel flange that moves with the patient
- High-volume, low-pressure cuff that provides a gentle and effective seal
- Jackson sizes 4, 6, 8, and 10

Latex free. Packaged individually, sterile.

† Refer to directions for use for cleaning instructions.

