

## Pal Seat

Colourful, comfortable, supportive seating system for children aged 1 - 12 years.

### Technical Details

Size	1	2	3	4
Age	1 - 3	2 - 5	4 - 8	7 - 12
Min User Height	700mm (27.6")	900mm (35.4")	1000mm (39.4")	1200mm (47")
Max User Height	950mm (37.4")	1050mm (41.3")	1300mm (51.2")	1600mm (63")
Max User Weight	30kgs (66lbs)	35kgs (77lbs)	40kgs (88lbs)	50kgs (110lbs)

### Product Features

Designed for children aged 1 - 12 years with mild postural needs.

Available in four bright colours, it gives improved stability which helps reduce fatigue, allowing the child longer periods of concentration and fine motor activity.

PAL: Posture+Attention=Learning.

Two-point lap belt and side pads for lateral stability. An optional ramped cushion can be added if children need a little extra help to keep pelvis and femurs in alignment.

Where a less upright angle is required, the backrest angle and armrests adjust to suit the child's needs, and along with the large tray or table, provide appropriate posterior and anterior support.

In the classroom environment, the matching table provides a large surface for anterior support, feeding, writing and other table top activities.

For small children, the optional high seat frame converts Pal size 1 to a high chair, ideal for feeding and inclusion at family meal times.

The seat depth is adjustable for maximum femoral support, and when children need additional help to maintain the alignment of their femurs, a pommel is available.

Lower leg and foot support is provided by the height adjustable legs and footplate, contributing to pelvic stability.

The armrests are height-adjustable.



Available in Green & Blue



## Age Ranges

PRODUCT NAME	SIZE	AGE OF USER																			
		0	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	18+
KIT SEAT	1																				
	2																				
MYGO SEAT	1																				
	2																				
SQUIGGLES SEAT																					
PAL	1																				
	2																				
	3																				
EVERYDAY ACTIVITY SEAT	1																				
	2																				
	3																				

## Feature Comparison

PRODUCT NAME	SIZE	MAX USER WEIGHT		GROSS MOTOR FUNCTION CLASSIFICATION SYSTEM <sup>1</sup>	CHAILEY LEVELS OF ABILITY <sup>2</sup>	4 POINT PELVIC HARNESS	PELVIC CRADLE	2 POINT PELVIC BELT	MULTI-ADJUSTABLE 3 PART BACKREST	SHOULDER PROTRACTION	COMPLEX THORACIC LATERAL SUPPORTS	STANDARD THORACIC LATERALS	HEAD SUPPORTS	CHEST HARNESS	INDEPENDENTLY ADJUSTABLE LUMBAR SUPPORT	ACCOMMODATE LEG LENGTH DISCREPANCY	ACCOMMODATE WINDSWEEPING	SPLIT SEAT BASE	ADDITION GUIDES / HIP PADS	MOUNTING FOR COMMUNICATION DEVICES	ABDUCTION GUIDES / POMMEL	MULTI ADJUSTABLE FOOTPLATES	INDEPENDENTLY ADJUSTABLE FOOTPLATES	STANDARD FOOTPLATE	KNEEBLOCKS / KNEEPADS	LONG SITTING	PRONE SITTING	MANUAL HEIGHT-ADJUSTABLE TILT IN SPACE BASE	MANUAL HI-LOW BASE WITH TILT IN SPACE	MACHINE WASHABLE COVERS	WIPE CLEAN COVERS	INFECTION CONTROL COVERS		
		kg	lbs																															
KIT SEAT	1	75	165	IV-V	1-3	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●									●	
	2	75	165	IV-V	1-2	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●									
MYGO SEAT	1	50	110	IV-V	1-3	●	●			●		●	●	●	●	●	●	●	●	●	●	●	●	●										
	2	60	132	IV-V	1-3	●	●			●		●	●	●	●	●	●	●	●	●	●	●	●	●	●									
SQUIGGLES SEAT		22	48	III-V	1-3	●					●	●	●	●					●	●	●	●	●	●										
PAL	1	30	66	II-III	4-7			●				●	●	●	●				●		●			●									●	
	2	35	77	II-III	4-7			●				●	●	●	●				●		●			●									●	
	3	40	88	II-III	4-7			●				●	●	●	●				●		●			●									●	
EVERYDAY ACTIVITY SEAT	1	30	66	III-IV	3-6	●				●		●	●	●	●				●		●	●												●
	2	50	110	III-IV	3-6	●				●		●	●	●	●				●		●	●												●
	3	100	220	III-IV	3-6	●				●		●	●	●	●				●		●	●												●