

Leglifter

The inflatable Leglifter raises, lifts and transfers legs into bed.

Fits on profiling beds



Technical Specifications

Max user weight	135 kg / 21 stone
Min height	Flat
Total length	124 cm / 48 ¾"
Width	57 cm / 22"
Total Product weight	1.5 kg / 3 ¼ lbs
Airflo Compressor	Mk3

Product code: MPCA140200

Easy to use

- Fits most beds including profiling beds
- Straight forward to fit or remove, and requires no permanent installation
- Attaches to the bed underneath the mattress
- One airline connects the Leglifter to the Airflo Compressor and inflation is activated by pushing one button

Safe to use

- Raises, lifts and transfers legs into the bed
- Leg break on profiling beds can still be used once leglifter is fitted.

Simple to maintain

- For a routine clean, sponge clean with a detergent solution and dry thoroughly.
- The covers can be machine washed
- Spare covers can be purchased

Improves wellbeing

- Minimises moving and handling so reducing the risk of musculoskeletal injury for the carer
- Fits discreetly against the side of the bed
- No permanent installation required

