


ELECTRO  
THERAPY  
UNITS

# MAX TENS

Transcutaneous Electrical Nerve Stimulation Unit

TENS UNITS



## Advantages of MAXTENS

### Fully adjustable

All parameters(Amplitude, Pulse Rate, Pulse Width, Timer, Mode) can be adjustable.

### Convenience

Easy operation, Easy setting, Easy control.

### Sturdy construction with an attractive appearance

### Protective cover to guard against accidental setting changed

Window cover for preventing accidental output adjustment and Slide cover for preventing accidental parameter setting changed.



# MAXTENS

## MAXTENS 1000

Type	Analog
Dimension	1" x 2 3/4" x 4 1/10"
Channel	Dual
Amplitude	Adjustable, 0-80mA peak into 500 ohm load each channel
Pulse Rate	Adjustable, 2 ~ 150 Hz
Pulse Width	Adjustable, 30 ~ 260 $\mu$ s
Timer	15, 30, Continue
Mode	Burst, Modulation, Normal
Power Supply	9V battery type 6F 22



## MAXTENS 2000

Type	Digital
Dimension	1" x 2 3/4" x 4 1/10"
Channel	Dual
Amplitude	Adjustable, 0-80mA peak into 500 ohm load each channel
Pulse Rate	Adjustable, 1 ~ 150 Hz
Pulse Width	Adjustable, 30 ~ 260 $\mu$ s
Timer	15, 30, 45, 60, 90 Continuous
Mode	Burst, Modulation, Normal
Power Supply	9V battery type 6F 22
LCD Monitoring	



### Product Configuration



TENS Unit X 1  
 TENS Electrodes X 4  
 Leadwires X 2  
 9V Battery, type 6F22 X 1  
 Instruction Manual X 1  
 Carrying Case X 1