

## Skin Closure Strips

McKesson skin closure strips are designed to provide superior wound closure and security. They are made of a porous non-woven polyamide material with pressure-sensitive adhesive. McKesson skin closure strips are ideal for primary closure, and help secure, close and support small cuts and wounds. Additionally, they provide wound support following early suture or staple removal. This product is available as non-reinforced or reinforced formulations, and comes in a variety of sizes.

### Features and Benefits

- Ideal for closing small cuts and wounds
- Unique porous construction is designed to permit skin ventilation and help prevent maceration
- Easy to apply and comfortable to wear
- Flexible strips are tan and reinforced strips are white
- Not made with natural rubber latex



### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them — making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

## ORDERING INFORMATION

MFR #	DESCRIPTION	QTY
3000	Skin Closure Strips, Non-Reinforced, 1/8" x 3", Sterile	5ea/pk, 50pk/bx, 4bx/cs
3001	Skin Closure Strips, Non-Reinforced, 1/8" x 3", Sterile	3ea/pk, 50pk/bx, 4bx/cs
3002	Skin Closure Strips, Non-Reinforced, 1/4" x 1 1/2", Sterile	6ea/pk, 50pk/bx, 4bx/cs
3003	Skin Closure Strips, Non-Reinforced, 1/4" x 4", Sterile	10ea/pk, 50pk/bx, 4bx/cs
3004	Skin Closure Strips, Non-Reinforced, 1/2" x 4", Sterile	6ea/pk, 50pk/bx, 4bx/cs
3005	Skin Closure Strips, Non-Reinforced, 1" x 5", Sterile	4ea/pk, 25pk/bx, 4bx/cs
3006	Skin Closure Strips, Reinforced, 1/8" x 3", Sterile	5ea/pk, 50pk/bx, 4bx/cs
3007	Skin Closure Strips, Reinforced, 1/4" x 3", Sterile	3ea/pk, 50pk/bx, 4bx/cs
3008	Skin Closure Strips, Reinforced, 1/4" x 1 1/2", Sterile	6ea/pk, 50pk/bx, 4bx/cs
3009	Skin Closure Strips, Reinforced, 1/4" x 4", Sterile	10ea/pk, 50pk/bx, 4bx/cs
3010	Skin Closure Strips, Reinforced, 1/2" x 4", Sterile	6ea/pk, 50pk/bx, 4bx/cs
3011	Skin Closure Strips, Reinforced, 1" x 5", Sterile	4ea/pk, 25pk/bx, 4bx/cs