

McKesson ECG Electrodes



Radio-Translucent Electrodes



Medi-Pak™ Performance Radio-Translucent Monitoring Foam Electrodes, 1-7/8"

- Silver/silver chloride sensor for quality traces
- Diaphoretic aggressive foam RT snap with solid adhesive gel
- For cath lab, x-ray and thallium procedures
- Easy removal, no gel residue or clean up

| Quantity | Reorder No. | EC No. |
|------------------------|-------------|----------|
| 5/card, 10 cards/pouch | 87-310-RT | 31012500 |



Medi-Pak™ Performance Radio-Translucent Monitoring Tape Electrodes, 1-7/8"

- Silver/silver chloride sensor for quality traces
- Diaphoretic breathable tape RT snap with solid adhesive gel
- For cath lab, x-ray and thallium procedures
- Easy removal, no gel residue or clean up

| Quantity | Reorder No. | EC No. |
|------------------------|-------------|----------|
| 5/card, 10 cards/pouch | 87-320-RT | 31022500 |

The product information contained in this catalog, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers, not McKesson Medical-Surgical or its affiliates ("McKesson") and have not been independently verified by McKesson. McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the catalog, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print catalogs by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this catalog is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson Electrodes are hypoallergenic, latex-free and single-use only.

McKesson also carries a complete line of recording chart paper and shave prep razors to meet your needs. Ask your McKesson Account Manager for more information.

Resting ECG Electrodes



Medi-Pak™ Performance Plus Resting ECG Tab Electrodes

- Larger surface area contact
- Silver/silver chloride and conductive adhesive gel for superior adhesion and quality traces
- Pliable backing conforms and is repositionable
- Easy removal, no gel residue or clean up
- Blue card allows for easier visibility of the liner and is designed to prevent accidental slipping if dropped
- Defib recoverable

| Quantity | Reorder No. | EC No. |
|------------------------------------|-------------|----------|
| 10/card, 10 cards/pouch, 10 pch/bx | 87-230 | 72302500 |



Medi-Pak™ Performance Resting ECG Tab Electrodes

- Specially designed strain relief tab
- Silver/silver chloride and conductive adhesive gel for superior adhesion and quality traces
- Pliable backing conforms and is repositionable
- Easy removal, no gel residue or clean up
- Blue card allows for easier visibility of the liner and is designed to prevent accidental slipping if dropped
- Defib recoverable

| Quantity | Reorder No. | EC No. |
|------------------------------------|-------------|----------|
| 10/card, 10 cards/pouch, 10 pch/bx | 87-220 | 87222500 |



Medi-Pak™ Performance Resting ECG and Echo Sound Foam Electrodes, 1-3/8"

- Silver/silver chloride sensor for quality traces
- Diaphoretic foam resting ECG snap
- Foam electrodes with 10% chloride wet gel

| Quantity | Reorder No. | EC No. |
|-------------------------|-------------|----------|
| 10/card, 10 cards/pouch | 87-90-100F | 79012500 |

Resting ECG Electrodes (cont.)



Medi-Pak™ Performance Resting ECG and Echo Sound Foam Electrodes, 7/8" x 1-3/4"

- Silver/silver chloride sensor for quality traces
- Diaphoretic foam resting ECG snap
- Foam electrodes with solid adhesive gel
- Easy removal, no gel residue or clean up

| Quantity | Reorder No. | EC No. |
|------------------------|-------------|----------|
| 5/card, 10 cards/pouch | 87-490 | 74902500 |



Medi-Pak™ Performance Resting ECG and Echo Sound Vinyl Electrodes, 1-3/8"

- Silver/silver chloride sensor for quality traces
- Vinyl resting ECG snap
- Vinyl electrodes with 10% chloride wet gel

| Quantity | Reorder No. | EC No. |
|-------------------------|-------------|----------|
| 10/card, 10 cards/pouch | 87-90-100V | 79102500 |

Stress, Holter and General Purpose Monitoring Electrodes



Medi-Pak™ Performance All Purpose Foam Electrodes

- Silver/silver chloride sensor for quality traces
- Diaphoretic aggressive foam snap with solid adhesive gel
- For stress test, long and short term monitoring or resting ECG
- Easy removal, no gel residue or clean up

| Size | Quantity | Reorder No. | EC No. |
|--------|------------------------|-------------|----------|
| 1-7/8" | 5/card, 10 cards/pouch | 87-310 | 73102500 |
| 1-1/2" | 5/card, 10 cards/pouch | 87-50SG | 87502500 |

Stress, Holter and General Purpose Monitoring Electrodes (cont.)



Medi-Pak™ Performance All Purpose Tape Electrodes, 1-7/8"

- Silver/silver chloride sensor for quality traces
- Diaphoretic breathable tape snap with solid adhesive gel
- For stress test, long and short term monitoring or resting ECG
- Easy removal, no gel residue or clean up


| Quantity | Reorder No. | EC No. |
|------------------------|-------------|----------|
| 5/card, 10 cards/pouch | 87-320 | 73202500 |



Medi-Pak™ Performance ECG Diagnostic Electrodes, 1-1/2"

- Silver/silver chloride and 10% chloride wet gel provide fast, clear traces
- Diaphoretic aggressive foam snap with high chloride gel
- For stress, holter and other diagnostic procedures
- Superior adhesion for reliability in stress testing

| Quantity | Reorder No. | EC No. |
|------------------------|-------------|----------|
| 5/card, 10 cards/pouch | 87-50 | 75002500 |



Medi-Pak™ Performance Stress Test Kits

- Silver/silver chloride
- Diaphoretic aggressive foam snap with 10% chloride wet gel
- Superior adhesion for reliability in stress testing
- 12-lead ECG kit with 10 electrodes included
- Each electrode designed for a specific location

| Quantity | Reorder No. | EC No. |
|----------|-------------|----------|
| 10/kit | 87-150-10S | 71502500 |



Medi-Pak™ Performance EMS 10-Pack Electrodes

- Silver/silver chloride
- Diaphoretic aggressive foam snap with 10% chloride wet gel
- Superior adhesion for reliability in stress testing
- 12-lead ECG kit with 10 electrodes included
- Each electrode designed for a specific location
- Defib recoverable

| Quantity | Reorder No. | EC No. |
|----------|-------------|----------|
| 10/kit | 87-520-10S | 75202500 |

Our philosophy

We help manage the business of healthcare, so you can manage the healthcare of patients.

McKesson Medical-Surgical

8741 Landmark Road
Richmond, VA 23228

www.mckesson.com
1.800.446.3008