

McKESSON



Sharps Containers

McKesson offers a full line of sharps containers and accessories that feature OSHA-compliant safety features and excellent performance. The line was created to help protect our customers by limiting their exposure to sharps and infectious waste.

We offer a full line of containers in a range of sizes—from 1 quart to 18 gallon—to accommodate different volumes, sizes and types of waste. Our large PG-11 containers comply with D.O.T. requirements to eliminate the need for secondary packaging, and are made of puncture- and leak-resistant plastic. They also feature a gasketed hinged lid. Our standard and multiuse sharps containers have OSHA-compliant design features and are suitable for use in a variety of healthcare settings.

Our containers are stackable to help save valuable storage space. They fit in Prevent® Wall Cabinets, but are also compatible with existing cabinets and brackets offered by leading manufacturers. Our sharps containers are not made with natural rubber latex, and each product in the line is designed to provide our customers with a combination of quality, performance and value.



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2016 McKesson Medical-Surgical Inc.

McKESSON

PREVENT®

Large PG-II Containers

McKesson Prevent® Premium Sharps Containers are rated D.O.T. PG-II, and require no secondary packaging when transporting regulated medical waste off-site for disposal. They allow for touchless, point-of-use, single-handed deposit of used sharps and biohazardous waste. Our large, premium containers are tamper-proof and have a sealing glue-gasket lid for permanent closure that resists leakage of fluids, as well as fully integrated side-handles for easy carrying of full, heavy containers.

The hermetically-sealed, duo-hinged-lid produces an audible click when locked correctly. It is designed to accommodate large volumes of waste of various types and sizes, while keeping the lid connected to the container to reduce the number of steps in the disposal process to save staff time while promoting OSHA-compliant sharps disposal best practices. These features make our large, premium containers ideal for the waste haulers as well as hospitals, labs, surgery centers, pharmacies and outpatient treatment centers.

FEATURES AND BENEFITS

Premium Biohazard Sharps Containers

- Made of puncture-resistant plastic
- Allows for touchless, one-hand deposit of used sharps
- Containers are offered in two colors (red and translucent) and three sizes (8-, 12- and 18-gallon) to accommodate various volumes, sizes, and types of sharps waste
- Compatible with existing wall cabinets and brackets for security and unwanted access prevention



MFR # 2266

Premium Chemotherapy Sharps Containers

- Made of solid, durable, plastic to provide high shock and leakage resistance
- Dual-opening, hinge-top lid; allows small or large side to remain open while other is locked
- Accommodates different types and sizes of trace chemotherapy waste and sharps
- Three sizes (8-, 12- and 18-gallon) available
- Color: yellow (with white lid)



MFR # 2258

MCKESSON



Premium Pharmaceutical Waste Containers

- Made of solid, durable, plastic to provide high shock and leakage resistance
- Dual-opening, hinge-top lid; accepts larger, bulkier, and liquid waste
- Ideal for disposal of non-hazardous pharmaceutical waste
- Closes temporarily between uses; locks for final disposal to help deter tampering
- Color: white (with blue lid)
- Easily distinguished from sharps, chemo and other disposal containers



MFR # 2256

Standard Biohazard Infectious Waste Sharps Containers

Made of puncture and leak-resistant plastic, containers are offered in multiple colors and sizes—from 1.7-quart to 8 gallon—to accommodate various volumes, sizes and types of sharps waste. We offer phlebotomy as well as multiuse and patient-room sharps containers. Most containers are nestable to maximize valuable storage space, and have many OSHA-compliant design features such as touchless deposit, visible fill-level indicators, and integrated handles for safe transportation. Our standard sharps containers are compatible with most existing wall cabinets and brackets for security and unwanted access prevention.

FEATURES AND BENEFITS

- Puncture resistant; strong polypropylene containers meet or exceed industry standards
- Horizontal drops allow full utilization of capacity to fill-line
- Locking lid helps ensure safe transport
- Translucent lid for quick visual identification of fill level
- 1-, 1.7- and 5-quart containers feature needle rewind capability
- Horizontal drop helps to maximize capacity
- Clear top allows for visual monitoring
- Stackable for compact storage prior to use
- For use on countertops or mounted to wall



MFR # 2269

McKESSON


PREVENT®

Sharps Container Accessories

McKesson Prevent® also offers a line of sharps container accessories, including wall brackets, locking cabinets and glove box holders, for user convenience and to help facilitate sharps disposal best practices in a variety of health care settings.

FEATURES AND BENEFITS

- Wall cabinets keep containers secured in locked enclosure
- Economical wall mounting brackets provide an unobtrusive method for mounting containers
- Ensure compliance; mount sharps collectors at OSHA recommended heights and locations
- Optional glove box holder keeps exam gloves at the ready
- Neutral cream color; suitable for any decor
- Keeps counter space clear
- Easy to clean and disinfect



MFR # 2240

DESCRIPTION	MFR#	QTY
PREMIUM BIOHAZARD SHARPS CONTAINERS		
8 Gallon, White Port Lid, 33 lbs. Weight Bearing Capacity, Pg-II, 13 L x 17.3 W x 13.5 H	2266	1ea, 9ea/cs
12 Gallon, White Port Lid, 55 lbs. Weight Bearing Capacity, Pg-II, 13 L x 17.3 W x 20.8 H	2267	1ea, 8ea/cs
18 Gallon, White Port Lid, 62 lbs. Weight Bearing Capacity, Pg-II, 13 L x 17.3 W x 24.68 H	2268	1ea, 7ea/cs
PREMIUM CHEMOTHERAPY SHARPS CONTAINERS		
8 Gallon, White Port Lid, 33 lbs. Weight Bearing Capacity, Pg-II, 13" L x 17.3" W x 13.5" H	2258	1ea, 9ea/cs
12 Gallon, White Port Lid, 52 lbs. Weight Bearing Capacity, Pg-II, 13" L x 17.3" W x 20.8" H	2259	1ea, 8ea/cs
18 Gallon, White Port Lid, 62 lbs. Weight Bearing Capacity, Pg-II, 13" L x 17.3" W x 24.68" H	2260	1ea, 7ea/cs
PREMIUM PHARMACEUTICAL WASTE CONTAINERS		
8 Gallon, White, Blue Port Lid, 33 lbs. Weight Bearing Capacity, Pg-II (Non-Hazardous), 13" L x 17.3" W x 13.5" H	2256	1ea, 9ea/cs
18 Gallon, White, Blue Port Lid, 62 lbs. Weight Bearing Capacity, Pg-II (Non-Hazardous), 13" L x 17.3" W x 24.68" H	2257	1ea, 7ea/cs

MCKESSON**DESCRIPTION****MFR#****QTY****STANDARD BIOHAZARD SHARPS CONTAINERS**

1.7 Quart, Red, Tall Tray	101-185	1ea, 20ea/cs
1 Gallon, Red	101-8703	1ea, 24ea/cs
3 Gallon, Red	101-8710	1ea, 12ea/cs
1 Quart, Red, 6.25" L x 4.25" W x 4.25" H	065	1ea, 80ea/cs
2 Gallon, Red, 10.25" L x 7" W x 10.5" H	047	1ea, 20ea/cs
8 Gallon, Red, 13.75" W x 13.75" D x 14" H	80-8705	1ea, 10ea/cs
5 Quart, Translucent, 10.75" x 10-1/5" x 4.75"	2261	1ea, 10ea/bx, 20bx/cs
5 Quart, Red, 10.75" x 10 1/5" x 4.75"	2262	1ea, 10ea/bx, 2bx/cs
5.4 Quart, Red, 11" x 12" x 4.75"	2269	1ea, 10ea/bx, 2bx/cs
5.4 Quart, Translucent, 11" x 12" x 4.75"	2270	1ea, 10ea/bx, 2bx/cs

ACCESSORIES

Sharps Wall Locking Cabinet, Cream, 5 Quart, 10" x 12.25" x 4.75", For Use With # 2261 and #2262	2263	1ea, 2ea/cs
Glove Box Holder, Cream, For Use With # 2263, 11" W x 6.5" H x 3-7/8" D	2264	1ea, 2ea/cs
Sharps Kit: Includes: #2263 (Wall Cabinet), #2265 (Glove Box Holder) and 5 Quart Sharps Containers (#2261/2262)	2265	1/cs
Sharps Wall Bracket, 5.4 Quart w/ Locking Swing Arm, 9.36" H x 10" W x 4.35" D, For Use With #2269 and #2270	2240	2/cs