



**AFFORDABLE STYLE
AND QUALITY.**

Medline Overbed Tables
A Complete Line-Up for Every Application



MATCH ANY ROOM. MEET EVERY NEED.

Perfect for comfortable dining, reading and other activities, Medline overbed tables are available in multiple finishes to match your facility's room case goods and offer a variety of features to meet different needs:

- » Extra space with a flip-up table for convenient eating
- » High weight capacity—no worrying about overbed table failure
- » Low overbed tables for low beds

ECONOMY HOME CARE OVERBED TABLES



COMPOSITE-TOP OVERBED TABLE

This durable table will not crack or peel.

- » Easy to clean
- » Impervious to moisture
- » Height adjustable with one hand, from 30–47" (76–119 cm)

Item Number	Description
MDS107015	Composite-Top Overbed Table



LAMINATE-TOP OVERBED TABLE

This table has a walnut, wood-grain and melamine laminate with a "T" mold edge.

- » Automatic spring-lift mechanism adjusts height from 28–45" (71–114 cm)

Item Number	Description
MDS104015	Walnut-Top Overbed Table

FINISH OPTIONS



LAMINATE TILT-TOP OVERBED TABLE

This table has separate handles for tilt, and raise and lower functions.

- » Three-position tilt with squeeze of the handle
- » Melamine laminate with a "T" mold edge
- » Height adjustment: 28–45" (71–114 cm)

Item Number	Description
MDS104950	Walnut Tilt-Top Overbed Table

FINISH OPTIONS

