

ML300



The ML300 is a flanged mounted pool lift that is ideal for a wide range of pool applications.

- Third party tested & verified ADA compliant
- Integrated armrests
- State of California compliant
- 300 lb/136kg lifting capacity
- Retrofit anchor jig is standard
- Flange mounted with 4 anchor bolts
- LiftOperator® Intelligent Controller
- Powder-coated stainless steel and aluminum construction



Lift Color



GRAY MIST

Due to printing technology actual color may differ.



New Construction Jig with Anchors
500-5000A

Model No.	Description	Shipping – (Class 85)			
		Weight	Length	Width	Height
575-3000	ML300	200 lbs 90kg	50" 127cm	24" 61cm	28" 71cm
575-3000N	ML300, no anchor				

Note:


- Includes battery, charger, battery console cover, water-resistant hand control, footrest, seat belt assembly, armrests and retrofit anchor jig.

Parts & Accessories

- 1001495 Battery
- 300-6900 ML300 Anchor Kit for existing deck
- 500-5200 Cover
- 300-6700A Anchors, set of 4
- 300-6800A Anchor Bolts, set of 4
- 900-1000 Seat Belt
- 970-5000 Seat Saver Cover
- 900-6000 Stability Strap
- 170-3000A Armrest Replacement (pair)

New Construction Guidelines

Scenario	Pool Lift	Anchor Jig
Order pool lift and new construction jig at same time	575-3000	500-5000 (no anchors)
Order new construction jig ahead of pool lift	575-3000N	500-5000A (comes with anchors)



A completed Deck Profile form is required with your pool lift order.

