



Introducing the new Moxie GT Gait Trainer!

The Moxie GT offers a tremendous amount of versatility to meet the needs of any user. Height adjustable and available in three sizes. The Moxie GT can be used in either a posterior and anterior position.



FEATURES

- Height adjustable in 1" increments
- Medium and Large frames fold for easier storage and transportation
- Swivel casters can be locked in non-swivel position
- Directional lever allows for both forward and reverse mobility, or forward mobility only
- Wheel locks on all four casters to keep gait trainer stationary during transfer
- Accessories install and adjust without tools

FRAME



Warrior Blue




Magic Green



Sword Gray
(available in Large Only)

ACCESSORIES



Hand Loops



Forearm Platforms




Handgrips



Soft Seat Harness



Trunk Support



Hip Positioner & Pad