

Folding 3-In-1 Commodes

- » Folding frame design easily opens and folds, and will fold flat for convenient storage and transportation
- » Clip-on seat can be removed for easy cleaning
- » Includes pail, lid and splashguard
- » Elongated seat commode has a deeper seat to allow more room for cleaning

Item No.	Description	Weight Capacity	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
G30213-4F	Basic, Steel	350 lbs	16.5–22.5"	16"	19"	21.25"	4/cs
G30210-4F	Premium, Aluminum	250 lbs	17.5–21.5"	16"	19"	21.25"	4/cs
MDS89664EL	Elongated Seat, Steel	400 lbs	17–23"	16"	19"	22"	4/cs



Non-Folding Commodes

- » Use in the bathroom as an elevated toilet seat or as toilet safety rails when combined with most standard toilets
- » Includes pail, lid and splashguard
- » Clip-on seat can be removed easily for cleaning
- » Gray powder-coated steel commode is rust-resistant and durable
- » Bariatric commodes are extra-wide for larger patients

Item No.	Description	Weight Capacity	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS89664	Basic, Steel	350 lbs	15–21"	15"	19"	21.25"	4/cs
G30216B	Bariatric, Steel	500 lbs	16.5–23.5"	18.5"	23.5"	27"	2/cs
MDS89664XW	Extra-Wide Bariatric, Steel	650 lbs	18–23.5"	15"	24"	28"	ea



Drop-Arm Commodes

- » Arms swing out of the way for safe lateral transfers on and off the commode
- » Durable, rust-resistant, powder coated reinforced frame
- » Includes pail, lids and splash guard

Item No.	Description	Weight Capacity	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Pkg.
MDS89668	Basic, Steel	250 lbs	17–21"	15"	17.5"	25.5"	ea
G98202	Standard Non-Padded, Steel	300 lbs	18.25–23.25"	18–20"	18.5"	22"	ea
G98204	Premium Padded, Steel	350 lbs	20–25"	18–20"	18.5"	22"	ea
MDS89668XW	Bariatric, Steel	850 lbs	20–25"	19"	28"	29"	ea

