


AirLife[®] oxygen therapy consumables product portfolio



AirLife® oxygen therapy consumables product portfolio

Our extensive AirLife portfolio of oxygen therapy solutions empowers you to provide the highest quality patient care. From the most basic nasal cannula to the clinically differentiated OxyMask™ from Southmedic, AirLife meets the needs of oxygen care patients both in the hospital and at home.

Why AirLife?

- Clinical differentiation through an industry-leading and unique product portfolio
- Standardization around a core respiratory portfolio with a leading set of products
- Quality care solutions that reduce variability in processes and products

Products for oxygen therapy patients*

** Not a comprehensive list; please contact your AirLife representative or refer to our ecatalog [carefusion.com/ALcatalog](https://www.carefusion.com/ALcatalog) if you seek a product not listed here.*

Table of contents

AirLife cannulas.....	1-2
Soft cushion cannulas.....	1
Foam cover cushion cannulas.....	1
Cushion cannulas.....	2
Standard nasal cannulas.....	2
AirLife masks.....	3-4
Disposable O ₂ masks.....	3
Venturi masks.....	4
Aerosol masks.....	4
Tracheostomy masks.....	4
OxyMask™ from Southmedic.....	5-6
OxyMask.....	5
OxyMulti-Mask™.....	6
OxyMask aerosol.....	6
OxyTrach™.....	6
Bubble humidifiers and large volume nebulizers.....	7-8
Bubble humidifier.....	7
Large volume nebulizer.....	7
Misty Ox™.....	8
Oxygen tubing.....	8

All AirLife cannulas come with crush-resistant tubing to minimize occlusion and kinking. Select cannulas come with U/Connect-It® oxygen tubing for easy connection to any oxygen source.

Soft cushion cannulas

AirLife soft cushion cannulas include Iariat tubing that is approximately 37% softer than our standard AirLife cannulas.** This soft tubing helps reduce chafing behind the ear to minimize skin irritation. They are made without DEHP.

Soft cushion cannula		
Catalog no.	Description	Qty. (cs.)
SFT2611	Neonate, 7' O ₂ tubing	50
SFT2601	Infant, 7' O ₂ tubing	50
SFT2602	Pediatric, 7' O ₂ tubing	50
SFT2692	Pediatric, 7' O ₂ U/Connect-It tubing	50
SFT2600	Adult, 7' O ₂ tubing	50
SFT2699	Adult, 7' O ₂ U/Connect-It tubing	50
SFT2614	Adult, 14' O ₂ tubing	50
SFT2614U	Adult, 14' O ₂ U/Connect-It tubing	50
SFT2625	Adult, 25' O ₂ tubing	25



Foam cover cushion cannulas

AirLife cushion cannulas with soft foam ear covers help reduce skin irritation around the ears. The soft, gray covers slide smoothly on the headset tubing, making them easy to position for enhanced comfort and fit. They are made without DEHP.

Foam cover cushion cannula		
Catalog no.	Description	Qty. (cs.)
FM2602	Pediatric, 7' O ₂ tubing	50
FM2692	Pediatric, 7' O ₂ U/Connect-It tubing	50
FM2600	Adult, 7' O ₂ tubing	50
FM2699	Adult, 7' O ₂ U/Connect-It tubing	50
FM2614	Adult, 14' O ₂ tubing	50
FM2625	Adult, 25' O ₂ tubing	25
002016	Cannula ear covers (<i>by pair</i>)	50




** Based on a material durometer difference of 30 using the Shore A hardness scale

Cushion cannulas

AirLife cushion cannulas have a soft, clear cushion design that supports the nasal prongs to optimize patient comfort. Cushion cannulas are available in sizes from neonate to adult with a variety of oxygen tubing lengths. They are made without DEHP.

Cushion cannula		
Catalog no.	Description	Qty. (cs.)
002611	Neonate, 7' O ₂ tubing	50
002601	Infant, 7' O ₂ tubing	50
002602	Pediatric, 7' O ₂ tubing	50
002692	Pediatric, 7' O ₂ U/Connect-It tubing	50
002600	Adult, 7' O ₂ tubing	50
002699	Adult, 7' O ₂ U/Connect-It tubing	50
002600-14	Adult, 14' O ₂ tubing	50
002600-14U	Adult, 14' O ₂ U/Connect-It tubing	50
002600-25	Adult, 25' O ₂ tubing	25
002600-50	Adult, 50' O ₂ tubing	20



Standard nasal cannulas

AirLife standard nasal cannulas offer an angulated, flexible lip plate, which provides comfort at the tip of the nasal region, with beveled tips to provide a softer edge. Disposable cannulas are curved and are available with flared or non-flared tips. They are packaged separately with varying lengths of oxygen tubing.

Standard cannula		
Catalog no.	Description	Qty. (cs.)
001320	Adult, flared tip, 7' O ₂ tubing	50
001365	Adult, flared tip, 7' O ₂ U/Connect-It tubing	50
001321	Adult, flared tip, 14' O ₂ tubing	50
001322	Adult, flared tip, 21' O ₂ tubing	25
001325	Adult, non-flared tip, 7' tubing	50
001366	Adult, non-flared tip, 7' O ₂ U/Connect-It tubing	50
001326	Adult, non-flared tip, 14' O ₂ tubing	50
001326U	Adult, non-flared tip, 14' O ₂ U/Connect-It tubing	50
001327	Adult, non-flared tip, 21' O ₂ tubing	25
001329	Adult, non-flared tip, 50' O ₂ tubing	15



Disposable O₂ masks

AirLife disposable oxygen masks are latex-free and made from soft, high-grade vinyl resins. The smooth, feathered edges are designed for patient comfort while reducing irritation points. AirLife oxygen masks are designed for anatomical compatibility. A full range is available to meet all requirements in infant, pediatric and adult sizes. Select cannulas come with U/Connect-It oxygen tubing for easy connection to any oxygen source. AirLife three-in-one masks are packaged as high-concentration, non-rebreathing masks complete with reservoir bags and one-way valves.



Oxygen mask		
Catalog no.	Description	Qty. (cs.)
2D4823	Infant, medium-concentration, 7' O ₂ tubing	50
001262	Pediatric, medium-concentration, 7' O ₂ tubing	50
001267	Pediatric, high-concentration, rebreather, 7' crush-resistant O ₂ tubing	50
001268	Pediatric, high-concentration, non-rebreather 3-in-1 with safety vent, 7' crush-resistant O ₂ tubing	50
001269	Pediatric, high-concentration, non-rebreather, 7' O ₂ tubing	50
001260	Pediatric, medium-concentration, short style, 7' O ₂ tubing	50
001265	Pediatric, high-concentration, short style, 7' crush-resistant O ₂ tubing	50
001201	Adult, medium-concentration, 7' crush-resistant O ₂ tubing	50
001361	Adult, medium-concentration, 7' crush-resistant O ₂ U/Connect-It tubing	50
001203	Adult, high-concentration, non-rebreather 3-in-1 with safety vent, 7' crush-resistant O ₂ tubing	50
001362	Adult, high-concentration, non-rebreather 3-in-1 with safety vent, 7' crush-resistant O ₂ U/Connect-it tubing	50
001210	Adult, high-concentration, non-rebreather 3-in-1, 7' crush-resistant O ₂ tubing	50
001211	Adult, high-concentration, non-rebreather 3-in-1 with safety vent, 7' O ₂ tubing	50
001212	Adult, high-concentration, non-rebreather 3-in-1, 7' O ₂ tubing	50

Venturi masks

AirLife offers a wide selection of both standard and Dual Dial™ Venturi-style masks that feature either six diluter jets or two oxygen concentration dials. All masks are under-the-chin style and supplied with a six-inch flextube, humidification hood and seven-foot crush-resistant O₂ tubing. We also offer a Venturi-style with a tracheostomy mask for use with a 50 psi oxygen source and flowmeter.

Venturi mask		
Catalog no.	Description	Qty. (cs.)
001250	Pediatric, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ tubing	50
001240	Adult, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ tubing	50
001248	Adult, tracheostomy, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ tubing	50
001363	Adult, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ U/Connect-It tubing	50
001259	Pediatric, 6" flextube, Dual Dial™, humidification hood, 7' crush-resistant O ₂ tubing	50
001255	Adult, 6" flextube, Dual Dial, humidification hood, 7' crush-resistant O ₂ tubing	50



Aerosol masks

These soft, high-quality AirLife disposable aerosol masks are made from high-grade resins. The smooth, feathered edges are designed for patient comfort by reducing irritation points. All AirLife masks are designed for anatomical compatibility. A comprehensive line is available to meet requirements in pediatric and adult sizes.

Aerosol mask		
Catalog no.	Description	Qty. (cs.)
001263	Pediatric	50
001261	Pediatric, short style	50
001206	Adult	50
001208	Adult, 6' corrugated tubing	50



Tracheostomy masks

The AirLife disposable tracheostomy masks are available in pediatric and adult sizes, and incorporate several patient features. The neckband is made from a comfortable, nonbiting material, and the swivel tubing connector permits access from either side of the patient. The special strap clips allow the mask to be removed with minimal disturbance to the patient.

Tracheostomy mask		
Catalog no.	Description	Qty. (cs.)
001225	Pediatric	50
001226	Adult	50
001227	Adult, large	50

OxyMask from Southmedic

The OxyMask delivers 24–90% oxygen at flows of one to 15 LPM or to flush. Its unique open design provides a comfortable, non-claustrophobic fit for the patient while facilitating clear speech communication and easy oral access. This helps to ensure uninterrupted therapy while reducing the possibility of rebreathing CO₂. The OxyKid™ and OxyTyke™ are for pediatric and infant patients, respectively. The OxyPlus™ large mask is 25% larger than the standard OxyMask, while the OxyTyke delivers 22–65% oxygen. OxyMask products do not contain natural rubber latex or phthalates such as DEHP.

OxyMask from Southmedic		
Catalog no.	Description	Qty. (cs.)
OM-1125-8	OxyMask, 7' O ₂ tubing	25
OM-1125-14	OxyMask, 14' O ₂ tubing	25
OM-1425-8	OxyMask with ear elastics, 7' O ₂ tubing	25
OM-7025-8	OxyMask, non-PVC, 7' O ₂ tubing	25
OP-1125-8	OxyPlus, large mask, 7' O ₂ tubing	25
OP-1125-14	OxyPlus, large mask, 14' O ₂ tubing	25
OP-1425-8	OxyPlus large mask with ear elastics, 7' O ₂ tubing	25
OP-7025-8	OxyPlus, non-PVC, 7' O ₂ tubing	25
OK-1125-8	OxyKid, pediatric mask, 7' O ₂ tubing	25
OK-1125-14	OxyKid, pediatric mask, 14' O ₂ tubing	25
OK-7025-8	OxyKid, non-PVC, 7' O ₂ tubing	25
OT-1025-8	OxyTyke, pediatric mask, 7' O ₂ tubing	25
OT-7025-8	OxyTyke, non-PVC, 7' O ₂ tubing	25




*AARC Clinical Practice Guideline reprinted from Respiratory Care (*Respir Care* 2002; 47(6): 717-720)

One OxyMask can replace multiple traditional oxygen therapy interfaces to help facilities drive:

- SKU standardization
- Inventory reduction
- Reduction of medical waste
- Reduction of complexity for bedside staff

OxyMask delivers oxygen concentrations from 24–90%* FiO₂, at flows ranging from 1–15+ LPM or to flush.





OxyMulti-Mask

The OxyMulti-Mask is a multipurpose mask with a 22 mm swivel adapter that delivers 23–83% oxygen at flows from one to 15 LPM or to flush. It can be used with or without high humidity or aerosol therapy. It is available in three sizes. The OxyMulti Plus mask is 25% larger than the OxyMulti-Mask. The OxyMulti Kid is for patients three to 10 years of age, 33–70 pounds. This mask does not contain natural rubber latex or phthalates such as DEHP.

OxyMulti-Mask		
Catalog no.	Description	Qty. (cs.)
OHH-1025-8	OxyMulti-Mask with adapter, 7' O ₂ tubing	25
OHH-1425-8	OxyMulti-Mask with ear elastics, adapter, 7' O ₂ tubing	25
OPHH-1025-8	OxyMulti-Mask Plus with adapter, 7' O ₂ tubing	25
OPHH-1425-8	OxyMulti-Mask Plus with ear elastics, adapter, 7' O ₂ tubing	25
OKHH-1025-8	OxyUltra Kid with adapter, 7' O ₂ tubing	25
OKHH-1425-8	OxyUltra Kid with ear elastics, adapter, 7' O ₂ tubing	25

OxyMask aerosol

The OxyMask aerosol mask promotes greater patient-use compliance for aerosol therapy because it is comfortable to wear and has a unique open design. The swivel adapter helps optimize aerosol delivery when the patient is lying down or in a vertical position. The OxyMask aerosol mask is available in four sizes.

The OxyPlus large aerosol mask is 25% larger than the OxyMask aerosol. The OxyKid aerosol mask is for patients three to 10 years of age, 16–33 pounds. The OxyTyke aerosol mask is for patients from six months to three years of age, 10–33 pounds. This mask does not contain natural rubber latex or phthalates such as DEHP.

OxyMask aerosol		
Catalog no.	Description	Qty. (cs.)
OMN-1025	OxyMask aerosol mask with 22 mm elbow	25
OPN-1025	OxyPlus large aerosol mask with 22 mm elbow	25
OKN-1025	OxyKid pediatric aerosol mask with 22 mm elbow	25
OTN-1025	OxyTyke pediatric aerosol mask with 22 mm elbow	25

OxyTrach

The OxyTrach delivers 23–83% FiO₂ on flows from one to 15+ LPM flush. Oxygen is delivered through a standard small bore tube using the 22 mm OD adapter attached to a universal swivel elbow. The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O₂ delivery. The mask can be used with or without high humidity or aerosol therapy. It does not contain natural rubber latex or phthalates such as DEHP.

OxyTrach		
Catalog no.	Description	Qty. (cs.)
OAT-1025	OxyTrach mask, 7' O ₂ tubing	25

Bubble humidifiers and large volume nebulizers

Bubble humidifier

Prefilled with sterile water for inhalation, these humidifier systems feature a dual-purpose reservoir with a snap-open trigger, which makes a round opening for secure attachment of a humidifier adapter (*available both with Humidifier Systems 002620 and AL2702, and also separately*). The oxygen pathway is designed to help prevent particulate from reaching oxygen delivery tubing. An audible relief valve warns if the patient's tubing becomes occluded. They are designed for convenient shelf storage.

Our latex-free AirLife empty humidifiers provide an economical alternative to prefilled humidifiers. These 370 mL capacity units have an audible-alarm, pressure-relief valve to warn patients and caregivers of obstructions to flow. They are available with either a three or six psi pressure-relief valve.

Bubble humidifier		
Catalog no.	Description	Qty. (cs.)
002620	Prefilled kit 500 mL	12
AL2702	Prefilled kit 750 mL	12
AL2500	Prefilled 500 mL (<i>replacement water</i>)	12
AL2700	Prefilled 750 mL (<i>replacement water</i>)	12
ALP2002	Humidifier adapter only	120
002003	Empty 370 mL with 3 psi valve	50
002006	Empty 370 mL with 6 psi valve	50



002620 Compact prefilled humidifier with adapter

Empty humidifiers left 002003, right 002006

Prefilled nebulizer kits left CK0010, right CK0005

Large volume nebulizer

Our prefilled nebulizers provide performance and ease of use. Obtain heated and cool high-output aerosols adjustable from 28–100% FiO₂. This wide selection of solutions includes: sterile water for inhalation; USP; 0.45% sodium chloride for inhalation, USP; and 0.9% sodium chloride for inhalation, USP. Easy setup includes an aluminum, heat-exchange system and external slip-on heater (*sold separately*), which heats the aerosol while reducing the possibility of contamination.

Large volume nebulizer		
Catalog no.	Description	Qty. (cs.)
CK0005	Prefilled kit 500 mL, sterile water for inhalation	12
CK0010	Prefilled kit 1000 mL, sterile water for inhalation	12
CK4505	Prefilled kit 500 mL, 0.45% sodium chloride for inhalation	12
CK4510	Prefilled kit 1000 mL, 0.45% sodium chloride for inhalation	12
CK9010	Prefilled kit 1000 mL, 0.9% sodium chloride for inhalation	12
CN0005	Prefilled 500 mL (<i>replacement water</i>)	18
CN0010	Prefilled 1000 mL (<i>replacement water</i>)	12
CN4505	Prefilled 500 mL 0.45% sodium chloride (<i>replacement water</i>)	18
CN4510	Prefilled 1000 mL 0.45% sodium chloride (<i>replacement water</i>)	12
CN9010	Prefilled 1000 mL 0.9% sodium chloride (<i>replacement water</i>)	12
CC10	Nebulizer cap only (<i>1000 mL bottles only</i>)	24
2M8021	Nebulizer heater (<i>used with CK0005–CK9010</i>)	1
002002	Empty 350 mL (<i>35–100% FiO₂</i>)	24
5007P	Empty 350 mL (<i>28–98% FiO₂</i>)	50
5207	Empty 500 mL (<i>28–98% FiO₂</i>)	50

Misty Ox

Misty Ox high-flow nebulizers are designed for highly accurate, continuous aerosol therapy at flow rates up to 110 LPM to meet even the most challenging patient's inspiratory demands. They are available in an air entrainment design (441H) as well as a closed system design (441G) for quieter operation. These nebulizer adapters are compatible with most standard, pour-bottle solutions.

Misty Ox		
Catalog no.	Description	Qty. (cs.)
441G	Misty Ox gas injection nebulizer (GIN)	24
441H	Misty Ox high-flow nebulizer	24



441G Misty Ox gas injection nebulizer (GIN)



Oxygen tubing



Corrugated tubing

Oxygen tubing

Our extensive line of disposable vinyl-tipped oxygen tubing comes with crush-resistant lumen to resist occlusion. Choose from a variety of lengths ranging from 7' (84") to 50' (600"). Tubing 001350 also comes with a U/Connect-It connector to attach directly to flowmeters without an adapter.

AirLife disposable corrugated tubing for aerosol and IPPB applications offers flexibility and convenience to meet therapeutic needs. Cuffs on all tube ends lock to 7/8" (22 mm) connections. *AirLife* corrugated tubing does not contain DEHP.

Oxygen tubing		
Catalog no.	Description	Qty. (cs.)
001302	7' crush resistant	50
001350	7' crush resistant with U/Connect-It	50
001303	14' crush resistant	50
001305	25' crush resistant	25
001308	40' crush resistant	20
001306	50' crush resistant	15
001303GRN	14' crush resistant, green	50
001305GRN	25' crush resistant, green	25
001306GRN	50' crush resistant, green	15

Corrugated tubing		
Catalog no.	Description	Qty. (cs.)
001401	5'	50
001400	6'	50
001410	6' in 6" segments	50
001404	100' in 6" segments	1
001405	100' in 6" segments with dispenser box	1
001425	6"	50
001424	18"	50
001421	5' blue in 6" segments	50
001420	6' blue in 6" segments	50
001426	100' blue in 6" segments	1
001427	100' blue in 6" segments with dispenser box	1



For more information, call Customer Service at **800.323.9088**.

CareFusion
Vernon Hills, IL

carefusion.com



© 2016 CareFusion Corporation or one of its subsidiaries. All rights reserved. AirLife, U/Connect-It, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. All other trademarks are property of their respective owners. AL6181 (0516/500)