



EMERGENCY OXYGEN UNITS & CPR ADMINISTRATION EQUIPMENT

[Emergency Oxygen Units](#)
[Cylinders](#)
[Regulators](#)
[Masks](#)
[Home](#)
[AHA Emergency Oxygen Guidelines](#)
[Contact LIFE®](#)
[About Us](#)
[Refill Information](#)
[Maintenance](#)
[Regulations](#)
[Warranty](#)
[Go To LIFE® StartSystem](#)
[Go To LIFE® SoftPac™](#)


LIFE is Time™

When minutes count
EMERGENCY OXYGEN
can be the difference between
LIFE OR DEATH.

- **Simple and EASY TO USE,**
External Two-step Instructions and Clear Cover
Inhalator: administer supplemental oxygen for a breathing victim
Resuscitator: oxygen enriched mouth-to-mask CPR for a non-breathing victim
- **Easily Portable** with Convenient Handle
20" x 6" x 9" - 16lbs - Aluminum Cylinder
- **Wall Mount** for Recognition and Accessibility, always available on display for easy access
- **Simple "On-Off" Control Lever**
- **Constant Reading Supply Gauge**
(other units must be turned on to read gauge)
- **90 Minute Supply @ 6LPM - 566 liters Oxygen**
- **Resuscitation LIFE® CPR Mask**
with One-Way Valve fits both adult and child
- **Cylinder is easily disposable/replaceable**
or refillable at local gas distributors.
- **Tight Seal Protects Unit from Dust and Moisture**
- **Units Come Ready-to-Use**
full of oxygen with wall brackets and screws
- **5 Year Warranty**



*This unit
could save
your life.*



**Make LIFE® OxygenPac
part of your first aid and
safety program.**

[DEMO VIDEO](#)


LIFE® CPR Mask

invertably fits both ADULT and CHILD

**IDEAL LOCATIONS FOR THE LIFE®
OxygenPac
...WHEREVER THERE ARE PEOPLE!**

For Employees & First Responders
INDUSTRIAL/FACILITY
PUBLIC FACILITIES
ENTERTAINMENT / RECREATION OFFICES



LIFE® OxygenPac is attractive. Prominently displayed,
it makes the simple but strong statement,
"This Organization Cares For Its Employees And Visitors."

*Uncompromisingly manufactured with quality
materials, technology, and workmanship*

MODEL FEATURES

LIFE-612 6&12 LPM - "**NORM**" & "**HIGH**" (2 settings)
provides to recommended 100% inspired O2

LIFE-025 0 to 25 LPM Variable Flow (only for trained EMTs)
12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25

LIFE-6FF 6 LPM (Liters Per Minute) Fixed-Flow

LIFE-101 Replacement Cylinder Assembly
complete with cylinder, valve, gauge & O2 fill

Regulators

FDA Regulation minimum capability is 6 LPM flow rate.
LIFE® 612 models deliver *both* the *minimum* of 6 LPM
and the AHA recommended 100% inspired oxygen at 12 LPM,
in just two simple settings which read "NORM" & "HIGH".
Only LIFE Corporation offers the simple 6&12 LPM flow rates.

American Heart Association Guidelines for CPR Oxygen

- ▶ "... give 100% inspired oxygen during Basic Life Support
... and Advanced Cardiac Life Support ... (BLS & ACLS)
... as soon as it becomes available."
- ▶ "... use supplemental oxygen...at a flow rate of 10-12 LPM
... to enable delivery of 100% inspired oxygen" (initially).

www.LIFECorporation.com

[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!