



Petrolatum Gauze Wound Dressing

Gauze Wound Dressing Impregnated With Petrolatum

Petrolatum Gauze Wound Dressing is a non-adherent dressing made with cotton gauze impregnated with white petrolatum. Maintains a moist wound bed to support healing. Clings and conforms to all body contours. Use beneath an absorbent secondary dressing. Sterile. Not made with natural rubber latex.

Indications:

For use as a primary dressing on low to non-exudating wounds, including:

- donor sites
- burns
- skin graft sites
- lacerations
- abrasions

Product	Size	Reorder #	Units
Petrolatum Gauze Dressing	1x8	23180	50/Box
Petrolatum Gauze Dressing	3x9	23390	50/Box
Petrolatum Gauze Dressing	3x18	23318	12/Box
Petrolatum Gauze Dressing	6x36	23636	12/Box

DermaRite®